A microcomputer-based software package for eye monitoring research

GEORGE W. McCONKIE University of Illinois at Urbana-Champaign, Urbana, Illinois

> CHARLES W. SCOUTEN Stoelting Company, Chicago, Illinois

> > and

PATRICK K. BRYANT and JERRY WILSON Reading Technology, Inc., Savoy, Illinois

A software package is described that collects and reduces eye behavior (eye position and pupil size) data using an IBM-compatible personal computer. Written in the C language for speed and portability, the package includes several unique features: data can be collected simultaneously from other sources (e.g., EEG, EMG), logically defined events can be detected in real time on any data channel, and either of two types of data matrix can be produced. Data reduction algorithms and data structures are described.

The monitoring of eye behavior, including eye position and pupil size, has proven to be useful in studying issues in a number of different fields, including reading (Rayner, 1983), music reading (Sloboda, 1978), X-ray examination (Nodine & Kundel, 1987), map design (Antes, Chang, & Mullis, 1985), mathematics education (Suppes, Cohen, Laddaga, Anliker, & Floyd, 1982), chess (Chase & Simon, 1973), art and design (Lochegr & Nodine, 1987), mental load (Kahneman & Beatty, 1966), emotional response (Hess, 1972), and oculomotor control (Zuber, 1981). Almost any area of research that studies either the acquisition and use of visually presented information or changes in emotional response or arousal can benefit from the recording of eye position or pupil size.

Increased interest in monitoring the eyes has stimulated the development of new devices for monitoring eye behavior. However, a major factor that has prevented a more widespread use of eye monitoring has been the amount of software development required before research is possible. Each laboratory typically writes the software necessary to collect and reduce data with the particular equipment configuration being used.

Recent widespread acceptance of the IBM-initiated standard among personal computers, together with the continually increasing speed of computers of this family, suggests the appropriateness of developing a set of application programs for this computer family that can be used with a variety of eye behavior monitoring devices. The present paper describes a set of programs that was developed specifically for the Stoelting Wide Angle Eye Movement Tracker and TV Pupillometer System, but which is compatible with other eye monitoring equipment.

SOME BASIC CONCEPTS

The Stoelting tracker/pupillometer, a device that monitors the eyes by means of a video camera, contains hardware that processes each video frame to determine the size and location of the pupil. It samples eye position 60 times per second and gives as output three digital values: 7 bits each for horizontal and vertical eye position and 8 bits for pupil size. In addition, it has a field-of-view camera that allows the stimulus field to be displayed on a monitor, with a fixation spot that moves in response to movements of the eyes. In these respects, the equipment is quite similar to video-camera-based eye monitoring equipment produced by other companies, and is distinguished primarily by its relatively low price.

Each type of eye monitoring device currently available has strengths and weaknesses (Young & Sheena, 1975). Video-camera-based equipment tends to be relatively easy to use, gives an indication of both horizontal and vertical eye position, typically produces pupil size information as well as eye position, and ordinarily provides output showing the direction of gaze superimposed on a television image of the stimulus field. At the same time, most equipment of this type has a relatively slow sampling rate (i.e., 60 samples per second), accuracy to only about 1° of visual angle, and a noise level that prevents the reliable detection of small eye movements. Furthermore, there are occasionally artifacts in the signal that result from blinks or failures in the pupil-tracking algorithms that may not reflect actual eye movements or pupil size changes. These

Address correspondence to George W. McConkie, Center for the Study of Reading, University of Illinois at Urbana-Champaign, 51 Gerty Dr., Champaign, IL 61820.

characteristics make it impossible for equipment of this type to give a precise specification of the fine-grained characteristics of eye movements. For example, with samples taken only every 16 msec, it is not possible to specify the precise time at which a 30-msec saccadic eye movement began and ended or to describe its dynamic characteristics. Also, what appears to be a single eye fixation in the data may actually be a series of eye fixations separated by short saccades. However, since most types of research can be carried out without a higher degree of precision, the equipment is quite adequate for many purposes.

In developing software for video-camera-based equipment, it is necessary to deal with the problems of noise and indeterminacy in the data and to provide the researcher with as much information as possible about the eyes' behavior without falsely interpreting noise and artifacts as eye movements. The effects of the choice of values for parameters in the reduction of noisy eye movement data has been vividly described by Karsh and Breitenbach (1983). The goal of the present development was to deal with the problems posed by video-camera-based eye movement data in a way that would maximize the valid information available and, at the same time, would handle appropriately the signal characteristics of other types of eye monitoring devices, including those that sample at other rates and those that have other noise levels and characteristics. The data collection software has the flexibility to collect data from a variety of input ports on the computer, depending on the hardware configuration of the system. The data reduction software has a number of parameters that allow it to be "tuned" to properly handle eye movement data having different characteristics. Finally, diagnostic information is provided in the data matrix to indicate places in the data where noise characteristics are such that the validity of decisions being made is lower than normal. In this way, the more questionable data can be excluded from use in data analysis, if desired.

When eye position is monitored by a tracking device that samples very rapidly (e.g., every millisecond) and with low noise levels, there is seldom a question about whether a saccade occurred; indeterminacy of when it began and ended is within a range of a few milliseconds. Although the software described here can collect and reduce data in which fixations and saccades are always clearly specified, we discuss in this article the more difficult case in which there is greater indeterminacy. In this latter case, it is not appropriate to speak of fixations and saccades, since these cannot be determined precisely. Rather, we speak of a gaze, which can be defined as the period of time between the onset of one stable period in the eye movement data and the next stable period in which the eyes are directed toward a new location. Each gaze, by definition, has an initial stable period. The gaze may also have a period at the end that is not stable but that shows a transition toward the next gaze location. The amount of variability in the data that is still classified as

stable and the amount of change in gaze location that is required in order to reliably detect the fact that the eyes have shifted to a new gaze location are controlled by parameters in the program.

THE COMPUTER PROGRAMS

The eye monitoring software package presently consists of four programs: one for data collection, two for data reduction, and one that allows the user to modify calibration values in the data file. They run on IBMcompatible personal computers with 256K of memory and a hard disk. The programs are written primarily in the Microsoft C language, using the Microsoft C compiler to achieve both speed and portability. However, this language choice presented certain problems that are discussed below. Flexibility is achieved through parameter files that specify a number of the programs' internal characteristics.

Data Collection Program

The Data Collection Program is designed to collect eye position and pupil size data, as well as data of other kinds, including buttonpresses and analog data. Up to eight 16bit input ports can be recorded. The rate of sampling is determined externally. The program samples on the leading edge of a positive-going pulse that is generated by the eye monitoring device or produced by a clock. The program also allows the user to display stimulus patterns on the computer screen, advancing pages with keypresses or buttonpresses, or to use external displays, perhaps projecting slides or physically placing photographs in the field of view.

The user runs the Data Collection Program by typing its name plus the names of three files: a Parameter File, a Display Control File, and a Data File. The Parameter File and Display Control File are ASCII files that can be prepared and modified using any text editor capable of producing pure ASCII output. One additional file that contains the hardware initialization code is required. This file is used to establish communication between the computer and the eye monitoring equipment. For an interface, we currently use the Qua Tech PXB-721 Parallel Expansion Board, for which A/D modules are also available for monitoring analog signals. If different equipment is used, then this file can be altered to reflect the characteristics of the equipment.

Upon being run, the Data Collection Program requests the subject's name and identification number and an experimental condition number. All are necessary parameters, but shell scripts (i.e., DOS batch files) can be written that automatically specify some of these parameters. This procedure reduces the likelihood of experimeter error.

The parameters that are contained in the Parameter File are listed in Table 1. These parameters define the hardware environment, the channels on which data are to be collected, and information to be used in the calibration

Table 1 Data Collection Program Parameters: Parameter File
Number of data collection channels defined.
For each channel (0-7): Channel title Address (I/O Port) Data width (8 or 16 bits) Read on interrupt (1=enabled; 0=disabled)
Interrupt request line that is active.
Type of display (Internal or External, where Internal is a display on the computer monitor).
Number of calibration locations.
For each calibration location: X value in display space Y value in display space X location on computer screen Y location on computer screen
Calibration threshold on X data.
Calibration threshold on Y data.
Dianlay colors:

Display colors: Text color Background color Emphasized text color

Type of subject interaction during calibration (Keyboard or Event).

Type of subject interaction during data collection (Keyboard or Event).

process. This file is read and the parameters are invoked as part of the initialization process each time the program is run.

The Display Control File defines the sequence of events that are to occur in the experiment. By means of a series of dot commands, shown in Table 2, displays can be put onto the screen, data collection channels can be enabled and disabled, and calibrations can be taken as desired. Alphanumeric display pages can be included in the Display Control File and can use the IBM extended character set, thus giving some limited graphics capability.

The Data Collection Program contains a calibration procedure, which, when called in the internal display mode (i.e., when using the computer screen to display stimulus material), displays a grid of Xs on the computer screen at locations defined in the Parameter File. The color of each X is modified, in sequence, as a signal for the subject to fixate that point and press a key or button. When the key is pressed, the computer samples the eye data values, saving the values in a calibration table. This is repeated for each location, and the values obtained are compared with the prior values. If eye position values

Table 2

Data Collection Program: Dot Commands in Display Control File

deviate by more than the X or Y calibration threshold values specified in the parameter file, additional samples are taken from that location. This sampling continues until two successive samples taken from each location fall within the threshold values specified or until the experimenter manually exits the calibration procedure. This calibration procedure can be repeated at any time during the data collection process.

When the Data Collection Program is initially called, a nine-point matrix of digits appears on the screen while the experimenter adjusts the eye monitoring equipment. During this time, the values of data channels that are defined as eye position and pupil size are displayed continuously on the bottom of the computer screen. Also during this time, the computer identifies the range of values obtained from the pupil size channel, if available. Thresholds then are calculated for later use in identifying blinks. A keyboard command next advances the program to the actual data collection procedures, where calibrations and data collection can occur.

In the Data File that is produced by the Data Collection Program, the first 2048 bytes are used to store complete information about the data collection session (see Table 3), including data from each occurrence of the calibration task. Following this, a data record of variable length is created to contain the data collected during each display. This record indicates the channels on which data were collected during that display, an experimental condition identifier for that display, and the actual data, which consists of a 2-byte integer per sample for each channel.

One additional feature of the Data Collection Program is the ability to detect events and to take action when they occur. An event is the occurrence of a logical condition in the data coming from a specified data channel. An event is defined by an identifying title, a channel number, a threshold value, and a logical condition, including equal, not equal, greater than, less than, bitwise AND, and bitwise OR. Each defined event is tested as each data sample is taken. If an event tests as "true," then immediate action can be taken. Up to four events can be specified. At present, this feature is used to respond to the press of buttons on one of the data channels, which allows the system to advance the display or to simply record a subject's response. However, this feature provides the basis for a large number of potential response-contingent stimulus manipulations, especially since the events can be defined on channels other than the eye behavior channels. More detail concerning the techniques used to allow a C language program to deal with real-time events is presented below.

Data Reduction Programs

Two data reduction programs have been developed. One program is for situations in which the primary interest is in eye position data that are of a fixation-saccade type, and the other is for situations in which the interest is in nonsaccade data, such as pursuit eye movements, pupil

	Data Collection Program	: Records in the Data File	
System Information Record:	2048 bytes	Calibration Record:	1024 bytes
General Information:	128 bytes	Record type identifier	1 byte = 1
Record type identifier	120 bytes = 0	Size of record	2 bytes
21	4 bytes	Number of calibration points	1 byte
Record size		For each calibration point:	
Program name	64 bytes	Value assigned to X location	2 bytes
Version number: major	1 byte	Value assigned to Y location	2 bytes
Version number: minor	1 byte	For each calibration:	
Date/time	26 bytes	Preceding display number	1 byte
Reserved	31 bytes	For each calibration point in	
Subject Information:	128 bytes	each calibration:	
Subject name	64 bytes	Obtained X value	2 bytes
Subject ID number	2 bytes	Obtained Y value	2 bytes
Subject condition code	1 byte	Obtained pupil size	2 bytes
Reserved	61 bytes	Status of next calibration	1 byte
		(0=good, 1=aborted,	
File Information:	384 bytes	-1 = nonexistent)	
Parameter File name	64 bytes	Event Record:	384 bytes
Display Control File name	64 bytes	Record type identifier	1 byte = 3
Display condition code	1 byte	Record size	4 bytes
Data File name	64 bytes	Number of Events defined (≤ 4)	1 byte
Reserved	191 bytes	Definition of Event 0	68 bytes
Data Channel Information:	768 bytes	Event title	64 bytes
Number of channels (≤ 8)	1 byte	Channel number	1 byte
Channel 0 information	68 bytes	Threshold value	2 bytes
Channel title	64 bytes	Logical condition	1 byte
I/O Port	2 bytes	$\tilde{0} = equal$	-
Data width (8/16)	1 byte	1 = not equal	
Read on interrupt	i oyu	2 = greater than	
flag (1=true)	1 byte	3 = less than	
Channel 1 information	68 bytes	4=bitwise AND	
Chamer 1 mormation	08 bytes	5=bitwise OR	
Channel 7 information	68 bytes	Definition of Event 1	68 bytes
	223 bytes	Definition of Event 2	68 bytes
Reserved	225 bytes	Definition of Event 3	68 bytes
Data Collection Information:	128 bytes		
Data collection rate	2 bytes	Data Record: (length unspecified)	1 1-4-0
(samples per second)		Record type identifier	1 byte = 2
Number of ASCII characters used	l	Record size	4 bytes
to define page condition	1 byte	Display number	1 byte
Reserved	125 bytes	Channel mask (Channels	1 6
- ·	•	on which data are collected)	1 byte
Reserved:	512 bytes	Page condition	80 bytes
		Data (per channel per sample)	2 bytes

 Table 3

 Data Collection Program: Records in the Data File

size, or EEG. The first program, the Gaze-Oriented Reduction Program, produces a matrix in which each data row corresponds to a gaze or fixation-saccade pair. The second program, the Stream Data Reduction Program, is for use when the data cannot appropriately be partitioned into gaze segments. This program produces a row in the data matrix for every data sample taken. However, if eye position data were collected, then the matrix can include columns indicating eye position, the state of the eyes (i.e., stable or in transition), and the beginning and ending times of gazes.

Gaze-Oriented Data Reduction Program. Table 4 contains definitions and parameters used in the Gaze-Oriented Data Reduction Program, and Table 5 presents the general algorithm used in identifying gazes. The algorithm identifies periods of stability in the data, then examines the periods of instability lying between them, attempting to resolve them as either periods of noise in the

data or transition periods (i.e., eye movements taking the eyes to a new location). If the period of instability has characteristics that make it inappropriate to interpret in either of these simple ways, then a complex transition analysis is carried out that attempts to break the period into additional gazes using recursion with new parameters. The program produces an ASCII data matrix that contains a row of data for each identified gaze. The columns in the matrix are listed in Table 6. The matrix is in a form typically used with statistical programs and contains information about the location and duration of each gaze, the duration of the stable part of the gaze, and indications of data quality, such as the percentage of samples with nonaberrant values, amount of drift during the stable part of a gaze (a possible indicator of small, undetected saccades), and abnormal pupil size (an indicator of a blink).

Once a gaze has been identified, averages of the stable X and Y channel values are calculated. These averages

Term or Parameter*	Definition	Term or Parameter*	Definition	
X Channel	The data channel that gives information about horizontal eye position.	Reduced Local Gaze	A specified number of samples, fewer than the Local Gaze Duration, that meet the Local Gaze	
Y Channel	The data channel that gives information about vertical eye position.		Threshold Definition for a Local Gaze. The Reduced Local Gaze is used in Complex Tran- sition Analysis.	
Gaze	The eye movement record is used to segment time into discrete periods, each termed a Gaze. Each Gaze begins when the eyes stabilize in a new position.	Reduced Local Gaze Duration	The number of samples in a Reduced Local Gaze. (3)	
Stable Period	The initial period of a Gaze, during which the eyes are relatively still.	Aberrant Value	A data value from either the X or Y channel that differs from both the preceding and following data values by at least a specified amount, and	
Gaze Location	The location in display space to which the eyes are directed during a Stable Period.		in the same direction. Thus, it is a local peak or valley in the data stream, and is assumed to	
Gaze Duration	The number of samples included in a Stable Period.		represent noise rather than actual change in eye position.	
Transition Period	The final period of a Gaze, during which the eyes move to a new location.	Aberrant Value Threshold	The threshold that defines an aberrant data value. (2)	
Transition Length	The distance between the Gaze Locations that precede and follow the Transition Period, mea- sured in units appropriate to the stimulus display.	Transition Analysis	A procedure for examining Transition Data in order to determine whether or not a Transition Period exists.	
Transition Direction	The angle of the line connecting two successive Gaze Locations, measured in degrees, where zero deg is a horizontal line extending directly rightward from a Gaze Location. Locations above that location have positive values, and lo- cations below have negative values.	Minimum Gaze Duration	The minimum number of samples that can ex- ist between the initial samples of two succes- sive Local Gazes, without collapsing these into a single Gaze. Notice that this parameter is used in dealing with the case of overlapping Local Gazes. (3)	
Transition Duration	The number of samples included in a Transi- tion Period.	Successive Gaze Thresholds –	If the locations of two successive Gazes differ	
Local Gaze	A specified number of samples in which the maximum and minimum values on both the X and the Y channels differ by no more than the	X & Y	by less than the Successive Gaze Thresholds on both X and Y channels, they are collapsed into a single Gaze. (1)	
Local Gaze Duration	values defined as thresholds for those channels. The number of samples in a Local Gaze. (5)	Complex Transition Period	Transition Data that (a) contain more samples than specified by Transition Duration Threshold,	
Local Gaze Threshold X	The maximum difference that can exist among a set of X data values and still be considered a Local Gaze. (1)		or (b) have fewer samples than this threshold, are inconsistent (see definition above), and with preceding and following Gaze Locations that differ by more than the Successive Gaze Threshold. Transition Data not meeting either criterion are dealt with in the normal Transition Analysis.	
Local Gaze Threshold Y	The maximum difference that can exist among a set of Y data values and still be considered a Local Gaze. (1)			
Transition Data	A string of data values that does not contain a region that satisfies the definition of a Local Gaze. Thus, a candidate for containing a Tran- sition Period.	Transition Duration Threshold	A threshold used in determining whether a set of Transition Data constitutes a simple or com- plex transition period. (3)	
Consistent Transition Data	Transition Data values from one channel that change in only one direction.	Complex Transition Analysis	A procedure for further analyzing a Complex Transition Period to determine whether it can be broken up into additional Gazes. This proce-	
Inconsistent Transition Data	Transition Data values from one channel that both increase and decrease.	-	dure uses the Reduced Local Gaze in its attempt to find Stable Periods.	

Table 4 Gaze-Oriented Data Reduction Program: Definitions and Parameters

Note-Default values for parameters used with Stoelting Pupillometer are given in parentheses. *Parameters are typed in italic.

are used to determine where the eyes were directed in the stimulus display during that gaze. The algorithm used in this mapping process produces a solution similar to that of the linear interpolation procedure suggested by McConkie, Zola, Wolverton, and Burns (1978), but in fact accomplishes this goal using recursion.

The mapping algorithm first divides the stimulus space into quadrangles, based on the locations of the grid of target locations used in the calibration task. The algorithm assigns to each calibration target location the X and Y channel values obtained during the last calibration prior to the display, or the values obtained during the first calibration following it, or an average of the values obtained on these two calibrations. The algorithm then takes the average location of the eyes during the stable part of the first gaze and identifies the quadrangle within which it lies. Two lines are then identified that bisect opposing sides of the quadrangle and that divide it into four smaller quadrangles. Again, the algorithm identifies which of these quadrangles contains the gaze location. This process conGiven: A Local Gaze has been identified.

- A. Find the next series of data values that qualifies as a Local Gaze and carry out Transition Analysis.
 - 1. If the new Local Gaze begins one data value after the beginning of the prior Local Gaze, extend the current Gaze to include new data point. Go to A.
 - 2. If the new Local Gaze begins no more than a MINIMUM GAZE DURATION of values following the beginning of the prior Local Gaze, extend the current Gaze to include these data. Go to A. (Note: This includes cases where Local Gazes overlap substantially but are not sequential, indicating the possibility of some movement during the gaze, but too little to identify as saccades.)
 - 3. Else, carry out Transition Analysis to examine the Transition Data.
 - a. If the difference between either the average X or Y data values from the preceding and following stable periods exceeds the GAZE MOVEMENT THRESHOLD, then set the Movement Flag = TRUE, else set it FALSE. (Set a flag if the eyes have changed location.)
 - b. Set the Transition Duration Code to indicate the duration of the Transition Period. If Transition Duration (TD)=0, then set Transition Duration Code (TDC)=0; if TD=1, then set TDC=1; if TD is greater than 1 but less than the Transition Duration Threshold, then set TDC=2; if TD is greater than or equal to the Transition Duration Threshold, then set TD=3.
 - c. If the Transition Data from either the X or the Y channel are Inconsistent, then set the Consistency Flag = FALSE, else set it TRUE. (Use a flag to indicate whether or not the data show the eyes proceeding consistently in one direction.)
 - d. If the Movement Flag = TRUE and the Consistency Flag = TRUE and the Transition Duration Code < 3, then classify this as a clean transition, close out the old Gaze, and go to A. (A short, consistent transition that takes the eyes to a new location are accepted as a saccadic movement.)
 - e. If the Transition Duration Code = 3 or if the Movement Flag = TRUE and the Transition Duration Code = 2 and the Con-

tinues until the gaze location is within a specified distance of one of the corners of a quadrangle, at which point that corner is taken as the gaze location. Since the algorithm knows the corresponding location of that corner in the stimulus field, the mapping of the data point onto a stimulus location is complete and the stimulus-based gaze location is entered into the matrix. The distance threshold used in this process is a parameter in the program. The gaze location, together with the distance and direction of this gaze from prior and following gazes, are entered into the data matrix.

Special provisions are made in the program to handle a gaze that lies outside the region in which calibration targets were used. This is done by extending the nearest quadrangle to include the data point, estimating X and Yvalues corresponding to the new corners, and then using the normal mapping algorithm to estimate the location of the gaze on the stimulus.

If an event occurred during a gaze, the matrix also contains information about which event occurred and the time of its occurrence. This can be used to record response sistency Flag = FALSE, then go to 4 and carry out Complex Transition Analysis. (If this is a long transition or if the direction of movement is not consistent, it needs additional analysis.)

- f. If the Transition Duration Code < 3 and the Movement Flag = FALSE and the Consistency Flag = FALSE, then mark all Aberrant Values, combine the preceding and following Stable Periods into one Gaze, and go to A. (If this is a period of inconsistent movement with the eyes ending up in the same place, conclude that it is simply a period of noise and include it as part of the gaze.)
- g. If the Transition Duration Code < 3 and the Movement Flag = FALSE and the Consistency Flag = TRUE, then combine the preceding and following Stable Periods into one Gaze, and go to A. (This is a condition in which there is a period of consistent movement that leaves the eyes in the same location. It is included for completeness, but should not normally occur.)
- h. If the Transition Duration Code < 3 and the Movement Flag = TRUE and the Consistency Flag = FALSE, then mark the Aberrant Values. If an Aberrant Value is adjacent to a Stable Period, then join it to that Period, removing it from the Transition Period. The remaining Transition Period is accepted as a legitimate transition. Go to A. (If the eyes have moved, but the transition period is noisy, extend the Stable Period if possible and accept the transition.)
- B. Carry out Complex Transition Analysis.
 - Mark all Aberrant Values in Transition Data. Consider them to be equal to their bounding values, and reapply Transition Analysis. If this eliminates all Transition Data, go to A. (The transition was apparently simply an extended noise period.)
 - 2. Else use Reduced Local Gaze definition and reapply Transition Analysis to remaining Transition Data. Indicate in the Problem Flag when data show unusual characteristics. Go to A. (This call will find shorter stable periods or periods during which there is more drift or movement, but which do not show the characteristics of saccadic movements. Thus, it has the capability of decomposing a long string of Transition Data into additional gazes.)

times and/or response selection if each event defined is the press of a particular button.

The Gaze-Oriented Data Reduction Program also puts summary information concerning data from non-eyeposition channels into the data matrix, including mean value on that channel during the eye fixation, as well as an indication of the direction and amount of change during the gaze period. Finally, the sample number of the first data sample in the gaze is included. This number provides a pointer from the gaze-based matrix into the original data matrix. Thus, it is possible to write programs that select the data collected on some non-eye-movement channel, starting with the onset of a particular gaze or event.

Stream Data Reduction Program. When the researcher's interest is primarily in changes in pupil size or some other continuous data, such as EEG or smooth pursuit eye movement data, the data matrix for analysis is prepared using the Stream Data Reduction Program. This program creates an ASCII matrix having one row for each sample taken. Subject and condition columns are

Table 6					
Gaze-Oriented Data	Reduction	Program:	Data	Matrix	

Col.	Bytes	Contents	
1.	4	Subject number	
2.	2	Subject condition code	
3.	3	Display number	
4.	As defined	Display condition code	
5.	3	Gaze number on display	
6.	6	Eye position on display: X	
7.	6	Eye position on display: Y	
8.	3	Percent of values that are nonaberrant: X	
9.	3	Percent of values that are nonaberrant: Y	
10.	5	Distance from prior Gaze	
11.	6	Direction from prior Gaze (degrees)	
12.	5	Distance to following Gaze	
13.	6	Direction to following Gaze (degrees)	
14.	4	Total number of samples in Gaze	
15.	4	Number of samples in stable part of	
		Gaze	
16.	7	Problem Flag	
17.	5	Sample number of first sample in Gaze	
18.	4	Number of aberrant samples	
19.	5	Average X value in Gaze	
20.	5	Average Y value in Gaze	
21.	4	Change in X during Gaze	
22.	4	Change in Y during Gaze	
For each Event, when it occurs:			
E1.	2	Event number	
E2.	4	Time since beginning of Gaze	
E3.	5	Time since beginning of Display	
For each Data Channel other than X and Y Eye Position:			
D1.	5	Mean data value during Gaze	
D2.	4	Maximum data value during Gaze	
D3.	4	Minimum data value during Gaze	
D4.	4	Mode data value during Gaze	
D5.	5	Slope of data during Gaze	

added to the data. If eye position has also been recorded, then columns can be created that indicate for each sample whether it occurred during a period of stability or transition, how many samples have been taken since the beginning of the gaze, where the eyes were directed during that gaze, and whether there was any problem in the eye position data at that time (e.g., blinks or aberrant values). The occurrence of events can also be noted. With this matrix it is easy to identify the data taken on one channel beginning at the onset of a fixation at a particular location in the display, or beginning at the point where the subject pressed a button or the pupil size exceeded a particular value. Columns in this matrix are listed in Table 7.

Calibration Data Editor

A fourth program included in the software package allows the experimenter to examine the calibration values contained in a data file. Furthermore, if it is found that a value is incorrect, then the program allows the experimenter to replace the original value with a new one. This is occasionally necessary when a value obtained in one calibration is very deviant from those obtained during other calibrations in the same session.

USING THE C LANGUAGE FOR REAL-TIME PROGRAMMING

The C Language that is available on the IBM PC was not designed for use in real-time data collection. In particular, it does not provide a method of responding to most interrupts. Attempting to use it to respond directly to interrupts can actually interfere with the proper execution of some of its standard commands. This fact created certain problems in developing the Data Collection Program, particularly the need for interrupt-driven data collection and the need to respond quickly to the occurrence of an event. This problem was dealt with by developing an input channel handler, written in assembly language. This routine receives commands from the C language program and, in turn, makes information available to it that can be read upon request.

The input channel handler receives instructions from the calling program concerning which channels to take data from, whether to read 8 or 16 bits of data, which port addresses to use, where the data are found, which interrupt request line to respond to, and what the definitions of events are. When activated, the handler then responds to each interrupt by sampling the data from each specified input port and storing it in a buffer. Furthermore, as each sample is taken, that sample is tested against any event definition specified for that channel. If the test succeeds, then a bit is set in a status word to indicate that the event occurred.

The C language program can read the 16-bit status word and a 16-bit error word generated by the input channel

Table 7 Stream Data Reduction Program: Data Matrix			
Col.	Bytes	Contents	
1.	4	Subject number	
2.	2	Subject condition code	
3.	3	Display number	
	As defined	Display condition code	
5.	5	Data sample number	
	F	for Each Channel (0-7):	
C1.	5	Channel data value	
		Gaze Information:	
G1.	3	Gaze number on display	
G2.	6	Eye position on display: X	
G3.	6	Eye position on display: Y	
G4.	4	Total number of samples in Gaze	
G5.	4	Number of samples in stable part of Gaze	
G6.	4	Sample number in the Gaze	
G7.	7	Problem flag for gaze	
G8.	4	Percent of samples that are aberrant	
G9.	3	Blink flag	
G10.	2	Aberrant value flag: X	
G11.	2	Aberrant value flag: Y	
	For	Each Defined Event (0-3):	
E1.	2	Event type	
E2.	4	Number of samples since Gaze beginning	
E3.	5	Number of samples since display be- ginning	

handler and can request values from the data buffer. For simply obtaining data, the program can request data from the handler at its own convenience, so long as this occurs frequently enough to avoid overrunning the buffer. However, to respond quickly to events, it is necessary for the C language program to check the status flag frequently to see if an event has occurred. Because the program was written so that it checks this status at least every millisecond, it can respond to an event within 1 msec of the event's occurrence.

CONCLUSION

Monitoring eye behavior can play any of three different roles in research. First, it provides data for studying oculomotor control, perception, cognition, emotional response, or arousal. Second, eye behavior data can serve as a basis for analyzing simultaneously collected data, thus allowing selection, elimination, or grouping of those data. Third, the eyes can be used to control aspects of the stimulus; for instance, the display can be changed when the eyes are directed at a certain location or a tone can be sounded when the size of the pupil exceeds a certain value.

The software package described above makes it possible to use an IBM-compatible personal computer to acquire and reduce eye behavior data for any of the three purposes just described. Although developed specifically for the Stoelting Wide Angle Eye Movement Tracker and TV Pupillometer System, the package was designed to be compatible with a variety of different eye monitoring devices that sample at different rates and have different noise characteristics. The availability of this software should substantially reduce the time and cost involved in establishing the facilities to obtain eye behavior data. One unique contribution of the software is the ability to simultaneously collect data of other types, such as other physiological measures, buttonpresses, or indications of changes taking place in the stimulus field. The second contribution is a general facility for detecting logically defined events in the data in real time. This capability can be used

for making response-contingent stimulus changes, although this currently requires making user-specific modifications to the existing software. These two characteristics open the door to new and innovative roles for eye monitoring in psychological research.

REFERENCES

- ANTES, J. R., CHANG, K., & MULLIS, C. (1985). The visual effect of map design: An eye-movement analysis. *The American Cartographer*, 12, 143-155.
- CHASE, W. G., & SIMON, H. A. (1973). The mind's eye in chess. In W. G. Chase (Ed.), *Visual information processing* (pp. 215-281). New York: Academic Press.
- HESS, E. H. (1972). Pupillometrics: A method of studying mental, emotional, and sensory processes. In N. S. Greenfield & R. A. Sternbach (Eds.), *Handbook of psychophysiology* (pp. 491-531). New York: Holt, Rinehart & Winston.
- KAHNEMAN, D., & BEATTY, J. (1966). Pupil diameter and load on memory. Science, 154, 1583-1585.
- KARSH, R., & BREITENBACH, F. W. (1983). Looking at looking: The amorphous fixation measure. In R. Groner, C. Menz, D. F. Fisher, & R. A. Monty (Eds.), *Eye movements and psychological functions: International views* (pp. 53-64). Hillsdale, NJ: Erlbaum.
- LOCHEGR, P., & NODINE, C. (1987). Symmetry catches the eye. In J. K. O'Regan & A. Levy-Schoen (Eds.), *Eye movements: From physiology to cognition* (pp. 353-361). Amsterdam: North Holland.
- MCCONKIE, G. W., ZOLA, D., WOLVERTON, G. S., & BURNS, D. D. (1978). Eye movement contingent display control in studying reading. Behavior Research Methods & Instrumentation, 10, 154-166.
- NODINE, C. F., & KUNDEL, H. L. (1987). The cognitive side of visual search in radiology. In J. K. O'Regan & A. Levy-Schoen (Eds.), *Eye movements: From physiology to cognition* (pp. 578-583). Amsterdam: North Holland.
- RAYNER, K. (Ed.). (1983). Eye movements in reading: Perceptual and language processes. New York: Academic Press.
- SLOBODA, J. A. (1978). The psychology of music reading. Psychology of Music, 6, 3-20.
- SUPPES, P., COHEN, M., LADDAGA, R., ANLIKER, J., & FLOYD, R. (1982). Research on eye movements in arithmetic performance. In R. Groner & P. Fraisse (Eds.), *Cognition and eye movements* (pp. 57-73). Amsterdam: North Holland.
- YOUNG, L. R., & SHEENA, D. (1975). Survey of eye movement recording techniques. Behavior Research Methods & Instrumentation, 7, 397-429.
- ZUBER, B. L. (Ed.). (1981). Models of oculomotor behavior and control. Boca Raton, FL: CRC Press.