

# A Statistically Validated Electrochemical (Conductometric) and Spectroscopical Study of Some Metal Ligand Complexes for Prospective Biological Action

Anindya Bagchi<sup>1\*</sup>, Anusree Raha<sup>1</sup>, Prosenjit Mukherjee<sup>1</sup>, Monit Pal<sup>1</sup>, Kunal Datta<sup>1</sup>, Arnab Goswami<sup>2</sup>, Paramita Karmakar<sup>1</sup>, Diptendu Sekhar Biswas<sup>2</sup>, Reshmi Mukherjee<sup>2</sup>

<sup>1</sup>Netaji Subhas Chandra Bose Institute of Pharmacy, Chakdaha, West Bengal, India

<sup>2</sup>Bengal College of Pharmaceutical Technology, Birbhum, West Bengal, India

DOI: [10.36348/sjimps.2023.v09i02.006](https://doi.org/10.36348/sjimps.2023.v09i02.006)

| Received: 06.01.2023 | Accepted: 11.02.2023 | Published: 19.02.2023

\*Corresponding author: Anindya Bagchi

Netaji Subhas Chandra Bose Institute of Pharmacy, Chakdaha, West Bengal, India

## Abstract

The present electrochemical study deals with the evaluation of a complexation process involving different metal and ligands and its physicochemical properties to find out the ligand- metal ratio of complex in solution. For the determination of complexing nature “Monovariation method” have been used to ascertain the ligand metal ratio in the complex. The stability constant of the formed complexes were calculated by conductance measurement using Modified Job’s method (Turner Anderson Method). The analysis had been carried out by using conductometric principle and the final outcome of the experiment may influenced by temperature and with the ionic concentration of metal and ligand solution. Free energy change values were determined to find out the feasibility of the complexation process. One of the method was validated statistically by using system and method precision parameters. Lastly the complexes was been assayed by using a specific chromogenic agent with a specific spectroscopic technique that may give biological action in future.

**Keywords:** Conductometry; Monovariation method; Modified Job’s method; Chelation; Metal ion.

**Copyright © 2023 The Author(s):** This is an open-access article distributed under the terms of the Creative Commons Attribution 4.0 International License (CC BY-NC 4.0) which permits unrestricted use, distribution, and reproduction in any medium for non-commercial use provided the original author and source are credited.


## INTRODUCTION

### Complexometric Titration

In complexometric titration the complexes consist of one or more molecules bound to a *central cation*. The molecules bound to the central cation are called *Ligand*. The maximum number of groups bound to an ion is called the *coordination number* of the ion. Coordination number concept comes because of the *coordinate bond* present between the Ligand and metal ion in the structure of complexes as the metal ion (Lewis Acid) accepts the lone pair of electron and the Ligand (Lewis Base) molecule donates it. Thus in a ligand molecule there is presence of at least one lone pair of electrons through which co-ordinate linkage takes place with the metal ion. Importantly electron has two types i.e. bonding pair of electrons and lone pair of electrons. Bonding pair of electrons is accountable for bond formation between the atoms (single bond)

whereas lone pair will be accountable for the formation of coordinate bond as it will be not accountable for the bond formation as the bonding pair of electrons. E.g. coordinate bond formation between  $\text{NH}_3$  and  $\text{H}^+$  ion. Ligands those are bound to the ion at only one point (because of having and attachment of a single functional group which are sharing only one lone pair of electron) are called *unidentate* (one toothed) two points are called *bidentate* (two teathed because a ligand is sharing two lone pair of electron) and so on. Many Ligand contains more than one group capable of bonding to the metal ion with more number of lone pair of electron as such ligands are called *multidentate*.

So denticity only depends on the number of lone pair of electron shared by the ligand with the metal ion. A single ligand can be unidentate or multidentate (e.g. bidentate, hexadentate etc.).



**Fig 1: Structure of Ammonia and  $\text{Cu}^{++}$  ion complex**

### Electrochemistry

It is a branch of chemistry that deals with the chemical changes produced by electricity and the production of electricity by chemical changes. These reactions involve electrical charge moving between electrodes and electrolytes.

### Principle of Conductometric Titration

The electrical current through a chemical cell is carried out by the ionic species in the solution conductometric. The ease with which current is conducted through a solution (under the influence of potential difference applied across two electrodes) depends up on the concentrations and kind of ions in the solution. If two suitable electrodes are present in a solution and potential difference is applied across those electrodes then current will flow through the solution. During progress of a titration changes in the conductivity of the solution usually occur and at the end point involving neutralization or precipitation reaction the conductivity of the solution will be minimum. The principle of conductometric titration is based on the fact that during the titration, one of the ions is replaced by the other and invariably these two ions differ in the ionic conductivity with the result that conductivity of the solution varies during the course of titration. The equivalence point may be located graphically by plotting the change in conductance as a function of the volume of titrant added.

**Conductance:** Solutions of electrolytes normally obey Ohm's law.

$$R = V/I.$$

R = resistance of the solution.

V = applied potential difference.

I = current through the solution.

The conductance (G) of a solution is reciprocal of resistance.

$G = 1/R$ . The unit of conductance is *Siemens (S)*. The conductivity of a solution is measured by an instrument known as conductometer.

**Conductivity:** If we think the solution in between two electrodes as a slab of material then the conductance of that material is found to be directly proportional to the area of cross section (A) of the solution and inversely proportional to the length (L) between the electrodes. Therefore  $G = K \times L/A$  where K = conductivity of the material = specific conductance of the material. Therefore  $K = G \times L/A$ .

So conductivity of a material is defined as the conductance of a cube of material having a surface area of  $1\text{m}^2$  and length of 1m. Unit of conductivity is S/m or S/cm or  $\mu\text{S/cm}$ .

Practically Conductivity = Conductance X Cell constant of that particular experimental conductometric cell or electrode which is constant for that particular electrode and independent for any electrolytic solution but it depends only on the temperature of the electrolytic solution whom the cell will measure the conductance.

### Validation Parameter

#### Precision

Precision is the measurement of how close the data values to each other for a number of measurements under the same analytical conditions. Precision may be considered at three levels according to ICH.

#### Repeatability

#### System Precision

Precision under same operative conditions (within a laboratory over a short period of time using the same analyst with the same equipment) was determined. Mean, SD and % RSD were calculated from data. The system precision is checked by using standard chemical substance to ensure that the

analytical system is working properly. In this retention time and area of six determinations is measured and % RSD should be calculated.

**Acceptance criteria:** % RSD should be in between 98%- 102%.

#### Method Precision

In method precision, a homogenous sample of single batch should be analysed 6 times. This indicates whether a method is giving consistent results for a single batch. In this analysis the sample has been analysed six times with the calculation of % RSD.

**Acceptance criteria:** % RSD should be in between 98%-102%.

#### Reproducibility

Precision between laboratories/intermediate precision can be considered during the standardization of a procedure before it is submitted to the pharmacopoeia. A simple logic behind this parameter was some degree of inconsistency (Occurrence of random error which is a type of error that is basically untraceable or beyond the control of analyst) was allowed for every analytical measurement. But, the extent depends on steps involved (Weighing, dilution etc.), technique used in other expected variables (Stability) and intended use of the procedure.

#### Intermediate Precision (Ruggedness)

Precision under different laboratory conditions (within-laboratory variation, as on different days, or with different analysts, or equipment within the same laboratory) has been carried out. This experiment is done to see whether the test results are coming under the acceptance criteria or not.

**Acceptance criteria:** % RSD should be in between 98%-102%.

#### Robustness

Here the closeness of the values are seen in small changes of different parameters like solvent, temperature, pH etc. Here the mean, SD, % RSD is calculated.

**Acceptance criteria:** % RSD should be in between 98%-102% [1].

#### Stability Constant

Stability or formation or binding constant is the type of equilibrium constant used for the formation of metal complexes in the solution. Acutely, stability constant is applicable to measure the strength of interactions between the ligands and metal ions that are involved in complex formation in the solution [2]. A generally stability constant equations are expressed as the following ways: equilibrium constants, and these equilibrium constants are known as overall stability

constants or overall formation [3]. Any metal complexes will be of greater stability if its stability constant has the higher value. Sometimes the  $1/k$  values are alternative values of stability constant, and now this is called as instability constant.

#### Thermodynamic stability

In a chemical reaction, chemical equilibrium is a state in which the concentration of reactants and products does not change over time. Often this condition occurs when the speed of forward reaction becomes the same as the speed of reverse reaction. It is worth noting that the velocities of the forward and backward reaction are not zero at this stage but are equal. For the formation of metal chelates, the thermodynamic technique provides a very significant information. Thermodynamics is a very useful technique in distinguishing between enthalpic effects and entropic effects. The bond strengths are totally effected by enthalpic effect, and this does not make any difference in the whole solution in order/disorder. Based on thermodynamics the chelate effect below can be best explained. The change of standard Gibbs free energy for equilibrium constant is response:

Free energy change:

$$\Delta G = -2.303 RT \log K$$

R = gas constant

T = absolute temperature

Log K = Stability Constant For metal complexes

Thermodynamic stability and kinetic stability are two interpretations of the stability constant in the solution. If reaction moves from reactants to products, it refers to a change in its energy as shown in the above equation. But for the reactivity, kinetic stability is responsible for this system, and this refers to ligand species [4].

#### Factors affecting the stability constant of metal complexes

Equilibrium concentration of the free metal ion

Temperature

Analytical concentration of the ligand.

pH of the solution

## MATERIAL AND METHODS

### Chemicals

All the chemicals used were of Analytical Grade. Ferric Chloride, Ferrous Chloride, Ferrous Sulphate (Lobachem, India, 0.1 M, 0.05M) and EDTA, Oxalic acid, PABA, (Merck, India, 0.1M, 0.05M) solution were prepared by using double distilled water in the 100ml volumetric flask for several times and has been used in subsequent times for different methods. Pure Curcumin (Lobachem, India, 0.1M, and 0.05M) and 1, 10-Phenanthroline (Lobachem, India) solution was also prepared for several types in this experimental process.

### Calibration of conductivity meter

A Systronics model 306 CONDUCTIVITY METER with Conductivity Cell and a simple weight machine from EAGLE and a BOROSILICATE GLASS. In electrode method conductance reading was noted which having the unit called Siemens.

0.7 gm KCl was dissolved in 100 ml of distilled water to prepare 0.1M solution and conductance reading was noted in Milli Siemens (mS) unit to calibrate the instrument.

**Table 1: Observation found out in calibration process**

| Concentration at 25 <sup>0</sup> C | Conductance at 25 <sup>0</sup> C (mS) |
|------------------------------------|---------------------------------------|
| 0.1 M KCl                          | 12.88                                 |
| 0.01 M KCl                         | 1.413                                 |
| 0.001 M KCl                        | 0.146                                 |

### Monovariation method (Mole ratio method/Yoe-Jones method)


This method was introduced by Yoe and Jones [5]. A series of solutions are prepared in which the total concentration of the metal is kept constant and the concentration of the ligand is varied under similar condition. A plot was prepared of conductance value as a function of the ratio of moles of ligand to moles of the metal. The corresponding point where end point was visualized as the occurrence of decrease or increase (lowest or highest) in conductance value followed by the increase or decrease in value gave the idea about molar metal: ligand complexation process.

### Conductometric titration for detection of Metal-Ligand ratio (Monovariant method) [6]

Metal salt solution having strength 0.1M was prepared with distilled water. Similarly, 0.1 M of ligand solution was prepared by using the same method as the previous one. Metal salt solution was taken in a beaker and this was titrated conductometrically against 0.1M ligand solution that was taken in a burette.

### Modified Job's method (Method of Continuous Variations)

The modification of the Job's [7] continuous variation method performed by Vesburgh and Cooper [8] was applied to find the stoichiometric and formation constant of the complex. The case of co-ordination may be described by:



A series of solutions are prepared in which the sum of total concentration of M and L is kept constant but their proportions are continuously varied. The conductance values of the series is plotted against the mole fraction of the ligand. The ratio of the stoichiometric coefficients is determined from the mole fraction at the point of decrease or increase (lowest or highest) in conductance value followed by increase or decrease takes place. That value corresponding to the mole metal: ligand ratio will be the endpoint of the titration.

Equimolar (0.1M) solutions of ligand and metal solutions were prepared and three series C1, C2,

C3 of solutions were made. In set C1 metal salt solution was filled with volume 0.0ml to 12.0ml and total volume was made to 12.0ml in each. Similarly, in C2 ligand solution was filled and set C3 was prepared by mixing metal salt solution from 0.0ml to 12.0ml and ligand solution from 12.0ml to 0.0ml. Conductance values was recorded for each solution.  $\Delta$  Conductance was calculated as "C1+C2-C3" [9]. Graphs was plotted between Conductance values and mole metal-ligand ratio. The composition and stability constant can be determined from the equivalence point in the graph [10, 11].

### Turner Anderson Method for determination of stability constant

Turner and Anderson [12] have modified Job's method and have successfully used for determination of stability constant. By plotting a continuous variation curve for a given range of compositions and then repeating the procedure for more dilute solutions.

Stability Constant equation according to Turner Anderson Method is-

$$x = \frac{a_2b_2 - a_1b_1}{(a_2 - a_1) + (b_2 - b_1)} \dots \dots \dots (1)$$

$$k = \frac{x}{(a-x)(b-x)} \dots \dots \dots (2)$$

k= stability constant, (a-x) and (b-x) are the concentration of (product-reactant) in case of ligand and metal if calculated from Modified Job's Method.

x= concentration of complex

a,b= initial concentration of metal and ligand (0.1M).

a<sub>1</sub>, b<sub>1</sub>= concentration of metal and complex for 0.1M metal and ligand solution.

a<sub>2</sub>, b<sub>2</sub>= concentration of metal and complex for 0.05M metal and ligand solution.

After evaluating the value of x from equation no. 2, it can be substituted in equation no. 1 to find out the stability constant value.

### Free energy change

Free energy change of all the complexes that were formed was evaluated by the following equation:

$$\Delta G = -2.303 RT \log K$$

R = gas constant

T = absolute temperature

Log K= Stability Constant For metal complexes

### Assay of Metal- Chromogenic agent

In the present experiment a chromogenic agent (1, 10- phenanthroline) had been added which will indicate the metal content present by imparting colour in free metal state. This concept can help for the experiment in-vivo in different disease state i.e. Thalassemia where excess haemoglobin staged up inside the body and excess ferrous ion has to be removed from the body to cure the disease. In that case a ligand molecule can bind with that ferrous ion of haemoglobin as ligand-metal complex and can be removed from the body. So that concept can be experimented in-vitro by plotting a standard curve of

metal-chromogenic agent complex followed by taking the metal-ligand complex solution after reaching the endpoint (Using Monovariation Method) and find out the concentration of the solution with the help of spectroscopic principle. Importantly that metal-ligand solution was added after addition of the chromogenic agent that will give the colour to the resulting solution and will indicate the metal content in the solution.

## RESULT AND DISCUSSION


M= Metal, S= Solvent, L= Ligand

**Ferric Chloride vs EDTA**  
Monovariation method

**Table 2: Monovariation Method**

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 0                      | 25.6        |
| 1                      | 26.2        |
| 2                      | 26.9        |
| 3                      | 27.5        |
| 4                      | 27.8        |
| 5                      | 28.8        |
| 6                      | 29.3        |
| 7                      | 30.9        |
| 8                      | 32          |
| 9                      | 32.6        |
| 10                     | 32.7        |
| 11                     | 33.3        |
| 12                     | 33.7        |
| 13                     | 35          |
| 14                     | 35.8        |
| 15                     | 36          |
| 16                     | 37.2        |
| 17                     | 37.8        |
| 18                     | 38.2        |
| 19                     | 39          |
| 20                     | 39.5        |
| 21                     | 40.4        |
| 22                     | 41          |
| 23                     | 41.8        |
| 24                     | 41.9        |
| 25                     | 42          |
| 26                     | 42.1        |
| 27                     | 42.2        |
| 28                     | 42.3        |
| 29                     | 42.7        |
| 30                     | 43.2        |
| 31                     | 43.6        |
| 32                     | 43.7        |
| 33                     | 43.9        |
| 34                     | 41.8        |
| 35                     | 41.5        |
| 36                     | 40.2        |
| 37                     | 38.5        |
| 38                     | 36.6        |
| 39                     | 34.9        |

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 40                     | 33.7        |
| 41                     | 32.9        |
| 42                     | 31.8        |
| 43                     | 31.3        |
| 44                     | 30.2        |
| 45                     | 28.6        |
| 46                     | 27.8        |
| 47                     | 25.9        |
| 48                     | 25.1        |
| 49                     | 24.8        |
| 50                     | 24.5        |



Graph 1: Monovariation Method

$$VIS1 = V2S2$$

V1S1 = Volume and Strength of EDTA

V2S2 = Volume and strength of Ferric Chloride

$$33XS1 = 50X0.1$$


$$S1 = 50X0.1/33 = 0.151 \text{ M}$$

So approx., Metal: Ligand Ratio = .1: .151.

Modified Job's Method

Table 3: Modified Job's Method

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 2.41    | 2.04    | 0.285   | 4.165    |
| 01:11 | 3.41    | 1.36    | 1.92    | 2.85     |
| 02:10 | 7.9     | 1.41    | 0.405   | 8.905    |
| 03:09 | 8.5     | 0.319   | 1.29    | 7.529    |
| 04:08 | 9.25    | 0.232   | 1.45    | 8.032    |
| 05:07 | 3.08    | 1.14    | 6.67    | -2.45    |
| 06:06 | 4.3     | 1.53    | 9.99    | -4.16    |
| 07:05 | 2.09    | 0.482   | 8.09    | -5.518   |
| 08:04 | 3.51    | 0.475   | 6.66    | -2.675   |
| 09:03 | 6.35    | 0.512   | 2.08    | 4.782    |
| 10:02 | 4.66    | 0.36    | 2.56    | 2.46     |
| 11:01 | 6.6     | 7.23    | 1.64    | 12.19    |
| 12:00 | 6.9     | 7.66    | 1.84    | 12.72    |



**Graph 2: Modified Job's Method**

The end point was found out at 7:5 of Metal: Ligand ratio.

Turner Anderson Method


**Table 4: Turner Anderson Method (0.1M)**

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 6.73    | 13.9    | 7.15    | 13.48    |
| 01:11 | 4.23    | 12.6    | 2.73    | 14.1     |
| 02:10 | 7.14    | 11.8    | 4.74    | 14.2     |
| 03:09 | 9.65    | 11.1    | 6.35    | 14.4     |
| 04:08 | 13.3    | 10      | 7.1     | 16.2     |
| 05:07 | 15.5    | 8.48    | 1.08    | 22.9     |
| 06:06 | 17.8    | 8.33    | -2.78   | 28.91    |
| 07:05 | 19.9    | 7.09    | -13.41  | 40.4     |
| 08:04 | 22.3    | 6.2     | -8.1    | 36.6     |
| 09:03 | 24.4    | 5.21    | -3.99   | 33.6     |
| 10:02 | 26.4    | 4.17    | 1.07    | 29.5     |
| 11:01 | 28.2    | 3.11    | 4.11    | 27.2     |
| 12:00 | 30.3    | 0.748   | 5.348   | 25.7     |

**Table 5: Turner Anderson Method (0.05M)**

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.399   | 7.5     | 7.55    | 0.349    |
| 01:11 | 7.18    | 7.01    | 3.19    | 11       |
| 02:10 | 6.8     | 7.17    | 4.49    | 9.48     |
| 03:09 | 6.42    | 7.97    | 4.65    | 9.74     |
| 04:08 | 5.93    | 8.81    | 4.83    | 9.91     |
| 05:07 | 4.86    | 10.8    | 3.01    | 12.65    |
| 06:06 | 3.97    | 14.2    | -0.03   | 18.2     |
| 07:05 | 2.43    | 22.7    | -8.07   | 33.2     |
| 08:04 | 2.2     | 21.9    | -7.2    | 31.3     |
| 09:03 | 1.98    | 20.6    | -4.62   | 27.2     |
| 10:02 | 1.47    | 18.3    | -1.73   | 21.5     |
| 11:01 | 1.12    | 17.5    | -0.18   | 18.8     |
| 12:00 | 0.41    | 16.7    | 1.01    | 16.1     |



**Graph 3: Turner Anderson Method**

End point were observed at 7:5, for 0.1M and for 0.05M also.


Ferric Chloride vs Oxalic acid  
Monovariation method

**Table 6: Monovariation method**

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 0                      | 25.2        |
| 1                      | 25.8        |
| 2                      | 26.4        |
| 3                      | 26.7        |
| 4                      | 26.7        |
| 5                      | 26.7        |
| 6                      | 26.7        |
| 7                      | 27.4        |
| 8                      | 27.4        |
| 9                      | 28          |
| 10                     | 28.5        |
| 11                     | 28.5        |
| 12                     | 29.3        |
| 13                     | 29.5        |
| 14                     | 30.4        |
| 15                     | 30.6        |
| 16                     | 31.7        |
| 17                     | 32          |
| 18                     | 32.3        |
| 19                     | 33.3        |
| 20                     | 33.4        |
| 21                     | 34.2        |
| 22                     | 35          |
| 23                     | 35.4        |
| 24                     | 36          |
| 25                     | 36.7        |
| 26                     | 36.8        |
| 27                     | 37          |
| 28                     | 37.8        |
| 29                     | 38.1        |
| 30                     | 38.8        |
| 31                     | 39.1        |
| 32                     | 39.6        |



| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 33                     | 40.4        |
| 34                     | 40.6        |
| 35                     | 41.1        |
| 37                     | 40.1        |
| 38                     | 39.8        |
| 39                     | 39.7        |
| 40                     | 39.6        |
| 41                     | 39          |
| 42                     | 38.3        |
| 43                     | 38.1        |
| 44                     | 37          |
| 45                     | 36.6        |
| 46                     | 36.1        |
| 47                     | 36          |
| 48                     | 35.5        |
| 49                     | 35.2        |
| 50                     | 35.1        |



**Graph 4: Monovariation method**

$$VIS1 = V2S2$$

V1S1= Volume and Strength of Oxalic Acid

V2S2= Volume and strength of Ferric Chloride

$$35XS1=50X0.1$$


$$S1 = 50X0.1/35 = 0.142 \text{ M}$$

So approx., Metal: Ligand Ratio = .1:.142.

Modified Job's Method

**Table 7: Modified Job's Method**

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.517   | 24      | 23.4    | 1.117    |
| 01:11 | 3.4     | 22.5    | 30.5    | -4.6     |
| 02:10 | 6.76    | 21.1    | 33.3    | -5.44    |
| 03:09 | 9.67    | 20.1    | 36.4    | -6.63    |
| 04:08 | 12.2    | 18.2    | 41.8    | -11.4    |
| 05:07 | 14.4    | 16.4    | 43.8    | -13      |
| 06:06 | 16.1    | 15.2    | 46.7    | -15.4    |
| 07:05 | 18      | 12.7    | 41.4    | -10.7    |
| 08:04 | 19.6    | 11.1    | 34.9    | -4.2     |
| 09:03 | 21.4    | 9.17    | 34.3    | -3.73    |
| 10:02 | 22.6    | 6.9     | 31      | -1.5     |
| 11:01 | 24.1    | 3.67    | 29.2    | -1.43    |
| 12:00 | 26.3    | 5.19    | 28.2    | 3.29     |



Graph 5: Modified Job's Method

The end point was found out at 6:6 of Metal: Ligand ratio


Turner Anderson Method

Table 8: Turner Anderson Method 0.1M

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.885   | 24.1    | 24.1    | 0.885    |
| 01:11 | 4.05    | 23.1    | 29.3    | -2.15    |
| 02:10 | 4.76    | 21.3    | 35      | -8.94    |
| 03:09 | 6.33    | 19.9    | 39.3    | -13.07   |
| 04:08 | 7.98    | 18.1    | 43.7    | -17.62   |
| 05:07 | 8.67    | 16.7    | 45.6    | -20.23   |
| 06:06 | 9.84    | 15.2    | 43.4    | -18.36   |
| 07:05 | 11      | 14.8    | 39.5    | -13.7    |
| 08:04 | 12      | 11      | 35      | -12      |
| 09:03 | 13.1    | 9.04    | 28.3    | -6.16    |
| 10:02 | 13.9    | 6.33    | 27      | -6.77    |
| 11:01 | 14.9    | 4.03    | 25.1    | -6.17    |
| 12:00 | 15.7    | 0.25    | 24      | -8.05    |

Table 9: Turner Anderson Method 0.05 M

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.46    | 15.2    | 17.4    | -1.74    |
| 01:11 | 2.67    | 18.2    | 23.6    | -2.73    |
| 02:10 | 4.56    | 11.2    | 20.9    | -5.14    |
| 03:09 | 6.22    | 10.3    | 23.4    | -6.88    |
| 04:08 | 6.83    | 9.36    | 24.7    | -8.51    |
| 05:07 | 9.3     | 8.56    | 24.6    | -6.74    |
| 06:06 | 10.5    | 7.67    | 22.7    | -4.53    |
| 07:05 | 11.7    | 6.7     | 20.7    | -2.3     |
| 08:04 | 13.1    | 5.28    | 19.5    | -1.12    |
| 09:03 | 14.1    | 4.31    | 18.1    | 0.31     |
| 10:02 | 15.8    | 3.14    | 17      | 1.94     |
| 11:01 | 16.1    | 2.14    | 16.2    | 2.04     |
| 12:00 | 17.6    | 0.24    | 16.3    | 1.54     |

**Graph 6: Turner Anderson Method**


End point were observed at 7:5, for 0.1M and for 0.05M also

Ferric Chloride vs Curcumin  
Monovariation Method

**Table 10: Monovariation Method**

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 0                      | 22.9        |
| 1                      | 23          |
| 2                      | 22.2        |
| 3                      | 21.2        |
| 4                      | 20.2        |
| 5                      | 19.5        |
| 6                      | 19          |
| 7                      | 18.3        |
| 8                      | 17.7        |
| 9                      | 17.1        |
| 10                     | 16.6        |
| 11                     | 16          |
| 12                     | 15.6        |
| 13                     | 15.1        |
| 14                     | 14.6        |
| 15                     | 14.6        |
| 16                     | 14.1        |
| 17                     | 14.1        |
| 18                     | 13.3        |
| 19                     | 13.3        |
| 20                     | 13.1        |
| 21                     | 12.7        |
| 22                     | 12.5        |
| 23                     | 12          |
| 24                     | 11.9        |
| 25                     | 12.7        |
| 26                     | 12.8        |
| 27                     | 13.2        |
| 28                     | 13.4        |
| 29                     | 13.8        |
| 30                     | 14.6        |
| 31                     | 14.7        |
| 32                     | 15.3        |
| 33                     | 15.4        |
| 34                     | 16          |

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 35                     | 16.2        |
| 36                     | 16.4        |
| 37                     | 16.6        |
| 38                     | 16.81       |
| 39                     | 16.91       |
| 40                     | 17.75       |
| 41                     | 18          |
| 42                     | 18.62       |
| 43                     | 19.5        |
| 44                     | 19.36       |
| 45                     | 19.43       |
| 46                     | 19.48       |
| 47                     | 19.57       |
| 48                     | 20.21       |
| 49                     | 20.22       |
| 50                     | 20.29       |



Graph 7: Monovariation Method

VIS1= V2S2

V1S1= Volume and Strength of Curcumin

V2S2= Volume and strength of Ferric Chloride

24XS1=50X0.1


S1= 50X0.1/24= 0.208 M

So approx., Metal: Ligand Ratio = .1:0.208.

Modified Job's Method

Table 11: Modified Job's Method

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.393   | 8.13    | 1.41    | 7.113    |
| 01:11 | 3.79    | 8.5     | 1.43    | 10.86    |
| 02:10 | 6.45    | 8.66    | 2.81    | 12.3     |
| 03:09 | 8.73    | 10.2    | 5.26    | 13.67    |
| 04:08 | 7.9     | 11.3    | 8.03    | 11.17    |
| 05:07 | 12.8    | 12.1    | 8.46    | 16.44    |
| 06:06 | 14.6    | 14.1    | 9.67    | 19.03    |
| 07:05 | 16.7    | 16.4    | 10.9    | 22.2     |
| 08:04 | 18.2    | 18.9    | 13.3    | 23.8     |
| 09:03 | 19.9    | 22.5    | 14.7    | 27.7     |
| 10:02 | 21.3    | 29.6    | 17.6    | 33.3     |
| 11:01 | 23      | 39.2    | 19.9    | 42.3     |
| 12:00 | 24.4    | 47.1    | 24.6    | 46.9     |



Graph 8: Modified Job's Method

The end point was found out at 4:8 of Metal: Ligand ratio

Turner Anderson Method


Table 12: Turner Anderson Method 0.1M

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.392   | 8.11    | 7.4     | 1.102    |
| 01:11 | 3.77    | 10.4    | 1.4     | 12.77    |
| 02:10 | 8.76    | 7.76    | 1.82    | 14.7     |
| 03:09 | 9.73    | 10.4    | 5.26    | 14.87    |
| 04:08 | 10.9    | 11.3    | 10.03   | 12.17    |
| 05:07 | 13.8    | 12.1    | 8.56    | 17.34    |
| 06:06 | 16.6    | 13.1    | 9.57    | 20.13    |
| 07:05 | 19.7    | 15.4    | 9.9     | 25.2     |
| 08:04 | 23.2    | 16.9    | 12.3    | 27.8     |
| 09:03 | 29.9    | 21.5    | 13.7    | 37.7     |
| 10:02 | 31.1    | 29.6    | 15.6    | 45.1     |
| 11:01 | 33      | 35.2    | 19.6    | 48.6     |
| 12:00 | 34.3    | 47      | 23.6    | 57.7     |

Turner Anderson Method

Table 13: Turner Anderson Method 0.05 M

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.387   | 8.1     | 7.3     | 1.187    |
| 01:11 | 3.65    | 9.23    | 6.1     | 6.78     |
| 02:10 | 5.93    | 10.21   | 9.21    | 6.93     |
| 03:09 | 7.86    | 12.1    | 11.2    | 8.76     |
| 04:08 | 6.34    | 11.3    | 10      | 7.64     |
| 05:07 | 10.3    | 11.9    | 10.2    | 12       |
| 06:06 | 11.1    | 12      | 10.3    | 12.8     |
| 07:05 | 12.2    | 13.5    | 10.9    | 14.8     |
| 08:04 | 12.9    | 14.1    | 11.3    | 15.7     |
| 09:03 | 16.3    | 15.3    | 11.9    | 19.7     |
| 10:02 | 19.8    | 16      | 12.1    | 23.7     |
| 11:01 | 23.6    | 17      | 12.6    | 28       |
| 12:00 | 29.3    | 18.6    | 13.9    | 34       |




Graph 9: Turner Anderson Method

PABA vs Ferric Chloride  
Monovariation Method

Table 14: Monovariation Method

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 0                      | 14.1        |
| 1                      | 14.4        |
| 2                      | 14.2        |
| 3                      | 14.3        |
| 4                      | 14.7        |
| 5                      | 15.2        |
| 6                      | 15.3        |
| 7                      | 15.9        |
| 8                      | 16          |
| 9                      | 16.3        |
| 10                     | 16.4        |
| 11                     | 16.6        |
| 12                     | 17.1        |
| 13                     | 17.8        |
| 14                     | 18.4        |
| 15                     | 18.6        |
| 16                     | 19.3        |
| 17                     | 19.8        |
| 18                     | 19.9        |
| 19                     | 21          |
| 20                     | 21.6        |
| 21                     | 21.8        |
| 22                     | 21.9        |
| 23                     | 21.2        |
| 24                     | 21          |
| 25                     | 19.95       |
| 26                     | 19.9        |
| 27                     | 19.78       |
| 28                     | 19.71       |
| 29                     | 19.64       |
| 30                     | 19.54       |
| 31                     | 19.39       |
| 32                     | 19.22       |

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 33                     | 19.14       |
| 34                     | 19.13       |
| 35                     | 19.03       |
| 37                     | 18.88       |
| 38                     | 18.83       |
| 39                     | 18.79       |
| 40                     | 18.7        |
| 41                     | 18.63       |
| 42                     | 18.54       |
| 43                     | 18.42       |
| 44                     | 18.4        |
| 45                     | 18.3        |
| 46                     | 18.32       |
| 47                     | 18.11       |
| 48                     | 18.1        |
| 49                     | 18.03       |
| 50                     | 17.98       |



Graph 10: Monovariation Method

$$VIS1 = V2S2$$

V1S1 = Volume and Strength of PABA

V2S2 = Volume and strength of Ferric Chloride

$$22XS1 = 50X0.1$$

$$S1 = 50X0.1/22 = 0.227 \text{ M}$$


So approx., Metal: Ligand Ratio = .1:0.227

Modified Job's Method

Table 15: Modified Job's Method

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.627   | 0.74    | 0.553   | 0.814    |
| 01:11 | 2.83    | 0.749   | 2.69    | 0.889    |
| 02:10 | 4.63    | 0.716   | 4.41    | 0.936    |
| 03:09 | 5.77    | 0.706   | 5.5     | 0.976    |
| 04:08 | 6.73    | 0.702   | 6.22    | 1.212    |
| 05:07 | 7.93    | 0.677   | 7.51    | 1.097    |
| 06:06 | 7.93    | 0.666   | 7.61    | 0.986    |
| 07:05 | 9.44    | 0.656   | 9.51    | 0.586    |
| 08:04 | 10.71   | 0.616   | 10.5    | 0.826    |
| 09:03 | 11      | 0.604   | 10.9    | 0.704    |
| 10:02 | 12.6    | 0.563   | 11.8    | 1.363    |
| 11:01 | 13      | 0.581   | 12      | 1.581    |
| 12:00 | 14.9    | 0.127   | 13.7    | 1.327    |





Graph 11: Modified Job's Method

The end point was found out at 4:8 of Metal: Ligand ratio

Turner Anderson Method


Table 16: Turner Anderson Method 0.1M

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.68    | 0.8     | 0.083   | 1.397    |
| 01:11 | 2.69    | 0.88    | 0.343   | 3.227    |
| 02:10 | 4.81    | 0.103   | 0.519   | 4.394    |
| 03:09 | 5.34    | 0.114   | 0.891   | 4.563    |
| 04:08 | 6.42    | 0.141   | 1.35    | 5.211    |
| 05:07 | 6.64    | 0.173   | 1.89    | 4.923    |
| 06:06 | 6.77    | 0.212   | 2.64    | 4.342    |
| 07:05 | 9.23    | 0.224   | 3.6     | 5.854    |
| 08:04 | 10.4    | 0.267   | 4.66    | 6.007    |
| 09:03 | 11.9    | 0.313   | 5.84    | 6.373    |
| 10:02 | 12.7    | 0.383   | 7.79    | 5.293    |
| 11:01 | 13.9    | 0.47    | 10.5    | 3.87     |
| 12:00 | 14.2    | 0.628   | 13.4    | 1.428    |

Turner Anderson Method

Table 17: Turner Anderson Method 0.05M

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.654   | 0.039   | 0.037   | 0.656    |
| 01:11 | 1.55    | 0.041   | 0.156   | 1.435    |
| 02:10 | 2.35    | 0.06    | 0.176   | 2.234    |
| 03:09 | 3.19    | 0.069   | 0.18    | 3.079    |
| 04:08 | 3.9     | 0.088   | 0.709   | 3.279    |
| 05:07 | 3.65    | 0.09    | 0.738   | 3.002    |
| 06:06 | 2.96    | 1.12    | 1.61    | 2.47     |
| 07:05 | 6.13    | 0.189   | 1.72    | 4.599    |
| 08:04 | 6.57    | 0.214   | 2.51    | 4.274    |
| 09:03 | 7.35    | 0.288   | 3.96    | 3.678    |
| 10:02 | 7.91    | 0.341   | 4.96    | 3.291    |
| 11:01 | 8.36    | 0.461   | 6.67    | 2.151    |
| 12:00 | 8.7     | 0.663   | 8.5     | 0.863    |


**Graph 12: Turner Anderson Method**

Ferrous Sulphate Vs EDTA  
Monovariation Method

**Table 18: Monovariation Method**

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 0                      | 3.61        |
| 1                      | 3.5         |
| 2                      | 3.45        |
| 3                      | 3.38        |
| 4                      | 3.32        |
| 5                      | 3.12        |
| 6                      | 3.09        |
| 7                      | 3           |
| 8                      | 2.98        |
| 9                      | 2.92        |
| 10                     | 2.9         |
| 11                     | 2.87        |
| 12                     | 2.81        |
| 13                     | 2.77        |
| 14                     | 2.74        |
| 15                     | 2.68        |
| 16                     | 2.63        |
| 17                     | 2.6         |
| 18                     | 2.55        |
| 19                     | 2.52        |
| 20                     | 2.47        |
| 21                     | 2.4         |
| 22                     | 2.38        |
| 23                     | 2.32        |
| 24                     | 2.31        |
| 25                     | 2.26        |
| 26                     | 2.22        |
| 27                     | 2.21        |
| 28                     | 2.18        |
| 29                     | 2.14        |
| 30                     | 2.1         |
| 31                     | 2.06        |
| 32                     | 2.02        |
| 33                     | 2.22        |
| 34                     | 2.27        |

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 35                     | 2.31        |
| 36                     | 2.45        |
| 37                     | 2.47        |
| 38                     | 2.51        |
| 39                     | 2.59        |
| 40                     | 2.65        |
| 41                     | 2.81        |
| 42                     | 3.8         |
| 43                     | 3.85        |
| 44                     | 3.92        |
| 45                     | 4.02        |
| 46                     | 4.06        |
| 47                     | 4.15        |
| 48                     | 4.22        |
| 49                     | 4.29        |
| 50                     | 4.4         |



Graph 13: Monovariation Method


 $VIS1 = V2S2$ 
 $V1S1 = \text{Volume and Strength of EDTA}$ 
 $V2S2 = \text{Volume and strength of Ferrous Sulphate}$ 
 $32XS1 = 50X0.1$ 
 $S1 = 50X0.1/32 = 0.156 \text{ M}$ 

So approx., Metal: Ligand Ratio = .1.:156.

Modified Job's Method

Table 19: Modified Job's Method

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.75    | 12.3    | 12.6    | 0.45     |
| 01:11 | 0.83    | 11.8    | 12.1    | 0.53     |
| 02:10 | 1.15    | 11.1    | 11.6    | 0.65     |
| 03:09 | 0.92    | 10.4    | 10.8    | 0.52     |
| 04:08 | 0.72    | 9.48    | 9.74    | 0.46     |
| 05:07 | 1.17    | 8.64    | 9.73    | 0.08     |
| 06:06 | 0.92    | 7.49    | 8.38    | 0.03     |
| 07:05 | 2.8     | 6.74    | 7.27    | 2.27     |
| 08:04 | 3.09    | 5.43    | 7.16    | 1.36     |
| 09:03 | 2.62    | 4.37    | 6.86    | 0.13     |
| 10:02 | 2.8     | 3.29    | 6.47    | -0.38    |
| 11:01 | 3.08    | 2.07    | 7.82    | -2.67    |
| 12:00 | 3.3     | 0.7     | 4.52    | -0.52    |



Graph 14: Modified Job's Method

The end point was found out at 6:6 of Metal: Ligand ratio Turner Anderson method


Table 20: Turner Anderson method 0.1M

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.743   | 11.6    | 12.3    | 0.043    |
| 01:11 | 0.81    | 11.71   | 12      | 0.52     |
| 02:10 | 0.796   | 12.43   | 12.8    | 0.426    |
| 03:09 | 0.979   | 9.66    | 10.6    | 0.039    |
| 04:08 | 1.2     | 8.71    | 10.2    | -0.29    |
| 05:07 | 1.51    | 7.62    | 9.31    | -0.18    |
| 06:06 | 1.71    | 6.93    | 8.17    | 0.47     |
| 07:05 | 1.95    | 6.22    | 7.36    | 0.81     |
| 08:04 | 2.11    | 5.11    | 7.21    | 0.01     |
| 09:03 | 2.46    | 4.4     | 6.69    | 0.17     |
| 10:02 | 2.61    | 2.86    | 6.41    | -0.94    |
| 11:01 | 2.77    | 1.91    | 7.71    | -3.03    |
| 12:00 | 2.93    | 0.8     | 4.43    | -0.7     |

Turner Anderson method

Table 21: Turner Anderson method 0.05M

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.757   | 12.1    | 11.1    | 1.757    |
| 01:11 | 0.707   | 11.7    | 10.8    | 1.607    |
| 02:10 | 0.731   | 11.3    | 10.4    | 1.631    |
| 03:09 | 0.865   | 10.6    | 9.8     | 1.665    |
| 04:08 | 0.885   | 9.8     | 9.3     | 1.385    |
| 05:07 | 0.903   | 8.7     | 8.6     | 1.003    |
| 06:06 | 0.924   | 8.3     | 8.5     | 0.724    |
| 07:05 | 1.16    | 7.9     | 7.8     | 1.26     |
| 08:04 | 1.29    | 7.4     | 7.1     | 1.59     |
| 09:03 | 1.43    | 6.9     | 5.4     | 2.93     |
| 10:02 | 1.58    | 6.6     | 5.1     | 3.08     |
| 11:01 | 1.73    | 4.19    | 5.1     | 0.82     |
| 12:00 | 1.87    | 3.16    | 4.86    | 0.17     |




Graph 15: Turner Anderson method

Oxalic acid vs Ferrous Chloride  
Monovariation Method

Table 22: Monovariation Method

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 0                      | 14.3        |
| 1                      | 14.2        |
| 2                      | 14.1        |
| 3                      | 13.8        |
| 4                      | 13.5        |
| 5                      | 13.2        |
| 6                      | 13          |
| 7                      | 12.8        |
| 8                      | 12.5        |
| 9                      | 12.3        |
| 10                     | 12.1        |
| 11                     | 11.8        |
| 12                     | 11.5        |
| 13                     | 11.3        |
| 14                     | 11.2        |
| 15                     | 11          |
| 16                     | 10.8        |
| 17                     | 10.6        |
| 18                     | 10.4        |
| 19                     | 10.2        |
| 20                     | 10          |
| 21                     | 9.8         |
| 22                     | 9.6         |
| 23                     | 9.4         |
| 24                     | 9.2         |
| 25                     | 9.1         |
| 26                     | 9           |
| 27                     | 8.7         |
| 28                     | 8.5         |
| 29                     | 8.4         |
| 30                     | 8.2         |
| 31                     | 9.6         |
| 32                     | 9.7         |
| 33                     | 10          |
| 34                     | 10.2        |
| 35                     | 10.3        |

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 36                     | 10.4        |
| 37                     | 10.6        |
| 38                     | 10.7        |
| 39                     | 11          |
| 40                     | 11.2        |
| 41                     | 11.6        |
| 42                     | 12          |
| 43                     | 12.1        |
| 44                     | 12.3        |
| 45                     | 12.5        |
| 46                     | 12.9        |
| 47                     | 13.1        |
| 48                     | 13.2        |
| 49                     | 13.5        |
| 50                     | 13.8        |



**Graph 16: Monovariation Method**

$VIS1 = V2S2$

V1S1= Volume and Strength of Oxalic Acid

V2S2= Volume and strength of Ferrous Chloride

$30XS1 = 50X0.1$


$S1 = 50X0.1/30 = 0.166 M$

So approx., Metal: Ligand Ratio = .1: .166.

Modified Job's Method

**Table 23: Modified Job's Method**

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.69    | 22.3    | 21      | 1.99     |
| 01:11 | 1.63    | 21.3    | 22.7    | 0.23     |
| 02:10 | 2.63    | 20.1    | 23.9    | -1.17    |
| 03:09 | 3.03    | 18.6    | 24.8    | -3.17    |
| 04:08 | 5.06    | 17.3    | 25.9    | -3.54    |
| 05:07 | 6.29    | 15.4    | 26      | -4.31    |
| 06:06 | 7.94    | 13.5    | 26.5    | -5.06    |
| 07:05 | 13.16   | 11.2    | 26.8    | -2.44    |
| 08:04 | 10.4    | 8.78    | 25.8    | -6.62    |
| 09:03 | 13      | 5.47    | 25.7    | -7.23    |
| 10:02 | 13.2    | 3.19    | 25.6    | -9.21    |
| 11:01 | 11.4    | 1.49    | 18.3    | -5.41    |
| 12:00 | 15.2    | 0.517   | 17.2    | -1.483   |



Graph 17: Modified Job's Method

The end point was found out at 6:6 of Metal: Ligand ratio  
Turner Anderson Method

Table 24: Turner Anderson Method 0.1M


| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.581   | 21.7    | 20.8    | 1.481    |
| 01:11 | 1.61    | 21      | 21.6    | 1.01     |
| 02:10 | 2.78    | 20.4    | 22.9    | 0.28     |
| 03:09 | 5.31    | 19.3    | 24.2    | 0.41     |
| 04:08 | 6.32    | 18.6    | 25.6    | -0.68    |
| 05:07 | 7.11    | 16.6    | 26.8    | -3.09    |
| 06:06 | 8       | 14.1    | 26.9    | -4.8     |
| 07:05 | 9.29    | 12.9    | 24.6    | -2.41    |
| 08:04 | 10.36   | 9.91    | 24.1    | -3.83    |
| 09:03 | 11.98   | 7.46    | 22.3    | -2.86    |
| 10:02 | 13.16   | 4.81    | 20.6    | -2.63    |
| 11:01 | 14.3    | 1.36    | 19.1    | -3.44    |
| 12:00 | 14.9    | 0.523   | 18.2    | -2.777   |

Turner Anderson Method

Table 25: Turner Anderson Method 0.05M

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 0.675   | 12.8    | 12.1    | 1.375    |
| 01:11 | 1.18    | 12.9    | 12.9    | 1.18     |
| 02:10 | 1.62    | 11.61   | 13.4    | -0.17    |
| 03:09 | 2.11    | 10      | 13.6    | -1.49    |
| 04:08 | 3.1     | 9.37    | 13.6    | -1.13    |
| 05:07 | 3.57    | 8.83    | 13.7    | -1.3     |
| 06:06 | 4.13    | 6.83    | 13.5    | -2.54    |
| 07:05 | 4.77    | 5.96    | 13.2    | -2.47    |
| 08:04 | 5.44    | 5.19    | 12.7    | -2.07    |
| 09:03 | 6.26    | 4.31    | 11.8    | -1.23    |
| 10:02 | 7.44    | 3.7     | 10.93   | 0.21     |
| 11:01 | 8.09    | 1.92    | 9.56    | 0.45     |
| 12:00 | 8.39    | 0.591   | 8.99    | -0.009   |





Graph 18: Turner Anderson Method

Stability Constant

Table 26: Stability constant Calculation

| Metal- Chelate Complex          | Stability Constant Value (logK) (at 29.8 <sup>0</sup> C) |
|---------------------------------|--|
| 1. Ferric Chloride- EDTA        | 1.28072381   |
| 2. Ferric Chloride- Oxalic Acid | 1.99215164   |
| 3. Ferric Chloride- Curcumin    | 1.17248713   |
| 4. Ferric Chloride- PABA        | 1.172487   |
| 5. Ferrous Sulphate-EDTA        | 1.17425422   |
| 6. Ferrous Sulphate-PABA        | 1.172487   |

Free energy Change

Free energy change:  $\Delta G = -2.303 RT \log K$ 

Table 27: Free energy change calculation

| Metal- Chelate Complex          | Free Energy Change (Kcal/mole) |
|---------------------------------|--------------------------------|
| 1. Ferric Chloride- EDTA        | -6694.96                       |
| 2. Ferric Chloride- Oxalic Acid | -10413.9                       |
| 3. Ferric Chloride- Curcumin    | -6129.16                       |
| 4. Ferric Chloride- PABA        | -6129.16                       |
| 5. Ferrous Sulphate-EDTA        | -6138.4                        |
| 6. Ferrous Sulphate-PABA        | -6129.16                       |

Assay of Metal-Chromogenic agent


Ferrous Sulphate- 1, 10- Phenanthroline

Monovariation method

Table 28: Monovariation method

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 0                      | 28.7        |
| 1                      | 28.3        |
| 2                      | 28.2        |
| 3                      | 28          |
| 4                      | 27.5        |
| 5                      | 27.4        |
| 6                      | 27.2        |
| 7                      | 27.1        |
| 8                      | 27          |
| 9                      | 26.8        |
| 10                     | 26.7        |

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 11                     | 26.4        |
| 12                     | 26.2        |
| 13                     | 26.1        |
| 14                     | 26          |
| 15                     | 25.4        |
| 16                     | 25.3        |
| 17                     | 25.1        |
| 18                     | 25          |
| 19                     | 24.3        |
| 20                     | 24.1        |
| 21                     | 24          |
| 22                     | 25.1        |
| 23                     | 25.4        |
| 24                     | 26.2        |
| 25                     | 26.3        |
| 26                     | 26.9        |
| 27                     | 27          |
| 28                     | 27.1        |
| 29                     | 27.2        |
| 30                     | 27.4        |
| 31                     | 27.5        |
| 32                     | 27.8        |
| 33                     | 27.9        |
| 34                     | 28.1        |
| 35                     | 28.2        |
| 37                     | 28.4        |
| 38                     | 28.5        |
| 39                     | 28.6        |
| 40                     | 29          |
| 41                     | 29.2        |
| 42                     | 29.3        |
| 43                     | 29.5        |
| 44                     | 29.7        |
| 45                     | 30          |
| 46                     | 30.2        |
| 47                     | 30.4        |
| 48                     | 30.7        |
| 49                     | 31          |
| 50                     | 31.2        |




**Graph 19: Monovariation method**

## Modified Job's Method

**Table 29: Modified Job's Method**

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 5.23    | 15.2    | 18.3    | 2.13     |
| 01:11 | 8.42    | 16.6    | 23.2    | 1.82     |
| 02:10 | 12.2    | 11.9    | 22.6    | 1.5      |
| 03:09 | 13.11   | 13.4    | 25.1    | 1.41     |
| 04:08 | 16      | 13.1    | 28      | 1.1      |
| 05:07 | 16.3    | 14      | 26.2    | 4.1      |
| 06:06 | 17.32   | 14.2    | 26.4    | 5.12     |
| 07:05 | 18.6    | 14.5    | 26      | 7.1      |
| 08:04 | 19.2    | 15      | 25.3    | 8.9      |
| 09:03 | 20.2    | 15.3    | 25.1    | 10.4     |
| 10:02 | 23.3    | 16      | 25      | 14.3     |
| 11:01 | 27.2    | 16.32   | 28.1    | 15.42    |
| 12:00 | 28      | 17      | 26.2    | 18.8     |

**Graph 20: Modified Job's Method**

The end point was found out at 4:8 of Metal: Ligand ratio

Turner Anderson Method


**Table 30: Turner Anderson Method 0.1M**

| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 5.23    | 15.2    | 18.3    | 2.13     |
| 01:11 | 8.42    | 16.6    | 23.2    | 1.82     |
| 02:10 | 12.2    | 11.9    | 22.6    | 1.5      |
| 03:09 | 13.11   | 13.4    | 25.1    | 1.41     |
| 04:08 | 16      | 13.4    | 28      | 1.4      |
| 05:07 | 16.3    | 14      | 26.2    | 4.1      |
| 06:06 | 17.32   | 14.2    | 26.4    | 5.12     |
| 07:05 | 18.6    | 14.5    | 26      | 7.1      |
| 08:04 | 19.2    | 15      | 25.3    | 8.9      |
| 09:03 | 20.2    | 15.3    | 25.1    | 10.4     |
| 10:02 | 23.3    | 16      | 25      | 14.3     |
| 11:01 | 27.2    | 16.32   | 28.1    | 15.42    |
| 12:00 | 28      | 17      | 26.2    | 18.8     |

## Turner Anderson Method

**Table 31: Turner Anderson Method 0.05M**


| RATIO | M:S(C1) | S:L(C2) | M:L(C3) | C1+C2-C3 |
|-------|---------|---------|---------|----------|
| 00:12 | 2.52    | 15.3    | 15.3    | 2.52     |
| 01:11 | 5.1     | 16.3    | 19.2    | 2.2      |
| 02:10 | 8.2     | 15.2    | 21.3    | 2.1      |
| 03:09 | 10.2    | 15.3    | 23.7    | 1.8      |
| 04:08 | 12.1    | 15.7    | 26.2    | 1.6      |
| 05:07 | 14      | 15      | 21.3    | 7.7      |
| 06:06 | 15.2    | 16.3    | 22.9    | 8.6      |
| 07:05 | 15.3    | 17.6    | 23.4    | 9.5      |
| 08:04 | 16.2    | 19.2    | 25.1    | 10.3     |
| 09:03 | 18.3    | 17      | 26.3    | 9        |
| 10:02 | 20.3    | 18.3    | 26.9    | 11.7     |
| 11:01 | 21.4    | 19      | 28      | 12.4     |
| 12:00 | 22.3    | 20.1    | 29.3    | 13.1     |

**Graph 21: Turner Anderson Method**

Spectroscopic estimation of Ferrous- 1, 10-Phenanthroline complex ( $\lambda_{\max}$ - 509.6 nm)

**Table 32: Spectroscopic data of Ferrous-1, 10-Phenanthroline complex**

| Conc. (ug/ml) | absorbance |
|---------------|------------|
| 0.025         | 0.0126     |
| 0.05          | 0.0114     |
| 0.1           | 0.0073     |
| 0.5           | 0.1982     |
| 1             | 0.4043     |
| 1.5           | 0.6143     |
| 2             | 0.7198     |
| 2.5           | 0.9175     |



**Graph 22: Standard Curve of Ferrous-1, 10-Phenanthroline complex**

From the standard curve it was seen that-  
Ferrous- EDTA complex having the concentration of 0.523129 ug/ml where absorbance of that solution was 0.582

Ferrous- PABA complex having the concentration of 0.337129 ug/ml where absorbance of that solution was 0.396.


Ferrous –Oxalic Acid complex having the concentration of 0.682129 ug/ml where absorbance of that solution was 0.741.

Ferrous Sulphate-EDTA + 1, 10 - Phenanthroline  
Monovariation Method

**Table 33: Monovariation Method**

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 0                      | 3.58        |
| 1                      | 3.51        |
| 2                      | 3.45        |
| 3                      | 3.41        |
| 4                      | 3.4         |
| 5                      | 3.38        |
| 6                      | 3.29        |
| 7                      | 3.21        |
| 8                      | 3.18        |
| 9                      | 3.11        |
| 10                     | 3.05        |
| 11                     | 2.99        |
| 12                     | 2.95        |
| 13                     | 2.93        |
| 14                     | 2.91        |
| 15                     | 2.87        |
| 16                     | 2.85        |
| 17                     | 2.8         |
| 18                     | 2.79        |
| 19                     | 2.72        |
| 20                     | 2.65        |
| 21                     | 2.62        |
| 22                     | 2.54        |
| 23                     | 2.5         |
| 24                     | 2.47        |
| 25                     | 2.45        |

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 26                     | 2.41        |
| 27                     | 2.38        |
| 28                     | 2.27        |
| 29                     | 2.24        |
| 30                     | 2.2         |
| 31                     | 2.18        |
| 32                     | 2.15        |
| 33                     | 2.12        |
| 34                     | 2.17        |
| 35                     | 2.19        |
| 36                     | 2.25        |
| 37                     | 2.27        |
| 38                     | 2.31        |
| 39                     | 2.4         |
| 40                     | 2.48        |
| 41                     | 2.5         |
| 42                     | 2.62        |
| 43                     | 2.65        |
| 44                     | 2.96        |
| 45                     | 3.02        |
| 46                     | 3.12        |
| 47                     | 3.19        |
| 48                     | 3.4         |
| 49                     | 3.42        |
| 50                     | 3.5         |



**Graph 23: Monovariation Method**


After addition of 33 ml of ligand solution to the metal solution where endpoint was seen 1 ml of 1, 10- phenanthroline was added and was estimated with spectroscopic technique.

Ferrous Sulphate-PABA +1, 10 –phenanthroline  
Monovariation Method

**Table 34: Monovariation Method**

| Volume of titrant (ml) | Conductance |
|------------------------|-------------|
| 0                      | 2.32        |
| 1                      | 2.39        |
| 2                      | 2.41        |
| 3                      | 2.43        |
| 4                      | 2.52        |
| 5                      | 2.55        |
| 6                      | 2.62        |
| 7                      | 2.68        |
| 8                      | 2.71        |
| 9                      | 2.74        |
| 10                     | 2.82        |
| 11                     | 2.87        |
| 12                     | 2.9         |
| 13                     | 2.91        |
| 14                     | 2.96        |
| 15                     | 3.02        |
| 16                     | 3.12        |
| 17                     | 3.23        |
| 18                     | 3.41        |
| 19                     | 3.52        |
| 20                     | 3.55        |
| 21                     | 3.61        |
| 22                     | 3.66        |
| 23                     | 3.54        |
| 24                     | 3.31        |
| 25                     | 3.08        |
| 26                     | 3.01        |
| 27                     | 2.98        |
| 28                     | 2.92        |
| 29                     | 2.87        |
| 30                     | 2.85        |
| 31                     | 2.81        |
| 32                     | 2.75        |
| 33                     | 2.72        |
| 34                     | 2.7         |
| 35                     | 2.68        |
| 36                     | 2.65        |
| 37                     | 2.62        |
| 38                     | 2.58        |
| 39                     | 2.52        |
| 40                     | 2.48        |
| 41                     | 2.45        |
| 42                     | 2.36        |
| 43                     | 2.3         |
| 44                     | 2.27        |
| 45                     | 2.25        |
| 46                     | 2.2         |
| 47                     | 2.18        |
| 48                     | 2.15        |
| 49                     | 2.1         |
| 50                     | 1.85        |





**Graph 24: Monovariation Method**


After addition of 22 ml of ligand solution to the metal solution where endpoint was seen 1 ml of 1, 10- phenanthroline was added and was estimated with spectroscopic technique.

Ferrous Sulphate- Oxalic Acid +1, 10- Phenanthroline Monovariation Method

**Table 35: Monovariation Method**

| Volume of Titrant added | Conductance |
|-------------------------|-------------|
| 0                       | 18.4        |
| 1                       | 18.1        |
| 2                       | 17.8        |
| 3                       | 17.6        |
| 4                       | 17.2        |
| 5                       | 17          |
| 6                       | 16.9        |
| 7                       | 16.8        |
| 8                       | 16.5        |
| 9                       | 16.4        |
| 10                      | 16.2        |
| 11                      | 16          |
| 12                      | 15.9        |
| 13                      | 15.7        |
| 14                      | 15.6        |
| 15                      | 15.3        |
| 16                      | 15.2        |
| 17                      | 15          |
| 18                      | 14.9        |
| 19                      | 14.8        |
| 20                      | 14.6        |
| 21                      | 14.5        |
| 22                      | 14.1        |
| 23                      | 14          |
| 24                      | 13.8        |
| 25                      | 13.6        |
| 26                      | 13.5        |
| 27                      | 13.2        |
| 28                      | 13.1        |
| 29                      | 12.9        |
| 30                      | 12.7        |
| 31                      | 12.5        |

| Volume of Titrant added | Conductance |
|-------------------------|-------------|
| 32                      | 12.3        |
| 33                      | 12.1        |
| 34                      | 12.6        |
| 35                      | 12.8        |
| 37                      | 13.1        |
| 38                      | 13.2        |
| 39                      | 13.4        |
| 40                      | 13.6        |
| 41                      | 14.1        |
| 42                      | 14.4        |
| 43                      | 14.5        |
| 44                      | 14.6        |
| 45                      | 15          |
| 46                      | 15.1        |
| 47                      | 15.2        |
| 48                      | 15.4        |
| 49                      | 15.6        |
| 50                      | 16          |



Graph 25: Monovariation Method

After addition of 33 ml of ligand solution to the metal solution where endpoint was seen 1 ml of 1, 10- phenanthroline was added and was estimated with spectroscopic technique.

**Precision**  
**System Precision**  
 Ferric Chloride-EDTA

Table 36: System Precision

| Volume of titrant (ml) | r1   | r2   | r3   | r4   | r5   | r6   | Mean        | SD       | %RSD     |
|------------------------|------|------|------|------|------|------|-------------|----------|----------|
| 1                      | 24.4 | 24.5 | 24.5 | 24.6 | 24.6 | 24.6 | 24.53333333 | 0.08165  | 0.332811 |
| 2                      | 24.1 | 24.2 | 24.2 | 24.3 | 24.3 | 24.3 | 24.23333333 | 0.08165  | 0.336931 |
| 3                      | 24.2 | 24.5 | 24.6 | 24.6 | 24.7 | 24.7 | 24.55       | 0.187083 | 0.762048 |
| 4                      | 24.7 | 24.8 | 24.7 | 24.7 | 24.7 | 24.7 | 24.71666667 | 0.040825 | 0.165171 |
| 5                      | 25   | 25   | 25   | 25.1 | 25   | 25   | 25.01666667 | 0.040825 | 0.163191 |
| 6                      | 25.4 | 25.4 | 25.4 | 25.3 | 25.4 | 25.4 | 25.38333333 | 0.040825 | 0.160833 |
| 7                      | 25.7 | 25.7 | 25.6 | 25.7 | 25.7 | 25.7 | 25.68333333 | 0.040825 | 0.158955 |
| 8                      | 25.5 | 25.9 | 26.3 | 26.3 | 26.3 | 26.3 | 26.1        | 0.334664 | 1.282238 |
| 9                      | 26.7 | 26.7 | 26.6 | 26.5 | 26.5 | 26.5 | 26.58333333 | 0.098319 | 0.369853 |
| 10                     | 27.1 | 27.1 | 27.1 | 27   | 27.1 | 27.1 | 27.08333333 | 0.040825 | 0.150738 |
| 11                     | 27.2 | 27.3 | 27.3 | 27.3 | 27.3 | 27.3 | 27.28333333 | 0.040825 | 0.149633 |

| Volume of titrant (ml) | r1   | r2   | r3   | r4   | r5   | r6   | Mean        | SD       | %RSD     |
|------------------------|------|------|------|------|------|------|-------------|----------|----------|
| 12                     | 27.9 | 28   | 28.1 | 28.1 | 28.2 | 28.2 | 28.08333333 | 0.116905 | 0.416277 |
| 13                     | 28.6 | 28.5 | 28.5 | 28.5 | 28.5 | 28.5 | 28.51666667 | 0.040825 | 0.143161 |
| 14                     | 28.9 | 29   | 29   | 29   | 29   | 29   | 28.98333333 | 0.040825 | 0.140856 |
| 15                     | 29.6 | 29.6 | 29.6 | 29.5 | 29.6 | 29.6 | 29.58333333 | 0.040825 | 0.137999 |
| 16                     | 30   | 30.1 | 30.1 | 30.1 | 30.1 | 30.1 | 30.08333333 | 0.040825 | 0.135706 |
| 17                     | 30.5 | 30.5 | 30.5 | 30.3 | 30.5 | 30.5 | 30.46666667 | 0.08165  | 0.167997 |
| 18                     | 30.9 | 31   | 31   | 31.1 | 31.1 | 31.1 | 31.03333333 | 0.08165  | 0.163103 |
| 19                     | 31.3 | 31.6 | 31.7 | 31.7 | 31.7 | 31.7 | 31.61666667 | 0.160208 | 0.506721 |
| 20                     | 31.9 | 31.9 | 31.9 | 31.8 | 31.9 | 31.9 | 31.88333333 | 0.040825 | 0.128044 |
| 21                     | 32.2 | 32.3 | 32.3 | 32.3 | 32.4 | 32.4 | 32.31666667 | 0.075277 | 0.132936 |
| 22                     | 32.8 | 32.8 | 32.9 | 32.9 | 32.9 | 32.9 | 32.86666667 | 0.05164  | 0.157119 |
| 23                     | 33.3 | 33.4 | 33.4 | 33.4 | 33.4 | 33.4 | 33.38333333 | 0.040825 | 0.122291 |
| 24                     | 33.4 | 33.5 | 33.6 | 33.6 | 33.6 | 33.6 | 33.55       | 0.083666 | 0.149377 |
| 25                     | 33.8 | 33.9 | 33.9 | 33.9 | 33.9 | 34   | 33.9        | 0.063246 | 0.186565 |
| 26                     | 35.1 | 35.1 | 35.1 | 35.2 | 35.2 | 35.2 | 35.15       | 0.054772 | 0.155824 |
| 27                     | 35.2 | 35.3 | 35.3 | 35.3 | 35.3 | 35.3 | 35.28333333 | 0.040825 | 0.115706 |
| 28                     | 35.5 | 35.5 | 35.5 | 35.6 | 35.5 | 35.5 | 35.51666667 | 0.040825 | 0.114946 |
| 29                     | 36.1 | 36.1 | 36.1 | 36.3 | 36.1 | 36.1 | 36.13333333 | 0.08165  | 0.225968 |
| 30                     | 36.5 | 36.5 | 36.6 | 36.6 | 36.6 | 36.6 | 36.56666667 | 0.05164  | 0.141221 |
| 31                     | 36.7 | 36.7 | 36.8 | 36.8 | 36.7 | 36.8 | 36.75       | 0.054772 | 0.14904  |
| 32                     | 37   | 37.2 | 37.2 | 37.2 | 37.2 | 37.2 | 37.16666667 | 0.08165  | 0.119685 |
| 33                     | 37.4 | 37.6 | 37.6 | 37.6 | 37.6 | 37.6 | 37.56666667 | 0.08165  | 0.117346 |
| 34                     | 37.7 | 37.7 | 37.7 | 37.7 | 37.7 | 37.8 | 37.71666667 | 0.040825 | 0.108241 |
| 35                     | 37.1 | 37.1 | 37.2 | 37   | 36.8 | 36.6 | 36.96666667 | 0.225093 | 0.608907 |
| 36                     | 35.7 | 35.6 | 35.6 | 35.6 | 35.7 | 35.7 | 35.65       | 0.054772 | 0.153639 |
| 37                     | 34.7 | 34.8 | 35   | 35   | 34.9 | 34.6 | 34.83333333 | 0.163299 | 0.468802 |
| 38                     | 32.6 | 32.7 | 32.9 | 33.1 | 32.9 | 32.9 | 32.85       | 0.176068 | 0.535976 |
| 39                     | 37.6 | 37.6 | 37.6 | 37.2 | 37.6 | 37.6 | 37.53333333 | 0.163299 | 0.435078 |
| 40                     | 37.7 | 37.7 | 37.6 | 37.6 | 37.6 | 37.7 | 37.65       | 0.054772 | 0.145477 |
| 41                     | 37.7 | 37.6 | 37.6 | 37.6 | 37.6 | 37.6 | 37.61666667 | 0.040825 | 0.108529 |
| 42                     | 37.6 | 37.7 | 37.7 | 37.7 | 37.7 | 37.7 | 37.68333333 | 0.040825 | 0.108337 |
| 43                     | 37.6 | 37.6 | 37.6 | 37.7 | 37.6 | 37.6 | 37.61666667 | 0.040825 | 0.108529 |
| 44                     | 37.7 | 37.6 | 37.6 | 37.6 | 37.6 | 37.6 | 37.61666667 | 0.040825 | 0.108529 |
| 45                     | 37.7 | 37.6 | 37.7 | 37.7 | 37.7 | 37.6 | 37.66666667 | 0.05164  | 0.137097 |
| 46                     | 37.7 | 37.7 | 37.7 | 37.6 | 37.7 | 37.2 | 37.6        | 0.2      | 0.531915 |
| 47                     | 37.7 | 37.6 | 37.6 | 37.7 | 37.7 | 37.7 | 37.66666667 | 0.05164  | 0.137097 |
| 48                     | 37.6 | 37.7 | 37.7 | 37.6 | 37.6 | 37.6 | 37.63333333 | 0.05164  | 0.137218 |
| 49                     | 37.7 | 37.6 | 37.6 | 37.7 | 37.7 | 37.6 | 37.65       | 0.054772 | 0.145477 |
| 50                     | 37.7 | 37.6 | 37.7 | 37.6 | 37.6 | 37.6 | 37.63333333 | 0.05164  | 0.137218 |

r= result of conductance.

All the %RSD value came within 2% so the above process is Precised Systematically.

### Method Precision

Ferric Ammonium Sulphate-Oxalic Acid

Table 37: Method Precision

| Volume of titrant (ml) | t1   | t2   | t3   | t4   | t5   | t6   | Mean        | SD       | %RSD     |
|------------------------|------|------|------|------|------|------|-------------|----------|----------|
| 1                      | 15.6 | 15.5 | 15.3 | 15.6 | 15.6 | 15.7 | 15.55       | 0.13784  | 0.886434 |
| 2                      | 15.4 | 15.7 | 15.7 | 15.8 | 15.7 | 15.9 | 15.7        | 0.167332 | 1.065809 |
| 3                      | 15.5 | 15.8 | 16.1 | 15.9 | 15.8 | 16   | 15.85       | 0.207364 | 1.308293 |
| 4                      | 15.9 | 16   | 16.4 | 16   | 15.9 | 16.2 | 16.06666667 | 0.196638 | 1.223891 |
| 5                      | 16.3 | 16   | 16.4 | 16.4 | 16.2 | 16.4 | 16.28333333 | 0.160208 | 0.983878 |
| 6                      | 16.4 | 16.5 | 16.7 | 16.4 | 16.3 | 16.6 | 16.48333333 | 0.147196 | 0.892999 |
| 7                      | 16.5 | 16.5 | 16.7 | 16.6 | 16.5 | 16.6 | 16.56666667 | 0.08165  | 0.492855 |
| 8                      | 16.8 | 16.5 | 17   | 16.9 | 16.7 | 16.8 | 16.78333333 | 0.17224  | 1.026257 |
| 9                      | 16.9 | 16.9 | 17.2 | 17.1 | 16.9 | 16.9 | 16.98333333 | 0.132916 | 0.782626 |
| 10                     | 17   | 17.2 | 17.4 | 17.2 | 17   | 17.2 | 17.16666667 | 0.150555 | 0.877017 |

| Volume of titrant (ml) | t1   | t2   | t3   | t4   | t5   | t6   | Mean        | SD       | %RSD     |
|------------------------|------|------|------|------|------|------|-------------|----------|----------|
| 11                     | 17.2 | 17.3 | 17.8 | 17.7 | 17.1 | 17.4 | 17.41666667 | 0.278687 | 1.600119 |
| 12                     | 17.6 | 17.6 | 18.1 | 17.7 | 17.2 | 17.6 | 17.63333333 | 0.287518 | 1.630538 |
| 13                     | 17.9 | 17.8 | 18.4 | 17.9 | 17.5 | 17.8 | 17.88333333 | 0.292689 | 1.636656 |
| 14                     | 18.2 | 18.2 | 18.7 | 18.4 | 18.2 | 18.1 | 18.3        | 0.219089 | 1.197208 |
| 15                     | 18.3 | 18.3 | 18.9 | 18.7 | 18.5 | 18.4 | 18.51666667 | 0.240139 | 1.296879 |
| 16                     | 18.6 | 18.7 | 19   | 18.9 | 18.6 | 18.9 | 18.78333333 | 0.17224  | 0.916984 |
| 17                     | 18.9 | 18.8 | 19.1 | 19.2 | 19   | 19.1 | 19.01666667 | 0.147196 | 0.774037 |
| 18                     | 19.2 | 19.2 | 19.2 | 19.4 | 19.1 | 19.3 | 19.23333333 | 0.10328  | 0.536982 |
| 19                     | 19.6 | 19.4 | 19.6 | 19.4 | 19.3 | 19.5 | 19.46666667 | 0.121106 | 0.62212  |
| 20                     | 19.5 | 19.8 | 19.8 | 20.1 | 19.5 | 19.8 | 19.75       | 0.225832 | 1.143452 |
| 21                     | 19.9 | 20   | 20.1 | 20.2 | 20.1 | 20   | 20.05       | 0.104881 | 0.523097 |
| 22                     | 20   | 19.6 | 20.2 | 20   | 20.2 | 20.4 | 20.06666667 | 0.273252 | 1.361721 |
| 23                     | 20.5 | 20.7 | 20.4 | 20.8 | 20.4 | 20.7 | 20.58333333 | 0.17224  | 0.836794 |
| 24                     | 20.6 | 20.9 | 20.8 | 21   | 20.7 | 20.9 | 20.81666667 | 0.147196 | 0.707107 |
| 25                     | 20.7 | 20.9 | 21.1 | 21.3 | 20.7 | 21   | 20.95       | 0.234521 | 1.119431 |
| 26                     | 21   | 21   | 21.4 | 21.6 | 21   | 21.4 | 21.23333333 | 0.265832 | 1.251956 |
| 27                     | 21.2 | 22   | 21.6 | 21.7 | 21.1 | 21.6 | 21.53333333 | 0.332666 | 1.544889 |
| 28                     | 21.7 | 22.1 | 21.8 | 21.7 | 21.3 | 21.8 | 21.73333333 | 0.258199 | 1.188032 |
| 29                     | 21.7 | 22.4 | 22.4 | 22.2 | 21.8 | 22.2 | 22.11666667 | 0.299444 | 1.353929 |
| 30                     | 22   | 22.7 | 22.5 | 22.5 | 22.1 | 22.3 | 22.35       | 0.266458 | 1.192207 |
| 31                     | 22   | 22.9 | 22.8 | 22.5 | 22.3 | 22.4 | 22.48333333 | 0.33116  | 1.472911 |
| 32                     | 22.4 | 23.6 | 23.2 | 23.2 | 22.9 | 22.6 | 22.98333333 | 0.440076 | 1.91476  |
| 33                     | 22.6 | 23.2 | 23.2 | 23.2 | 22.6 | 23.2 | 23          | 0.309839 | 1.347125 |
| 34                     | 22.8 | 23.9 | 23.7 | 23.6 | 23.1 | 23.6 | 23.45       | 0.413521 | 1.763418 |
| 35                     | 23.2 | 24.3 | 24.2 | 24.3 | 24.4 | 24.3 | 24.11666667 | 0.453505 | 1.880462 |
| 36                     | 24   | 24.5 | 24.4 | 24   | 23.6 | 24.3 | 24.13333333 | 0.332666 | 1.37845  |
| 37                     | 23.5 | 24.7 | 24.7 | 24.2 | 24   | 24.8 | 24.31666667 | 0.511534 | 2.103634 |
| 38                     | 23.8 | 24.9 | 25   | 24.6 | 24.2 | 24.9 | 24.56666667 | 0.476095 | 1.937972 |
| 39                     | 24   | 25.2 | 25.1 | 24.7 | 24.4 | 25.2 | 24.76666667 | 0.492612 | 1.989012 |
| 40                     | 25.1 | 25.5 | 25.4 | 24.7 | 24.7 | 25.4 | 25.13333333 | 0.361478 | 1.438243 |
| 41                     | 25   | 25.8 | 25.6 | 25.3 | 24.9 | 25.7 | 25.38333333 | 0.376386 | 1.482809 |
| 42                     | 25.2 | 26.1 | 25.2 | 25.3 | 25   | 25.9 | 25.45       | 0.441588 | 1.73512  |
| 43                     | 26   | 26.4 | 26.1 | 25.8 | 25.1 | 26.1 | 25.91666667 | 0.444597 | 1.715488 |
| 44                     | 26   | 26.9 | 26.5 | 25.9 | 26   | 26.3 | 26.26666667 | 0.382971 | 1.458011 |
| 45                     | 25.2 | 26.3 | 26.6 | 26.1 | 26   | 26.6 | 26.13333333 | 0.520256 | 1.990777 |
| 46                     | 27   | 27.5 | 26.9 | 26.5 | 27   | 26.6 | 26.91666667 | 0.354495 | 1.317009 |
| 47                     | 27   | 27.7 | 27.3 | 26.8 | 26.2 | 26.8 | 26.96666667 | 0.508593 | 1.886006 |
| 48                     | 28   | 28   | 27.4 | 27   | 27   | 26.9 | 27.38333333 | 0.507609 | 1.853714 |
| 49                     | 28   | 28.2 | 27.5 | 27.3 | 27.1 | 27.1 | 27.53333333 | 0.467618 | 1.698371 |
| 50                     | 28.4 | 28.4 | 27.9 | 27.4 | 27.3 | 27.3 | 27.78333333 | 0.526941 | 1.896608 |

t= number of titration with its conductance value.

This process is precised methodically as %RSD value comes within 2%.

### Difference between the observed value and expected value

Table 38: Difference between the observed value and expected value

| Method  | Complex            | Expected value                    | Observed Value (Approx.)           | Reason for deviation                                |
|---|--------------------|-----------------------------------|------------------------------------|---|
| Monovariation method<br>Modified Job's method<br>Turner Anderson method | Ferric-EDTA        | 1:1                               | .1:.151                            | Temperature/Ionic concentration of metal or ligand. |
|   |                    | 6:6                               | 7:5                                |   |
|   |                    | For 0.05M- 6:6<br>For 0.1M- : 6:6 | For 0.05M- 7:5<br>For 0.1M- : 7:5. |   |
| Monovariation method<br>Modified Job's method<br>Turner Anderson method | Ferric-Oxalic acid | 1:1                               | .1:.142                            | Temperature/Ionic concentration of metal or ligand. |
|   |                    | 6:6                               | 6:6                                |   |

| Method  | Complex                     | Expected value                                  | Observed Value (Approx.)                            | Reason for deviation                                |
|---|-----------------------------|---|---|---|
|   |                             | For 0.05M- 6:6<br>For 0.1M- : 6:6               | For 0.05M- 4:8<br>For 0.1M- : 5:7                   |   |
| Monovariation method<br>Modified Job's method<br>Turner Anderson method | Ferric- Curcumin            | 1:2<br>4:8<br>For 0.05M- 4:8<br>For 0.1M- : 4:8 | .1:.208<br>4:8<br>For 0.05M- 4:8<br>For 0.1M- : 4:8 | Temperature/Ionic concentration of metal or ligand. |
| Monovariation method<br>Modified Job's method<br>Turner Anderson method | Ferric-PABA                 | 1:2<br>4:8<br>For 0.05M- 4:8<br>For 0.1M- : 4:8 | .1:.227<br>4:8<br>For 0.05M- 4:8<br>For 0.1M- : 4:8 | Temperature/Ionic concentration of metal or ligand. |
| Monovariation method<br>Modified Job's method<br>Turner Anderson method | Ferrous-Oxalic acid         | 1:2<br>4:8<br>For 0.05M- 4:8<br>For 0.1M- : 4:8 | .1:.166<br>6:6<br>For 0.05M- 6:6<br>For 0.1M- : 6:6 | Temperature/Ionic concentration of metal or ligand. |
| Monovariation method<br>Modified Job's method<br>Turner Anderson method | Ferrous-EDTA                | 1:1<br>6:6<br>For 0.05M- 6:6<br>For 0.1M- : 6:6 | .1:.156<br>6:6<br>For 0.05M- 6:6<br>For 0.1M- : 4:8 | Temperature/Ionic concentration of metal or ligand. |
| Monovariation method<br>Modified Job's method<br>Turner Anderson method | Ferrous-1,10-Phenanthroline | 1:1<br>6:6<br>For 0.05M- 6:6<br>For 0.1M- : 6:6 | .1:.07<br>4:8<br>For 0.05M- 4:8<br>For 0.1M- : 4:8  | Temperature/Ionic concentration of metal or ligand  |

## CONCLUSION

Through this analysis, it has been observed that the formation of complex of Ferric and Ferrous salt with EDTA was taken place approximately at 1:1 ratio while using 0.1M metal and ligand solution for experiment (Monovariant method). Similarly in the case of PABA and Curcumin it was seen that the complexation takes place approximately at 1:2 ratio as metal and ligand ratio. The modified Job's method of continuous variation (Turner Anderson Method) was used to calculate the stability constant of the complex and the free energy change. The value of free energy change is negative showing the feasibility of complex formation.

1, 10- Phenanthroline was used as the chromogenic agent which essentially confirmed the presence of free metal content in the solution as EDTA, PABA and Oxalic Acid were used as ligand in the assay process of metal- chromogenic agent complex. The spectroscopic data shows the concentration from the

standard curve that where free metal was exist in the titration solution.

In future this method can also be used to find out the stoichiometric ratio of other metal and ligands along with the evaluation of different biological activities like in the treatment of Thalassemia when conjugation with ferrous ion takes place.

## ACKNOWLEDGEMENT

The authors are thankful to the respected Principal Sir, Dr. Arnab Samanta, Netaji Subhas Chandra Bose Institute of Pharmacy, West Bengal for providing necessary facilities for the completion of research work.

## REFERENCE

1. ICH Topic Q2B. (1996b). Validation of Analytical procedures: Methodology, 68.
2. Ingri, N., & Sillen, L. G. (1964). High-speed computers as a supplement to graphical methods. *Arkivor Kemi*, 23, 97-121.

3. Sayce, I. G. (1968). Computer calculations of equilibrium constants of species present in mixtures of metal ions and complexing reagents. *Talanta*, 15(12), 1397-1421.
4. Pearson, R. G. (2005). *Chemical Hardness: Applications from Molecules to Solids*. Manhattan, New York City: Springer-VCH; p. 210. ISBN: 978-3-527-60617-7.
5. Yoe, J. R., & Jones, A. L. (1944). *Ind. Eng. Chem., Anal. Ed.*, 16,111.
6. Sengupta, N. R. (1966). *Indian Journal of Chemistry*; 29, 33.
7. Job, P. (1928). Formation and stability of inorganic complexes in solution. *Ann. chim*, 9(10), 113-134.
8. Vosburgh, W. C., & Cooper, G. R. (1941). Complex ions. I. The identification of complex ions in solution by spectrophotometric measurements. *Journal of the American Chemical Society*, 63(2), 437-442.
9. Huhey, J. E. (1973). Proceedings of the Regional Annual Chemistry Teaching Symposium; *React.* 52-78.
10. Caradonna, J. P., & Lippard, S. J. (1988). Synthesis and characterization of [d (ApGpGpCpCpT)] 2 and its adduct with the anticancer drug cis-diamminedichloroplatinum (II). *Inorganic Chemistry*, 27(8), 1454-1466.
11. *Qualitative Inorganic Analysis*, Longmann Green and Co. London Vogel-I, II; (1959).
12. Anderson, W. F. (1980). *Inorganic Chemistry in Biology and Medicine*, ACS Symposium series, American Chemical Society, Washington DC.