

Open access • Posted Content • DOI:10.1101/2020.04.28.20083048

An epidemic CC1-MRSA-IV clone yields false negative test results in molecular MRSA identification assays: a note of caution — Source link [2]

Stefan Monecke, Stefan Monecke, Elisabeth König, Megan R. Earls ...+17 more authors


Institutions: Dresden University of Technology, Leibniz Institute of Photonic Technology, Medical University of Graz, University College Dublin ...+2 more institutions

Published on: 05 May 2020 - medRxiv (Cold Spring Harbor Laboratory Press)

Topics: SCCmec

Related papers:

- · Staphylococcal chromosome cassette evolution in Staphylococcus aureus inferred from ccr gene complex sequence typing analysis
- Extensive Genetic Diversity Identified among Sporadic Methicillin-Resistant Staphylococcus aureus Isolates Recovered in Irish Hospitals between 2000 and 2012
- · Genetic diagnosis of community-acquired MRSA: a multiplex real-time PCR method for Staphylococcal cassette chromosome mec typing and detecting toxin genes.
- Unexpected Diversity of Staphylococcal Cassette Chromosome mec Type IV in Methicillin-Resistant Staphylococcus aureus Strains
- Distinction of Staphylococcal Cassette Chromosome mec type V elements from Staphylococcus aureus ST398.











An epidemic CC1-MRSA-IV clone yields false negative test results in molecular MRSA identification assays: a note of caution

(gabriel.wagner-lichtenegger@medunigraz.at)

Stefan Monecke 1.2.3* (stefan.monecke@leibniz-ipht.de)
Elisabeth König 4* (elisabeth.ullrich@medunigraz.at)
Megan R. Earls 5* (megan.earls@dental.tcd.ie)
Eva Leitner 4** (eva.leitner@medunigraz.at)

Elke Müller ^{1,3} (elke.mueller@leibniz-ipht.de),

David M. Poitz ⁶ (David.Poitz@uniklinikum-dresden.de)

Lutz Jatzwauk ⁷ (Lutz.Jatzwauk@uniklinikum-dresden.de)

Teodora Vremeră ⁸ (doravre@yahoo.com)

Olivia S. Dorneanu ⁸ (odorneanu@yahoo.com)

Gabriel Wagner 4

Alexandra Simbeck ⁹ (alexandra.simbeck@klinik.uni-regensburg.de)

Andreas Ambrosch ¹⁰ (Andreas.ambrosch@barmherzige-regensburg.de)

Ines Zollner-Schwetz ¹¹ (ines.schwetz@medunigraz.at)

Robert Krause ¹¹ (robert.krause@medunigraz.at)

Werner Ruppitsch ¹² (werner.ruppitsch@ages.at)

Wulf Schneider-Brachert ⁹ (wulf.schneider@klinik.uni-regensburg.de)

David C. Coleman ⁵ (david.coleman@dental.tcd.ie)

Ivo Steinmetz ⁴ (ivo.steinmetz@medunigraz.at)

Ralf Ehricht ^{1,3,1,3} (ralf.ehricht@leibniz-ipht.de),

- 1) Leibniz Institute of Photonic Technology (IPHT), Jena, Germany
- 2) Institute for Medical Microbiology and Hygiene, Medical Faculty "Carl Gustav Carus", Technische Universität Dresden, Dresden, Germany
- 3) InfectoGnostics Research Campus Jena, Jena, Germany
- 4) Diagnostic and Research Institute for Hygiene, Microbiology and Environmental Medicine, Medical University of Graz, Graz, Austria
- 5) Microbiology Research Unit, Division of Oral Biosciences, Dublin Dental University Hospital, Trinity College, University of Dublin, Dublin, Ireland
- 6). Institute of Clinical Chemistry and Laboratory Medicine, University Hospital "Carl Gustav Carus", Technische Universität Dresden, Dresden, Germany
- 7) Department of Hospital Infection Control , University Hospital "Carl Gustav Carus", Technische Universität Dresden, Dresden, Germany
- 8) Microbiology Unit, Department of Preventive and Interdisciplinary Medicine, University of Medicine & Pharmacy "Grigore T Popa", Iași, Romania
- 9) Department of Infection Control and Infectious Diseases, University Hospital Regensburg, Regensburg, Germany
- 10) Institute of Laboratory Medicine, Microbiology and Hygiene, Barmherzige Brüder Hospital, Regensburg, Germany
- 11) Section of Infectious Diseases and Tropical Medicine, Department of Internal Medicine, Medical University of Graz, Austria
- 12) Institute of Medical Microbiology and Hygiene, Austrian Agency for Health and Food Safety, 1090 Vienna, Austria
- 13) Friedrich- Schiller University; Institute of Physical Chemistry, Jena, Germany

^{*}Contributed equally ** Corresponding author

ABSTRACT

Background: A variety of rapid molecular PCR tests has been developed and commercialised

that interrogate the junction site between the staphylococcal core genome, and the mobile

genetic element (SCCmec) which harbours the gene responsible for methicillin-/beta-lactam-

resistance.

Aim: The purpose of the present study was to investigate why a clinical MRSA isolate yielded

false negative test results with widely used, commercial orfX/SCCmec junction site PCR tests.

Methods: A collection of isolates that belonged to the same epidemic strain as the index

isolate were investigated with commercial MRSA assays and all isolates were sequenced in

order to explain this observation.

Results: It was found that isolates of the epidemic "European CC1-MRSA-IV" clone, which

likely originated in South-Eastern Europe and subsequently spread to Western Europe,

generally exhibit this behaviour. The failure of the assays was attributable to a characteristic

large insertion in the orfX/SCCmec integration site presumably targeted by such tests. In

contrast to MW2 (GenBank BA000033.2, a CC1 "USA400" strain which also harbours

SCCmec IVa), the European CC1 clone harbours an insertion of ca. 5,350 nucleotides

adjacent to orfX. This sequence starts with a novel SCC terminal sequence alternate to dcs and

encodes six different hypothetical proteins (E7MHX1, ydiL2, C5QAP8, A8YYX4, npd-SCC,

H4AYD7; nucleotide positions 280,690–286,024 of GenBank RBVO000005.1). An SCCmec

element with the same insertion was previously found in a Staphylococcus epidermidis isolate

(GenBank MH188467.1) suggesting transfer between staphylococcal species.

Conclusion: In order to ensure the reliability of molecular MRSA tests, it is vital to monitor

the emergence of new SCCmec junction sites, not only in Staphylococcus aureus but also in

coagulase-negative staphylococci.

INTRODUCTION

Staphylococcus aureus is a Gram-positive bacterium that frequently colonises both humans and animals, mostly as a harmless commensal of the skin and mucous membranes, especially the nares. However, *S. aureus* can also be pathogenic, causing skin and soft tissue infections, pneumonia, meningitis, endocarditis, mastitis, osteomyelitis, toxic shock syndrome and sepsis. It can harbour numerous virulence-associated genes, many of which reside on mobile genetic elements (MGEs). Infections caused by *S. aureus* can be treated with a variety of antibiotics including flucloxacillin, cefazolin, co-trimoxazole, clindamycin, vancomycin, doxycycline amongst others, depending on the clinical presentation, severity and the susceptibility profile of the infecting strain.

Methicillin-resistant *S. aureus* (MRSA) were first identified almost 60 years ago and now constitute a major obstacle to therapy. Methicillin resistance is mediated by a modified penicillin binding protein (PBP2′ or PBP2a) that confers resistance to virtually all beta-lactam antibiotics, with the exception of some recently developed compounds such as ceftobiprole and ceftaroline. The genes *mecA* or *mecC* typically encode this property. The *mecA* gene is most commonly responsible for methicillin resistance in *S. aureus* and is localised on a variety of large and complex mobile genetic elements (MGEs), called staphylococcal cassette chromosome elements (SCC*mec*). The *mecC* gene is localised on an uncommon SCC*mec* element (type XI; (1, 2)) and its presence is largely restricted to Western Europe where a distinct zoonotic connection to wildlife, especially hedgehogs, has been observed (3-5). Finally, there have been sporadic reports of plasmid-encoded genes, *mecB* (6) and *mecD* (7), which possibly originated in micrococci.

As therapy options for MRSA infections are limited, emphasis must be placed on infection prevention and control (IPC) measures including hand hygiene, the use of protective clothing and equipment and isolation of patients. Effective IPC implementation requires the rapid

detection of MRSA. Traditionally, this has been achieved using selective growth media containing methicillin or cefoxitin, salt and a chromogenic substrate (for an overview, see (8)). The latter facilitates the quick differentiation of MRSA from other, clinically less relevant staphylococcal species that may also harbour mec genes. While the use of chromogenic media is relatively inexpensive and convenient, plates must be incubated for 18-24 h before results are obtained. Considering that clinical judgements on whether to isolate a patient or to commence presumptive therapy with antibiotics that covers MRSA commonly require a much quicker decision, a variety of rapid molecular PCR tests has been developed and commercialised (9). These need to identify the presence of mecA and mecC and/or the presence of SCCmec elements together with a S. aureus-specific genetic marker(s). In addition, molecular tests should remain negative if the mec gene and/or the SCCmec element are present in another staphylococcal species that also might occur in the same sample. This can be achieved by placing one primer for the amplification reaction into the core chromosomal gene adjacent to the SCCmec element (orfX), which is present in all staphylococci, the sequence of which differs between different species. Importantly, the second primer cannot be placed into the mecA or mecC sequence as the distance between orfX and the mec gene is simply too large to be covered by a standard PCR. A "downstream constant segment" (dcs) of SCC is widely distributed and highly conserved among MRSA strains, even among strains harbouring different types of SCCmec elements, making it an obvious choice for the placement of a second amplification primer. However, this approach can have drawbacks. Firstly, methicillin-susceptible S. aureus might yield a positive test result if they carry an SCC element without a mec gene, such as a dysfunctional remnant of SCCmec (10) or an element associated with fusidic acid resistance. Secondly, an MRSA strain with a non-dcs downstream sequence might yield a negative result in such a test. In a central European setting, ca. 5% of MRSA strains might carry one of as many as fourteen known different SCC terminal sequences alternate to dcs (11). In order to ensure the reliability of

PCR-based MRSA detection assays, newly obtained sequences of the *orfX/SCCmec* junction site should be monitored, and the emergence of new SCC terminal sequences should prompt an update to these tests. Thirdly, there are strains that are phenotypically susceptible to beta-lactam antibiotics (12) despite being positive in molecular tests targeting *mecA* and/or the *orfX/SCCmec* junction site.

The purpose of the present study was to investigate why a clinical isolate yielded false negative test results with a widely used, commercial *orfX*/SCC*mec* junction site PCR tests (Cepheid GeneXpert MRSA/SA BC and BD MAX Staph SR Assay), although it was confirmed to be MRSA by other means. It was found that isolates of the epidemic "European CC1-MRSA-IV" clone, which likely originated in South-Eastern Europe and subsequently spread to Western Europe ((13, 14); for a detailed description, see also below), generally exhibit this behaviour. The failure of the assay was attributed to a characteristic large insertion in the *orfX*/SCC*mec* integration site presumably targeted by this test.

MATERIALS AND METHODS

Index case

The index case of the present study was a 62-year-old man with bacteraemia caused by CC1-MRSA which yielded a negative result using commercially available MRSA PCR systems in widespread use.

The index patient, who had metastasised lung cancer, was admitted to the Department of Internal Medicine, Medical University of Graz (Austria) with suspected pneumonia, and amoxicillin/clavulanic acid therapy was started. Blood cultures, incubated in an automated BACTEC system (Becton Dickinson; Heidelberg, Germany/Franklin Lakes, NJ, USA) became positive within one day. The organism was identified as S. aureus directly from positive blood culture bottles by peptide nucleic acid fluorescence in situ hybridisation (PNA FISH, AdvanDx, Woburn, MA, USA). The S. aureus isolate from one blood culture bottle (Graz_511421-19) was investigated by using the GeneXpert MRSA/SA BC PCR assay (Cepheid, Sunnyvale, CA, USA) in order to identify S. aureus/MRSA. Simultaneously, rapid antimicrobial susceptibility testing (RAST) was performed according to EUCAST guidelines (15, 16). The GeneXpert test was negative for MRSA but RAST revealed cefoxitin resistance after six hours of incubation. Additional antimicrobial susceptibility testing (Vitek2, bioMérieux, Marcy-l'Étoile, France) was performed confirming that it was indeed an MRSA isolate. Another molecular MRSA detection assay, BD MAXTM StaphSR assay (Becton Dickinson) was performed; it also yielded a negative result. Microarray-based characterisation with the S. aureus Genotyping Kit 2.0 (Abbott [Alere Technologies GmbH], Jena, Germany; see below) detected mecA, confirming the disk diffusion result for MRSA, and assigned the isolate to CC1-MRSA-IV.

Initially, before there was evidence for an MRSA infection, the patient's condition improved under beta-lactam treatment, and therapy was changed to flucloxacillin and fosfomycin. One day later, after antimicrobial susceptibility testing confirmed MRSA with susceptibility to fosfomycin and daptomycin, the treatment was modified and flucloxacillin replaced by daptomycin in addition to fosfomycin. Despite the fact that daptomycin has no pulmonary activity, alternative substances like vancomycin or linezolid were contraindicated due to nephro- and myelotoxicity. Fifth generation cephalosporins were not available. The patient died a few days later because of tumour progression.

The index isolate was identified as European CC1-MRSA-IV using DNA microarray profiling and whole-genome sequencing (see below).

Additional CC1-MRSA-IV isolates investigated

Beside the index case, another 38 additional isolates of the European CC1-MRSA-IV clone were investigated.

Ten Romanian CC1-MRSA-IV isolates were recovered between 2007 and 2010, during routine microbiological diagnostics in the Sfanta Parascheva Hospital for Infectious Diseases in Iasi, a 300-bed tertiary centre in North-Eastern Romania, where European CC1-MRSA-IV appears to be common (17).

Four isolates were recovered from diagnostic or screening procedures at the Medical Faculty of Dresden, a 1,200-bed tertiary care hospital in South-Eastern Germany. In Dresden, European CC1-MRSA-IV is very rare, with only six isolates being identified at the University Hospital among 1,752 isolates characterised by microarray profiling between 2000 and 2020 ((11) and unpublished data). In addition to European CC1-MRSA-IV isolates, four control strains were included in this study (see Table 1). Controls were assigned by microarray profiling (see below) to epidemic strains that are or were prevalent locally, or that belonged to

another CC1-MRSA strain (isolate #220663). In addition, fully sequenced reference strains MU50, MW2, N315 and USA300-FPR3757 were also included as controls.

Nineteen Bavarian isolates recovered at the University Medical Centre in Regensburg between 2010 and 2019 were included, along with five isolates from other tertiary care centres in Bavaria. In Bavaria, the European CC1-MRSA-IV strain appears to be common. At the University Medical Centre in Regensburg, retrospective microarray profiling-based typing of 3,067 isolates revealed a clear increase in the prevalence of this strain. It accounted for <1% of typed MRSA between 2010 and 2013, but reached a prevalence of 9.4% (26 out of 277) in 2019, being the third most common strain (behind CC5-MRSA-II and CC398-MRSA-V/VT) at that centre.

All study isolates were preserved at -80°C in cryovials following identification as MRSA, and molecular typing was performed retrospectively by microarray profiling.

Table 1:MRSA strains and isolates investigated in the present study and their detection using commercially available *orfX*/SCC*mec* junction site assays.

Isolate	Strain affiliation (according to microarray)	SCCmec element	Origin	BD MAX results	GeneXpert results
Graz_511421-19 (index case)	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Austria, 2019	NEGATIVE (G)	NEGATIVE (2 x BC; G)
Dresden-94757	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Saxony, 2010	NEGATIVE (D, G)	N/A
Dresden-94758	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Saxony, 2014	NEGATIVE (D, G)	N/A
Dresden-94759	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Saxony, 2009	NEGATIVE (D, G)	N/A
Dresden-94760	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Saxony, 2010	NEGATIVE (D, G)	N/A
Iasi-95033	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Romania, 2009	NEGATIVE (D, G)	N/A
Iasi-95034	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Romania, 2009	NEG. (D)/AMBIGUOUS (G)*	N/A
Iasi-95035	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Romania, 2009	NEG. (D)/AMBIGUOUS (G)*	N/A
Iasi-95037	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Romania, 2009	NEGATIVE (D, G)	N/A
Iasi-95038	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Romania, 2009	NEGATIVE(D, G)	N/A
Iasi-95039	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Romania, 2009	NEG. (D)/AMBIGUOUS (G)*	N/A
Iasi-95040	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Romania, 2009	NEGATIVE (D, G)	N/A
Iasi-95041	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Romania, 2009	NEGATIVE (D, G)	N/A
Iasi-174752	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Romania, 2010	NEGATIVE (D, G)	N/A
Iasi-176047	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Romania, 2009	NEGATIVE (D, G)	N/A
Bavaria-0643	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2018	N/A	POSITIVE (SSTI; R)
Bavaria-0824	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2015	N/A	POSITIVE (SSTI; R)
Bavaria-1185	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2018	N/A	POSITIVE (SSTI; R)
Bavaria-1274	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2014	N/A	POSITIVE (SSTI; R)
Bavaria-1537	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2013	N/A	POSITIVE (SSTI; R)
Bavaria-1780	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2013	N/A	POSITIVE (SSTI; R)
Bavaria-1962	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria	N/A	POSITIVE (SSTI; R)
Bavaria-2102	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2019	N/A	POSITIVE (SSTI; R)
Bavaria-2220	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria	N/A	POSITIVE (SSTI; R)
	1				

Isolate	Strain affiliation (according to microarray)	SCCmec element	Origin	BD MAX results	GeneXpert results
Bavaria-2312	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria	N/A	POSITIVE (SSTI; R)
Bavaria-2360	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria	N/A	POSITIVE (SSTI; R)
Bavaria-2391	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2018	N/A	POSITIVE (SSTI; R)
Bavaria-2483	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2019	N/A	POSITIVE (SSTI; R)
Bavaria-2535	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2019	N/A	POSITIVE (SSTI; R)
Bavaria-2584	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2019	N/A	POSITIVE (SSTI; R)
Bavaria-2585	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2019	N/A	POSITIVE (SSTI; R)
Bavaria-2588	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2019	N/A	POSITIVE (SSTI; R)
Bavaria-2596	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2019	N/A	POSITIVE (SSTI; R)
Bavaria-2618	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2012	N/A	POSITIVE (SSTI; R)
Bavaria-3012	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2011	N/A	POSITIVE (SSTI; R)
Bavaria-3254	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2010	N/A	POSITIVE (SSTI; R)
Bavaria-3702	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2019	N/A	POSITIVE (SSTI; R)
Bavaria-3741	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2019	N/A	POSITIVE (SSTI; R)
Bavaria-3784	CC1-MRSA-IV (PVL-neg., aphA3/sat-positive)	SCCmec IVa with insertion	Bavaria, 2019	N/A	POSITIVE (SSTI; R)
Dresden-220663	CC1-MRSA-IV, (PVL-neg., aphA3/sat-negative)	SCCmec IVa (as in MW2)	Saxony, 2007	POSITIVE (D)	N/A
Dresden-124288	CC22-MRSA-IV (Barnim/UK EMRSA-15)	SCCmec IVh/j	Saxony	POSITIVE (D)	N/A
Dresden-124289	CC22-MRSA-IV (Barnim/UK EMRSA-15)	SCCmec IVh/j	Saxony	POSITIVE (D)	N/A
Dresden-124281	CC45-MRSA-IV (Berlin EMRSA)	SCCmec IVa	Saxony	POSITIVE (D)	N/A
MU50	CC5-MRSA-II (New York/Japan Clone)	SCCmec II	Japan (sequenced reference strain)	POSITIVE (D)	N/A
MW2	CC1-MRSA-IV (PVL-pos. USA400)	SCCmec IVa	USA (sequenced reference strain)	POSITIVE (D)	N/A
N315	CC5-MRSA-II (New York/Japan Clone)	SCCmec II	Japan (sequenced reference strain)	POSITIVE (D)	N/A
USA300-FPR3757	CC8-MRSA-[IVa+ACME1] (PVL+), USA300	SCC [mec IVa+ACME1+Cu]	USA (GenBank CP000255.1)	POSITIVE (D)	N/A

D, G, and R indicate that assays were performed in Dresden, Graz and Regensburg, respectively. BC, GeneXpert MRSA/SA BC. SSTI, GeneXpert MRSA/SA SSTI. N/A, not available/not performed.

^{*}AMBIGUOUS, weak signal observed at Ct>35

Cepheid GeneXpert

The index specimen from a positive blood culture was tested with the with the Cepheid

GeneXpert system using the MRSA/SA BC assay (Lot. No. 1000148707). It was later re-

tested with another lot of the same assay (Lot no. 1000179462). Further investigations on

GeneXpert MRSA/SA BC were not possible due to restricted laboratory capacities during the

influenza season and the SARS-CoV-2 pandemic.

Twenty-four Bavarian isolates were tested in Regensburg with the Cepheid GeneXpert

MRSA/SA SSTI assay (Lot no. 1000180532). It was used according to manufacturer's

instructions starting from pure overnight cultures on CBA.

BD MAX

The index isolate was tested in Graz using the BD MAX Staph SR assay (Ref.: 443419; Lot

no. 9303156). Additional isolates as well as control strains (see above and Table 1) were

grown overnight on Columbia Blood Agar (CBA) and sent for processing to the laboratories

in Graz and Dresden where the BD MAX Staph SR Assay (Graz, Lot no. 9303156; Dresden,

Lot no. K55928980720210312) was used according to the manufacturer's instructions.

Microarray-based molecular characterisation

Genotyping of all strains was performed using the S. aureus Genotyping Kit 2.0 (Abbott

[Alere Technologies GmbH]). The microarray-based assay covers 333 different targets related

to approximately 170 different genes and their allelic variants. The list of target genes as well

as sequences of probes and primers have been described in detail previously along with all

relevant protocols (18).

Briefly, *S. aureus* was cultivated and harvested on CBA. Subsequently, DNA extraction was performed with an enzymatic combination of lysis enzymes and buffer from the *S. aureus* Genotyping Kit 2.0 kit and Qiagen DNA extraction columns (Qiagen, Hilden, Germany), according to the manufacturers' instructions. Linear amplification was performed using one specific primer per target. During amplification, biotin-16-dUTP was incorporated into amplicons. After incubation and several washing steps, hybridisation to the probes on the array was detected using streptavidin horseradish peroxidase that caused local precipitation at the spot where the specific amplicon was bound to the probe. Microarrays were photographed and analysed with a designated reader and software (IconoClust, Abbott (Alere Technologies GmbH)). Array profiling permits rapid isolate characterisation, detects the presence or absence of certain genes or alleles, automatically compares query isolates to a database of existing patterns, and assigns isolates to clonal complexes, known strains and SCC*mec* types.

Whole-genome sequencing

All isolates underwent whole-genome sequencing (WGS). Genomic DNA was extracted as detailed above for microarray-based molecular characterisation, and its quality was assessed as previously described (19). The Nextera DNA Flex Library Preparation Kit (Illumina, Eindhoven, The Netherlands) was used according to the manufacturer's instructions and libraries underwent paired-end sequencing using the 500-cycle MiSeq Reagent Kit v2 (Illumina). Libraries were scaled to exhibit at least 50-fold coverage and the quality of each sequencing run was assured following cluster density and Q30 assessment. Raw sequence reads were trimmed using fastp version 0.19.11 (20). Trimmed sequence read sets underwent *de novo* assembly using SPAdes version v3.9.1 (21) and all contigs under 1,000 bp were removed.

RESULTS

Cepheid GeneXpert tests and BD MAX

Test results are provided in detail in Table 1.

Testing of the index MRSA isolate with two different lots of the GeneXpert MRSA/SA BC assay yielded negative results. However, the GeneXpert MRSA/SA SSTI assay correctly identified 24 European CC1-MRSA-IV isolates as MRSA.

The BD MAX assay system failed to identify isolates (n = 15) of the European CC1-MRSA-IV clone although control reference MRSA strains handled in parallel were correctly identified.

Description of the European CC1 strain

Microarray profiling and WGS analysis showed that the index MRSA isolate (Graz_511421-19) belonged to a PVL-negative CC1-MRSA-IV clone previously described as "European CC1-MRSA-IV" (13, 14, 17). This clone has been identified in Romania, Ireland and Germany (13, 14, 17) and is represented by a published genome sequence (RBVO000000001; (14)).

European CC1-MRSA-IV isolates typically exhibit sequence type (ST) 1 (1-1-1-1-1-1) or its single locus variant ST4110 (1-1-1-1-1-558), and *spa* type t127 (07-23-21-16-34-33-13) or to related types such as t386 (07-23-13) or t13790 (07-23-21-16-34-33-34-34-33-34).

This clone is resistant to beta-lactam antimicrobials due to a functional *mecA* encoded on an SCC*mec* IV element (see below). It is usually resistant to erythromycin and clindamycin (*ermC*), tetracycline (*tetK*), as well as to kanamycin and neomycin (*aphA3*, accompanied by

aadE and sat). Some isolates of this clone also carry aacA-aphD and are resistant to gentamicin. This was not the case for the index isolate from Graz. A particular variant of this clone from Ireland is resistant to mupirocin (mupR) and quaternary ammonium compounds such as chlorhexidine due to the presence of a plasmid encoding iles2/mupR and qacA, as previously described (14). Fusidic acid resistance has not yet been detected in this epidemic clone. In CC1, this property is usually encoded by SCCmec-borne fusC, which indicates affiliation to other CC1 strains (of Middle Eastern or Southern Hemisphere background); and fusB/far-1 was not detected in isolates tested hitherto.

European CC1-MRSA-IV is PVL-negative and only very rarely carries supplementary enterotoxin genes (*sek/seq*) in addition to *seh* that is ubiquitously present in CC1.

Description of the SCCmec element of the European CC1 strain

The SCC*mec* element sequences of the European CC1-MRSA-IV clone were essentially identical in all study isolates as well as in the previously sequenced Irish isolate (GenBank RBVO). Minor polymorphisms were noted within *dru* and within the region between Q2FKL7 and Q8VUV8. These were not relevant for the present study, and can be considered random duplications/deletions, or possibly even as assembly artefacts.


This CC1-MRSA strain has an SCC*mec* IVa element. In contrast to MW2-like SCC*mec* IVa however, it harbours an insertion of ca. 5,350 nucleotides, adjacent to *orfX*. This sequence starts with a novel SCC terminal sequence alternate to *dcs* ("SCCterm 15"; see below) and encodes six different hypothetical proteins (E7MHX1, *ydiL*2, C5QAP8, A8YYX4, *npd*, H4AYD7; nucleotide positions 280,690–286,024 of GenBank RBVO000005.1). Compared to MW2 (BA000033.2), this insertion replaces *dcs* and Q9XB68-*dcs*, and removes the larger part (212 of 240 nt) of the gene encoding hypothetical protein, Q7A213.

The entire gene content of the SCC*mec* IV element is listed in Table 2, and Figure 1 provides a schematic representation of the MGE, including a comparison to MW2 (BA000033.2) and to the sequence of an *S. epidermidis* composite SCC*mec* element (see below).

Table 2: Genes in the variant SCC*mec* IVa element of the European CC1-MRSA-IV strain.

Gene ID	Gene ID/definition of gene product and comments (see also annotation of MH188467.1)	Orientation	Locus tag in MW2 (BA000033.2)	Coordinates in GenBank RBVO	Coordinates in Iasi- 95037 (MT380478)
orfX	23S rRNA methyltransferase		MW0024	RBVO01000005.1; pos. 280209 to 280689	
sRNA6	Antisense RNA associated with <i>orfX</i>		N/A	RBVO01000005.1; pos. 280389 to 280673	
DR_SCC	Direct repeat of SCC, to 19 nt of the 3' end of the coding sequence of orfX		N/A	RBVO01000005.1; pos. 280670 to 280689	pos. 1 to 19
sccterm15	SCC-terminal sequence adjacent to orfX, and alternate to dcs, see Discussion		Not present	RBVO01000005.1; pos. 280689 to 280912	pos. 20 to 242
E7MHX1	Transcription regulator	FORWARD	Not present	RBVO01000005.1; pos. 280912 to 281239	pos. 243 to 569
ydiL2	Hypothetical protein/putative membrane peptidase, associated with SCC elements	FORWARD	Not present	RBVO01000005.1; pos. 281275 to 282109	pos. 606 to 1439
C5QAP8-M299	Hypothetical protein	FORWARD	Not present	RBVO01000005.1; pos. 282794 to 283568	pos. 2125 to 2898
A8YYX4	Hypothetical protein	REVERSE	Not present	RBVO01000005.1; pos. 283805 to 284144	pos. 3136 to 3474
npd-SCC	Enoyl-[acyl-carrier-protein] reductase-like protein	REVERSE	Not present	RBVO01000005.1; pos. 284329 to 285400	pos. 3660 to 4730
H4AYD7-trunc	Transcriptional regulator, LysR family	TRUNCATED	Not present	RBVO01000005.1; pos. 285412 to 286024	pos. 4743 to 5354
Q7A213-trunc	Putative protein; it comprises the inverted repeat of IS431. In MW2 it is not truncated and comprises 240 nt.	TRUNCATED	MW0026	RBVO01000005.1; pos. 286024 to 286052	pos. 5355 to 5382
IR_IS431	Inverted repeat of IS431	TRUNCATED	N/A	RBVO01000005.1; pos. 286024 to 286040	pos. 5355 to 5370
tnpIS431	Transposase for IS431	REVERSE	MW0027	RBVO01000005.1; pos. 286083 to end of contig (pos. 286184) (partial)	pos. 5414 to 6088
Teg143	Trans-encoded RNA associated with tnpIS431(22)		N/A	RBVO01000003.1; pos. 203 to 237	pos. 6119 to 6152
IR_IS431	Inverted repeat of IS431	TRUNCATED	N/A	RBVO01000003.1; pos. 213 to 229	pos. 6129 to 6144
mvaS-SCC	Truncated HMG-CoA synthase	FORWARD	MW0028	RBVO01000003.1; pos. 245 to 598	pos. 6161 to 6513
Q5HJW6	Hypothetical protein	FORWARD	N/A	RBVO01000003.1; pos. 695 to 1046	pos. 6611 to 6841
dru	SCC direct repeat units	TRUNCATED	N/A	RBVO01000003.1; pos. 835 to 1273	pos. 6751 to 7148
ugpQ	Glycerophosphoryl diester phosphodiesterase-like protein	FORWARD	MW0029	RBVO01000003.1; pos. 1474 to 2218	pos. 7350 to 8093
ydeM	Acyl dehydratase MaoC	FORWARD	MW0030	RBVO01000003.1; pos. 2314 to 2743	pos. 8190 to 8618
mecA	Encodes penicillin binding protein 2 prime, defining MRSA	REVERSE	MW0031	RBVO01000003.1; pos. 2812 to 4795	pos. 8664 to 10670
mecR1-trunc	Methicillin resistance operon repressor 1, signal transducer protein, truncated in SCCmec IV	TRUNCATED	MW0032	RBVO01000003.1; pos. 4894 to 5862	pos. 10770 to 10816
hsdR2-IS1272	Type I site-specific deoxyribonuclease restriction subunit	TRUNCATED	MW0033	RBVO01000003.1; pos. 5869 to 6103	pos. 11745 to 11978
tnpIS1272	Transposase	REVERSE	MW0034	RBVO01000003.1; pos. 6103 to 7627	pos. 11979 to 13502
Q9KX75	Hypothetical protein	REVERSE	MW0035	RBVO01000003.1; pos. 7762 to 8269	pos. 13638 to 14144
Q7A207	Hypothetical protein	REVERSE	MW0036	RBVO01000003.1; pos. 8283 to 8595	pos. 14159 to 14470
Q7A206-trunc	Hypothetical protein, truncated	TRUNCATED	N/A	RBVO01000003.1; pos. 8596 to 8683	pos. 14472 to 14558
Q7A206	Hypothetical protein	REVERSE	MW0037	RBVO01000003.1; pos. 8681 to 9032	pos. 14557 to 14907
UTR_ccrB-2	Highly conserved 3'-untranslated region of <i>ccrB</i>		N/A	RBVO01000003.1; pos. 9032 to 9553	pos. 14908 to 15428
ccrB-2	Cassette chromosome recombinase B2	REVERSE	MW0038	RBVO01000003.1; pos. 9553 to 11182	pos. 15429 to 17057
ccrA-2	Cassette chromosome recombinase A2	REVERSE	MW0039	RBVO01000003.1; pos. 11203 to 12553	pos. 17079 to 18428
cch-2	Hypothetical protein/cassette chromosome helicase	REVERSE	MW0040	RBVO01000003.1; pos. 12786 to 14574	pos. 18662 to 20449
DUF1413	Hypothetical protein, associated with cch	REVERSE	MW0041	RBVO01000003.1; pos. 14573 to 14864	pos. 20449 to 20739
Q2FKL7	Putative membrane protein	FORWARD	MW0042	RBVO01000003.1; pos. 15002 to 16052	pos. 20878 to 21927
Q8VUV8	Putative transcriptional regulator	FORWARD	MW0043	RBVO01000003.1; pos. 16504 to 17995	pos. 22380 to 23870
cstB-scc2	Includes a putative beta-lactamase; marker for SCCmec IVa	TRUNCATED	MW0045	RBVO01000003.1; pos. 18367 to 19685	pos. 24243 to 25560
Q2FKL3	HNH endonuclease family protein	FORWARD	MW0046	RBVO01000003.1; pos. 19875 to 20247	pos. 25751 to 26122
Q8VUW0	Putative membrane protein	FORWARD	MW0047	RBVO01000003.1; pos. 20375 to 20996	pos. 26251 to 26871
DR_SCC	Direct repeat of SCC		N/A	RBVO01000003.1; pos. 21300 to 21319	pos. 27176 to 27194

Figure 1: Graphic representation of the SCC*mec* elements in the CC1 reference sequence MW2 (GenBank BA000033.23), in the European CC1-MRSA-IV isolate Iasi-95037 (GenBank submission pending) and in the *S. epidermidis* isolate P11PPP12 (GenBank MH188467.1).



DISCUSSION

The study demonstrates that a CC1-MRSA-IV strain prevalent in Europe can yield falsenegative results with frequently used commercial MRSA assays (GenXpert MRSA/SA BC
and BD MAX Staph SR). Interestingly, CC1-MRSA-IV strains yielded correctly positive
results when using the GeneXpert MRSA/SA SSTI, indicating that the tests are not equipped
with the same primers and/or that there is already an updated version for one MRSA test
available from Cepheid. Unfortunately, this is not the case for the blood culture (MRSA/SA
BC) version that is primarily used for diagnosis in invasive infections. The primers used in
these commercial tests are not disclosed, however, the absence of *dcs* and the presence of a
large insert in the region downstream of *orfX* correlate with negative GenXpert MRSA/SA
BC test and BDmax results. It can therefore be assumed that the diverging results of the
assays are related to their coverage of a sequence fragment alternate to *dcs* (SCCterm 15; see
Table 2 and Figure 1).

This sequence (SCCterm 15) has previously been observed in at least three different *S. aureus* strains (for discussion of *S. epidermidis*, see below). In CC152-MRSA-XIII, SCCterm 15 is located adjacent to *orfX* (GenBank MG674089, CP024998; (23)) and thus, may affect MSRA assays spanning the *orfX*/SCC*mec* junction. For the remaining two strains, this sequence is situated away from *orfX*, within complex composite SCC*mec* elements. One of these strains is a Danish CC8-MRSA strain (24) with an SCC*mec* IVa/ACME composite element (*e.g.*, Strain M1, GenBank HM030720.1). In this strain, the same insert as in the "European CC1" strain is localised at the boundary between an ACME-II element and an SCC*mec* element (25). Here, SCCterm 15 is not adjacent to *orfX*, and the negative junction site PCR results observed can be attributed to yet another SCCterm sequence. Another, CC22-MRSA, strain (CMFT3119, GenBank HF569105.1) from Saudi Arabia also harboured SCCterm 15, E7MHX1, *ydiL*2,

IR_IS431 and tnpIS431 localised between a copper resistance element and a composite element encompassing ACME-II and SCC*mec* IVh/j (26).

SCCterm 15, E7MHX1, ydiL2, C5QAP8, A8YYX4, npd and H4AYD7 are also present in S. epidermidis P11PPP12 (GenBank MH188467.1), suggesting a transfer of this region between coagulase-negative staphylococci and CC1 S. aureus. Beyond these genes, the entire SCCmec IV element of P11PPP12 is identical to that in the CC1-MRSA strain. Most significantly, the fault line cutting short Q7A213 is conserved in both strains (position 5818/5819 in Supplemental file S1), suggesting it unlikely that the insert alone was transferred from a P11PP12-like ancestor into S. aureus CC1. Indeed, it is more likely that the entire variant SCCmec IVa cassette (including the insert) was transferred between the ancestors of the two strains, although the direction of this transfer remains unknown. P11PPP12 also harbours an ACME-II element downstream of SCCmec IVa, which comprises another SCC terminal sequence (as in FR753166.1, pos. 481-568), arc genes, another copy of tnpIS431, and heavy metal resistance genes. ACME-II is absent from the European CC1-MRSA-IV clone. Therefore, either it must have been lost during or after transfer of the SCC element between S. epidermidis and CC1, or it was later acquired by ancestors of S. epidermidis P11PPP12. Most likely, the MRSA emerged following transfer of an "atypical" (dcs-negative) SCCmec IVa element from coagulase-negative staphylococci into a CC1-MSSA strain that is prevalent in South-Eastern Europe.

The presence of MRSA that are negative in *orfX/SCCmec* junction site PCRs provides reasons for concern. Firstly, molecular assays are used to predict a presence of MRSA in positive blood cultures and, hence, to prompt administration antibiotics that cover MRSA. A false negative result would harm the patient by delaying effective therapy until conventional susceptibility test results become available. This could easily have fatal consequences. The use of molecular assays for the guidance of therapy in life-threatening circumstances is

beneficial for a vast majority of patients due to speed. However, these assays can only detect target sequences that were considered upon primer design. Conventional susceptibility tests are slower but they are not constrained by the choice of primers or by the presence of unknown or unexpected genotypes. Thus, they must not be neglected and should be carried out in parallel to molecular assays.

Secondly, molecular assays are used to guide IPC. A false negative result may result in failure to isolate an infected patient and further spread of the strain. Regarding IPC, it is also noteworthy that this strain can harbour plasmid-borne genes *qacA* and *ileS2/mupR* (14), encoding resistance to quaternary ammonium compounds, biguanides and mupirocin, with the latter compounds commonly being recommended for MRSA decolonisation (27). Thus, it may not only escape detection by molecular screening, but also circumvent decolonisation.

A third and less obvious consequence is a potential shift in the clonal structure of local MRSA population with molecular assays themselves providing a selective pressure favouring the false-negative strain. Other strains will get "penalised" for being PCR-positive as subsequent interventions would hinder their proliferation and transmission. In the long term, this would lead to a higher prevalence of the false-negative strain, and increasingly more therapy and IPC failures.

As previously discussed, the European CC1-MRSA-IV clone may have emerged in South-Eastern Europe. A putative methicillin-susceptible ancestor is common in Romania where the MRSA strain was early and frequently observed (14, 17). Significant case numbers have been reported from Bavaria (14), North Rhine-Westphalia (28), Italy (29) and Ireland (14). Sporadic cases have been described in Austria (30) and Saxony (11).

A strategy to contain the CC1 strain must rely on conventional susceptibility tests or updated molecular tests. We therefore propose to screen, at least, patients with conspicuous travel histories as well as medical or nursing staff recruited in epidemic regions, not only by

molecular means but also by culture of nasal swabs on selective growth media. Sequence information and isolates described herein can be used by test manufacturers to adapt future versions of their assays in order to identify this strain.

The problem of unknown or unexpected genotypes also represents a considerable economic risk for the manufacturer of molecular point-of-care tests. The staging and approval of CE—marked in vitro diagnostic tests is very expensive and time-consuming, while there is an absolute necessity to adapt tests to constant evolution of targeted pathogens. It is therefore important for the manufacturers to be able to integrate new markers or new sequence variants into existing tests and platforms with minimal regulatory and economic effort. Unfortunately, the European Commission's new In Vitro Diagnostic Regulation (IVDR 2017/746) will further raise the regulatory burden (from 2022 on), which may discourage companies from frequently updating assays. Thus, this regulation might, paradoxically, endanger the patients it was intended to protect.

Finally, it is crucial to monitor the emergence of new SCC*mec* junction sites in staphylococci, both in *S. aureus* and coagulase-negative staphylococci, as mobile SCC*mec* elements can be readily transmitted between different strains and species, as was the case in the strain described herein. It is imperative that companies that offer rapid tests should react to these changes and adapt their products to the market as quickly as possible.

AUTHORS' CONTRIBUTIONS

Stefan Monecke: Conception and design of the study, data visualisation, analysis and interpretation of data, drafting of the manuscript

Elisabeth König: Collection of strains and data generation, analysis and interpretation of data

Megan R. Earls: Data generation, analysis and interpretation of data; drafting of the manuscript

Eva Leitner: Conception and design of the study, data generation, analysis and interpretation of data, drafting of the manuscript

Elke Müller: Data generation, analysis and interpretation of data

Gabriel Wagner: Data generation, analysis and interpretation of data

David Poitz: Data generation, analysis and interpretation of data

Lutz Jatzwauk: Data generation, analysis and interpretation of data

Teodora Vremeră: Collection of strains and data generation **Olivia S. Dorneanu:** Collection of strains and data generation

Alexandra Simbeck: Collection of strains and data generation, analysis and interpretation of data

Andreas Ambrosch: Data generation, analysis and interpretation of data

Ines Zollner-Schwetz: Data generation, analysis and interpretation of data

Robert Krause: Data generation, analysis and interpretation of data

Werner Ruppitsch: Data generation, analysis and interpretation of data
Wulf Schneider: Data generation, analysis and interpretation of data

David C. Coleman: Analysis and interpretation of data; drafting of the manuscript **Ivo Steinmetz:** Analysis and interpretation of data; drafting of the manuscript

Ralf Ehricht: Conception and design of the study, analysis and interpretation of data; drafting of the manuscript

ACKNOWLEDGMENTS

DC and ME wish to acknowledge the support of the staff of the Irish National MRSA Reference Laboratory at St. James's Hospital, Dublin, Ireland.

COMPETING INTEREST

Not applicable.

ETHICS APPROVAL AND CONSENT TO PARTICIPATE

Not applicable.

FUNDING

There was no external funding for this study.

REFERENCES

- 1. Shore AC, Deasy EC, Slickers P, Brennan G, O'Connell B, Monecke S, et al. Detection of staphylococcal cassette chromosome *mec* type XI carrying highly divergent mecA, mecI, mecR1, blaZ, and ccr genes in human clinical isolates of clonal complex 130 methicillin-resistant *Staphylococcus aureus*. Antimicrob Agents Chemother. 2011;55(8):3765-73.
- 2. Garcia-Alvarez L, Holden MT, Lindsay H, Webb CR, Brown DF, Curran MD, et al. Meticillin-resistant *Staphylococcus aureus* with a novel mecA homologue in human and bovine populations in the UK and Denmark: a descriptive study. Lancet Infect Dis. 2011;11(8):595-603.
- 3. Monecke S, Gavier-Widen D, Mattsson R, Rangstrup-Christensen L, Lazaris A, Coleman DC, et al. Detection of mecC-positive *Staphylococcus aureus* (CC130-MRSA-XI) in diseased European hedgehogs (*Erinaceus europaeus*) in Sweden. PLoS ONE. 2013;8(6):e66166.
- 4. Rasmussen SL, Larsen J, van Wijk RE, Jones OR, Berg TB, Angen Ø, et al. European hedgehogs (*Erinaceus europaeus*) as a natural reservoir of methicillin-resistant *Staphylococcus aureus* carrying *mecC* in Denmark. PLoS ONE. 2019;14(9):e0222031.
- 5. Bengtsson B, Persson L, Ekström K, Unnerstad HE, Uhlhorn H, Börjesson S. High occurrence of *mecC*-MRSA in wild hedgehogs (*Erinaceus europaeus*) in Sweden. Vet Microbiol. 2017;207:103-7.
- 6. Becker K, van Alen S, Idelevich EA, Schleimer N, Seggewiss J, Mellmann A, et al. Plasmid-Encoded Transferable *mecB*-Mediated Methicillin Resistance in *Staphylococcus aureus*. Emerg Infect Dis. 2018;24(2):242-8.
- 7. Schwendener S, Cotting K, Perreten V. Novel methicillin resistance gene *mecD* in clinical *Macrococcus caseolyticus s*trains from bovine and canine sources. Scientific Reports. 2017;7(1):43797.
- 8. Xu Z, Hou Y, Peters BM, Chen D, Li B, Li L, et al. Chromogenic media for MRSA diagnostics. Molecular Biology Reports. 2016;43(11):1205-12.
- 9. van Belkum A, Rochas O. Laboratory-Based and Point-of-Care Testing for MSSA/MRSA Detection in the Age of Whole Genome Sequencing. Frontiers in Microbiology. 2018;9(1437).
- 10. Shore AC, Rossney AS, O'Connell B, Herra CM, Sullivan DJ, Humphreys H, et al. Detection of staphylococcal cassette chromosome *mec*-associated DNA segments in multiresistant methicillin-susceptible *Staphylococcus aureus* (MSSA) and identification of *Staphylococcus epidermidis ccrAB*4 in both methicillin-resistant *S. aureus* and MSSA. Antimicrob Agents Chemother. 2008;52(12):4407-19.
- 11. Monecke S, Jatzwauk L, Muller E, Nitschke H, Pfohl K, Slickers P, et al. Diversity of SCC*mec* elements in *Staphylococcus aureus* as observed in South-Eastern Germany. PLoS ONE. 2016;11(9):e0162654.
- 12. Sharff KA, Monecke S, Slaughter S, Forrest G, Pfeiffer C, Ehricht R, et al. Genotypic Resistance Testing Creates New Treatment Challenges: Two Cases of Oxacillin-Susceptible Methicillin-Resistant *Staphylococcus aureus*. Journal of Clinical Microbiology. 2012;50(12):4151-3.
- 13. Earls MR, Kinnevey PM, Brennan GI, Lazaris A, Skally M, O'Connell B, et al. The recent emergence in hospitals of multidrug-resistant community-associated sequence type 1 and *spa* type t127 methicillin-resistant *Staphylococcus aureus* investigated by whole-genome sequencing: Implications for screening. PLoS ONE. 2017;12(4):e0175542.
- 14. Earls MR, Shore AC, Brennan GI, Simbeck A, Schneider-Brachert W, Vremerå T, et al. A novel multidrug-resistant PVL-negative CC1-MRSA-IV clone emerging in Ireland and Germany likely originated in South-Eastern Europe. Infection, Genetics and Evolution. 2019;69:117-26.
- 15. EUCAST. EUCAST rapid antimicrobial susceptibility testing (RAST) directly from positive blood culture bottles. . 2018.
- 16. EUCAST. Zone diameter breakpoints for rapid antimicrobial susceptibility testing (RAST) directly from blood culture bottles. Version 1.1. 2019.
- 17. Monecke S, Müller E, Dorneanu OS, Vremeră T, Ehricht R. Molecular typing of MRSA and of clinical *Staphylococcus aureus* isolates from Ia □i, Romania. PLoS ONE. 2014;9(5):e97833.
- 18. Monecke S, Coombs G, Shore AC, Coleman DC, Akpaka P, Borg M, et al. A field guide to pandemic, epidemic and sporadic clones of methicillin-resistant *Sta*phylococcus aureus. PLoS ONE. 2011;6(4):e17936.
- 19. Earls MR, Coleman DC, Brennan GI, Fleming T, Monecke S, Slickers P, et al. Intra-Hospital, Inter-Hospital and Intercontinental Spread of ST78 MRSA From Two Neonatal Intensive Care Unit Outbreaks Established Using Whole-Genome Sequencing. Frontiers in Microbiology. 2018;9(1485).
- 20. Chen S, Zhou Y, Chen Y, Gu J. fastp: an ultra-fast all-in-one FASTQ preprocessor. Bioinformatics. 2018;34(17):i884-i90.
- 21. Bankevich A, Nurk S, Antipov D, Gurevich AA, Dvorkin M, Kulikov AS, et al. SPAdes: a new genome assembly algorithm and its applications to single-cell sequencing. J Comput Biol. 2012;19(5):455-77.
- 22. Howden BP, Beaume M, Harrison PF, Hernandez D, Schrenzel J, Seemann T, et al. Analysis of the small RNA transcriptional response in multidrug-resistant *Staphylococcus aureus* after antimicrobial exposure. Antimicrobial agents and chemotherapy. 2013;57(8):3864-74.

- 23. Baig S, Johannesen TB, Overballe-Petersen S, Larsen J, Larsen AR, Stegger M. Novel SCC*mec* type XIII (9A) identified in an ST152 methicillin-resistant *Staphylococcus aureus*. Infection, Genetics and Evolution. 2018;61:74-6.
- 24. Bartels MD, Boye K, Rohde SM, Larsen AR, Torfs H, Bouchy P, et al. A common variant of staphylococcal cassette chromosome *mec* type IVa in isolates from Copenhagen, Denmark, is not detected by the BD GeneOhm methicillin-resistant *Staphylococcus aureus* assay. J Clin Microbiol. 2009;47(5):1524-7.
- 25. Bartels MD, Hansen LH, Boye K, Sørensen SJ, Westh H. An Unexpected Location of the Arginine Catabolic Mobile Element (ACME) in a USA300-Related MRSA Strain. PLoS ONE. 2011;6(1):e16193.
- 26. Hill-Cawthorne GA, Hudson LO, El Ghany MFA, Piepenburg O, Nair M, Dodgson A, et al. Recombinations in staphylococcal cassette chromosome mec elements compromise the molecular detection of methicillin resistance in *Staphylococcus aureus*. PLoS ONE. 2014:9(6):e101419.
- 27. Poovelikunnel T, Gethin G, Humphreys H. Mupirocin resistance: clinical implications and potential alternatives for the eradication of MRSA. Journal of Antimicrobial Chemotherapy. 2015;70(10):2681-92.
- 28. Scheithauer S, Trepels-Kottek S, Häfner H, Keller D, Ittel T, Wagner N, et al. Healthcare worker-related MRSA cluster in a German neonatology level III ICU: A true European story. International Journal of Hygiene and Environmental Health. 2014;217(2):307-11.
- 29. Manara S, Pasolli E, Dolce D, Ravenni N, Campana S, Armanini F, et al. Whole-genome epidemiology, characterisation, and phylogenetic reconstruction of *Staphylococcus aureus* strains in a paediatric hospital. Genome Medicine. 2018;10(1):82.
- 30. Loncaric I, Stalder GL, Mehinagic K, Rosengarten R, Hoelzl F, Knauer F, et al. Comparison of ESBL--and AmpC producing Enterobacteriaceae and methicillin-resistant *Staphylococcus aureus* (MRSA) isolated from migratory and resident population of rooks (*Corvus frugilegus*) in Austria. PLoS ONE. 2013;8(12):e84048.

SUPPLEMENTS

Supplemental File S1: Alignment of the SCC*mec* elements of study isolates, reference strains and the individual SCC-associated genes.

Supplemental File S2: Sequence of the SCC*mec* element of study isolate Iasi-95037 (GenBank MT380478, submission being processed)