

# **ORAL PRESENTATION**



# Free-breathing 3D whole-heart coronary mra using respiratory motion-resolved sparse reconstruction

Davide Piccini<sup>1,2\*</sup>, Li Feng<sup>3</sup>, Gabriele Bonanno<sup>1,4</sup>, Simone Coppo<sup>2,4</sup>, Jérôme Yerly<sup>2,4</sup>, Ruth P Lim<sup>5</sup>, Juerg Schwitter<sup>6</sup>, Daniel K Sodickson<sup>3</sup>, Ricardo Otazo<sup>3</sup>, Matthias Stuber<sup>2,4</sup>

*From* 19th Annual SCMR Scientific Sessions Los Angeles, CA, USA. 27-30 January 2016

## Background

Navigator gating is commonly used to minimize respiratory motion in free-breathing whole-heart coronary MRA [1]. However, lengthy and unpredictable acquisition times remain a drawback. Respiratory self-navigation (SN) [2-3], conversely, enables 100% scan efficiency, but performs motion correction over a broad range of respiratory displacements, which can result in image artifacts. Here, we propose an alternative respiratory motion-resolved approach based on 3D radial phyllotaxis sampling, respiratory motion sorting and sparse reconstruction.

## Methods

Examinations in N = 11 healthy volunteers (9 male, age:  $29 \pm 4$  y) were performed on a 1.5T clinical MRI scanner (MAGNETOM Aera, Siemens Healthcare) with a prototype 3D radial phyllotaxis bSSFP sequence [4]: TR/TE 3.1/ 1.56 ms, FOV (220 mm)<sup>3</sup>, matrix 192<sup>3</sup>, voxel (1.15 mm)<sup>3</sup>, RF angle 115°, and receiver BW 898 Hz/Px. Using a respiratory signal directly extracted from the modulations of the k-space center amplitude within the radial imaging data [5], signal-readouts were grouped according to the respiratory state at which they were acquired (Fig. 1). The resulting series of undersampled respiratory states were reconstructed using an eXtra-Dimensional Golden-angle RAdial Sparse Parallel imaging (XD-GRASP) [6] algorithm, which exploits sparsity along the newly created respiratory dimension. Datasets for 4 respiratory states were reconstructed. Image quality of the end-expiratory phase was compared to 1D respiratory self-navigation in terms of vessel sharpness (VS) [7], visible length (VL) and diagnostic quality on a scale from 0 (non-visible) to 2 (diagnostic).

### Results

Respiratory-resolved XD-GRASP reconstruction effectively suppresses respiratory motion artifacts (Fig. 1). Average VS and VL were always superior for the respiratory-resolved datasets, reaching statistical significance (p < 0.05) for the left main (LM), for the proximal and mid left anterior descending artery (LAD) (e.g. VS of mid LAD 40.8  $\pm$  9.1% vs 34.9  $\pm$  10.2%) and for the mid right coronary artery (RCA). Visualized length of LM+LAD was significantly increased as well. A total of 41/88 coronary segments were graded as diagnostic for 1D SN, while this ratio increased to 61/88 for the XD-GRASP reconstruction (Tab.1). The XD-GRASP reconstruction reached 100% diagnostic quality for LM, proximal-LAD, and proximal-RCA.

## Conclusions

Instead of discarding data or enforcing motion models for motion correction, XD-GRASP makes constructive use of all respiratory phases to improve image quality, and achieves superior quality compared to 1D respiratory SN without the need for breath-holding, navigators, or complex 3D respiratory motion correction schemes. The phyllotaxis trajectory and XD-GRASP reconstruction provide a synergistic combination that may lead routine coronary MRA closer to clinical practice.

Full list of author information is available at the end of the article



© 2016 Piccini et al. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (http:// creativecommons.org/licenses/by/4.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/ zero/1.0/) applies to the data made available in this article, unless otherwise stated.

<sup>&</sup>lt;sup>1</sup>Advanced Clinical Imaging Technology, Siemens Healthcare, Lausanne, Switzerland



#### Table 1 Diagnostic guality grading of all coronary segments

Coronary Segment	1D Respiratory Self-Navigation	4-Phase X-D GRASP (End-exp)
Left Main	1.8 ± 0.4	2.0 ± 0.0*
LAD Prox.	1.6 ± 0.5	2.0 ± 0.0*
LAD Mid	1.3 ± 0.6	$1.4 \pm 0.5$
LAD Dist.	$0.9 \pm 0.5$	$1.3 \pm 0.5$
LCX Prox.	1.4 ± 0.7	$1.4 \pm 0.7$
RCA Prox.	$1.8 \pm 0.4$	$2.0 \pm 0.0$
RCA Mid	$1.3 \pm 0.5$	$1.7 \pm 0.5$
RCA Dist.	$1.4 \pm 0.7$	$1.7 \pm 0.5$
Total Diagnostic Segments	41/88 (47%)	61/88 (70%)

All values are expressed as mean ± one standard deviation

\* Indicates statistical significance compared to 1D Respiratory Self-Navigation.

Diagnostic Grading: 0 = non-visible, 1 = visible but non diagnostic and 2 = visible and diagnostic

#### Authors' details

<sup>1</sup>Advanced Clinical Imaging Technology, Siemens Healthcare, Lausanne, Switzerland. <sup>2</sup>Department of Radiology, University Hospital (CHUV) and University of Lausanne (UNIL), Lausanne, Switzerland. <sup>3</sup>Department of Radiology, Center for Advanced Imaging Innovation and Researc (CAI2R) and Bernard and Irene Schwartz Center for Biomedical Imaging, New York, NY, USA. <sup>4</sup>Center for Biomedical Imaging (CIBM), Lausanne, Switzerland. <sup>5</sup>Department of Radiology, Austin Health and The University of Melbourne, Melbourne, VIC, Australia. <sup>6</sup>Division of Cardiology and Cardiac MR Center, University Hospital of Lausanne (CHUV), Lausanne, Switzerland.

#### Published: 27 January 2016

#### doi:10.1186/1532-429X-18-S1-O105

**Cite this article as:** Piccini *et al.*: **Free-breathing 3D whole-heart coronary mra using respiratory motion-resolved sparse reconstruction.** *Journal of Cardiovascular Magnetic Resonance* 2016 **18**(Suppl 1):O105.

## Submit your next manuscript to BioMed Central and take full advantage of: • Convenient online submission • Thorough peer review

BioMed Central

- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at www.biomedcentral.com/submit