



Cold Spring Harbor Protocols

Preparation of Samples for Single-Worm Tracking

Eviatar Yemini, Rex A. Kerr and William R. Schafer

Cold Spring Harb Protoc 2011; doi: 10.1101/pdb.prot066993

Email Alerting Service

Receive free email alerts when new articles cite this article - [click here](#).

Subject Categories

Browse articles on similar topics from *Cold Spring Harbor Protocols*.

[Behavioral Assays](#) (38 articles)
[C. elegans](#) (30 articles)
[Image Analysis](#) (61 articles)
[Imaging for Neuroscience](#) (185 articles)
[Imaging/Microscopy, general](#) (499 articles)
[In Vivo Imaging](#) (215 articles)
[Video Imaging / Time Lapse Imaging](#) (143 articles)

To subscribe to *Cold Spring Harbor Protocols* go to:
<http://cshprotocols.cshlp.org/subscriptions>

Protocol

Preparation of Samples for Single-Worm Tracking

Eviatar Yemini, Rex A. Kerr, and William R. Schafer

Neurobiological research in genetically tractable organisms relies heavily on robust assays for behavioral phenotypes. The simple body plan of the nematode *Caenorhabditis elegans* makes it particularly amenable to the use of automated microscopy and image analysis to describe behavioral patterns quantitatively. This protocol first describes the preparation and use of media for growing and maintaining worms for tracking. The second part of the protocol describes how to prepare a single young adult worm for recording during video analysis. Although the protocol was developed for use in a single-worm tracker, it addresses factors important for the generation of reproducible, standardized images in all systems.

MATERIALS

It is essential that you consult the appropriate Material Safety Data Sheets and your institution's Environmental Health and Safety Office for proper handling of equipment and hazardous materials used in this protocol.

RECIPE: Please see the end of this article for recipes indicated by <R>. Additional recipes can be found online at <http://cshprotocols.cshlp.org/site/recipes>.

Reagents

Nematode adult worms

Nematode growth medium (NGM) <R>

NGM, low-peptone

Follow the recipe for NGM with the following modifications: increase the agar to 20 g and decrease the peptone to 0.13 g.

OP50 *Escherichia coli*

Equipment

Imaging setup, including stage with platform to support the sample

The hardware setup, including stage, camera, and light source, is shown in Figure 1.

Petri dishes, ≥ 5 cm

Worm tracking hardware and software package (e.g., Worm Tracker 2.0, available at <http://www.mrc-lmb.cam.ac.uk/wormtracker/>)

Adapted from *Imaging in Neuroscience* (ed. Helmchen and Konnerth). CSHL Press, Cold Spring Harbor, NY, USA, 2011.

© 2011 Cold Spring Harbor Laboratory Press

Cite this article as *Cold Spring Harbor Protoc*; 2011; doi:10.1101/pdb.prot066993




FIGURE 1. Hardware: (A) The stage, composed of x - and y -axis actuators and linear translation stages; (B) the camera; (C) the sample, a worm on an agar Petri dish; (D) a diffuser above a collimating Fresnel lens, to provide uniform illumination; (E) a high-intensity red light-emitting diode (LED) for illumination; (F) course and fine focusing knobs to focus the sample in view of the camera; (G) a rigid cage to center the illumination above the camera and eliminate vibrations caused by moving the stage.

METHOD

Maintaining Strains

1. Prepare and pour fresh NGM into Petri dishes, and allow them to dry for 2 d with the lid on. (After drying, these NGM Petri dishes can be stored at 4°C for several days to retain their moisture; although, before use, they must be left out to reach room temperature.)
2. Seed each dish with three drops (~150 μ L) of fresh OP50 bacteria. Leave them to dry, for 2 d, with the lid on.
3. Select six adult worms and transfer onto a dish.
4. Incubate the plate at 22°C.
5. Repeat Steps 1–4 every 3 d to maintain your strains.

Tracking

6. Prepare and pour fresh NGM into Petri dishes, and allow them to dry for 2 d with the lid on. (After drying, these NGM Petri dishes can be stored at 4°C for several days to retain their moisture; although, before use, they must be left out to reach room temperature.)
7. Seed each dish with three drops (~150 μ L) of fresh OP50 bacteria. Leave them to dry, for 2 d, at room temperature with the lid on.
8. Approximately 16–18 h before tracking, select L4 worms from plates prepared in Steps 1–5 and transfer them to an NGM dish.
9. Store the plates at 22°C.
10. At 2 d before the experiment, prepare freshly poured low-peptone NGM dishes (these produce thin food lawns that produce less distortion in worm images). Leave them to dry, for 2 d, with the lid on.
11. At 30 min before each tracking experiment, seed each low-peptone dish with 20 μ L of fresh OP50 bacteria. Leave them to dry with the lid on (30 min).
12. Select and transfer one worm to a plate (prepared in Step 11). Move the plate under the tracker. Wait 30 min for the worm to habituate, and record for 15 min at room temperature.
13. Record at least 20 worms (5 h of data). When performing multiple experiments on multiple trackers, they should be randomized in time and in choice of tracker. This helps avoid inconsistent spatial and temporal cues influencing behavior by driving their statistical effects into noise rather than amplifying them as signal. For similar reasons, it is best to space experiments out over several days using different experimenters.


DISCUSSION

Background

Although most early work on *C. elegans* behavior was based on real-time observation, in recent years there has been increasing use of automated microscopy and video imaging analysis systems (known as “worm trackers”) for behavioral assays. These approaches offer a number of advantages over nonautomated methods. For example, many behavioral phenotypes have initially been described only in qualitative, imprecise terms (e.g., “weak kinker,” “loopy,” “coiler”). More quantitative assays lend precision to these descriptions and make it feasible to break behaviors down into components that can be related to particular neural functions. Another advantage is increased sensitivity, especially with regard to behaviors that occur over long timescales. There is also the potential for standardization between laboratories; if the same algorithms are used to score a behavior, data collected at different locations can in principle be compared as long as recording conditions are reproducibly controlled. Finally, an automated system (in particular a multiworm tracking microscope) offers the possibility of conducting high-throughput screens for mutants with complex behavioral phenotypes.

Worm trackers fall into two general types. Single-worm trackers follow individual animals at relatively high magnification, keeping them in the field of view by moving the camera or a motorized stage. Such systems have been used to study aspects of locomotion such as turning and speed (Pierce-Shimomura et al. 1999) as well as body posture and morphology (Karbowski et al. 2008). Multiworm trackers track multiple individuals in populations, by necessity at lower magnification, and usually focus on locomotion features. Generally speaking, single-worm trackers allow one to measure a greater number of parameters, whereas multiworm trackers offer much higher efficiency and throughput. Examples of single-worm trackers can be found in several published studies (e.g., Waggoner et al. 1998; Pierce-Shimomura et al. 1999; Hardaker et al. 2001; Feng et al. 2004; Geng et al. 2004; Cronin et al. 2005; Stephens et al. 2008). Another single-worm system, Worm Tracker 2.0, can be found on our website <http://www.mrc-lmb.cam.ac.uk/wormtracker/>. Imaging output from Worm Tracker 2.0.3.1 is shown in Figure 2.

A protocol is also available for **Illumination for Worm Tracking and Behavioral Imaging** (Yemini et al. 2011a). A discussion of various methods for recording many worms simultaneously, including a protocol, can be found in **Tracking Movement Behavior of Multiple Worms on Food** (Yemini et al. 2011b). Other multiworm trackers have also been described (Ramot et al. 2008; Tsechenakis et al. 2008).



Principles of Single-Worm Tracking

The benefit of single-worm over multiworm tracking is that it delivers videos with much higher spatial resolution—resolutions that enable neuronal tracking as well. Moreover, single-worm tracking opens up the door for stimulus targeting to specific locations of a worm on specific behavioral conditions. As such, image processing and computer vision to discover the worm must run as fast as possible. Among the simplest methods for worm localization, and the one seen most frequently in tracking software, is the following.

1. Incoming video is reduced to grayscale as worms present no additional information within color channels.
2. Video frames are dealt with individually and, using either a fixed or adaptive pixel value threshold, they are converted to a binary image of foreground and background pixels.
3. A connected-components algorithm (also known as particle or blob analysis) is run to label each group of connected foreground pixels with a separate numerical descriptor. The term “connected” usually refers to either four-connected pixels that touch only adjacently or eight-connected pixels that are permitted to touch diagonally as well.
4. Once the connected foreground pixels have been labeled, the “worm” is chosen to be either the largest group of connected pixels or a group that has an appropriate size and/or shape.
5. This “worm” is then tracked by moving the stage to keep it in view.

In general, refinements on the discovery algorithm described in point 3 include subtracting frames to determine motion and moving the stage to keep this motion in view. Motion tracking is necessary when illumination or other occlusions prevent accurate worm discovery using the more traditional approach of connected-components analysis. Additional refinements involve image processing to determine the worm’s contour, skeleton, or other descriptors that allow one to determine worm body part locations and/or behavior. These refinements are currently ad hoc as they can quickly approach the limit of what modern computers can perform in real time.



RECIPE

Nematode Growth Medium (NGM)

Reagent	Amount to add (for 1 L)	Final concentration
Agar	17.0 g	1.7% (w/v)
NaCl	2.9 g	50 mM
Peptone	2.5 g	0.25% (w/v)
CaCl ₂ (1 M)	1 mL	1 mM
Cholesterol (5 mg/mL)	1 mL	5 µg/mL
KH ₂ PO ₄ (1 M)	25 mL	25 mM
MgSO ₄ (1 M)	1 mL	1 mM

Mix the first three reagents and autoclave. After the mixture is cool, add the last four reagents.

REFERENCES

- Cronin CJ, Mendel JE, Mukhtar S, Kim YM, Stirbl RC, Bruck J, Sternberg PW. 2005. An automated system for measuring parameters of nematode sinusoidal movement. *BMC Genet* 6: 5.
- Feng Z, Cronin CJ, Wittig JH Jr, Sternberg PW, Schafer WR. 2004. An imaging system for standardized quantitative analysis of *C. elegans* behavior. *BMC Bioinformatics* 5: 115.
- Geng W, Cosman P, Berry CC, Feng Z, Schafer WR. 2004. Automatic tracking, feature extraction and classification of *C. elegans* phenotypes. *IEEE Trans Biomed Eng* 51: 1811–1820.
- Hardaker LA, Singer E, Kerr R, Zhou G, Schafer WR. 2001. Serotonin modulates locomotory behavior and coordinates egg-laying and movement in *Caenorhabditis elegans*. *J Neurobiol* 49: 303–313.

- Karbowsky J, Schindelman G, Cronin CJ, Seah A, Sternberg PW. 2008. Systems level circuit model of *C. elegans* undulatory locomotion: Mathematical modeling and molecular genetics. *J Comput Neurosci* 24: 253–276.
- Pierce-Shimomura JT, Morse TM, Lockery SR. 1999. The fundamental role of pirouettes in *Caenorhabditis elegans* chemotaxis. *J Neurosci* 19: 9557–9569.
- Ramot D, Johnson BE, Berry TL Jr, Carnell L, Goodman MB. 2008. The Parallel Worm Tracker: A platform for measuring average speed and drug-induced paralysis in nematodes. *PLoS One* 3: e2208.
- Stephens GJ, Johnson-Kerner B, Bialek W, Ryu WS. 2008. Dimensionality and dynamics in the behavior of *C. elegans*. *PLoS Comput Biol* 4: e1000028.
- Tsechpenakis G, Bianchi L, Metaxas D, Driscoll M. 2008. A novel computational approach for simultaneous tracking and feature extraction of *C. elegans* populations in fluid environments. *IEEE Trans Biomed Eng* 55: 1539–1549.
- Waggoner LE, Zhou GT, Schafer RW, Schafer WR. 1998. Control of alternative behavioral states by serotonin in *Caenorhabditis elegans*. *Neuron* 21: 203–214.
- Yemini E, Kerr RA, Schafer WR. 2011a. Illumination for worm tracking and behavioral imaging. *Cold Spring Harb Protoc* doi: 10.1101/pdb.prot067009.
- Yemini E, Kerr RA, Schafer WR. 2011b. Tracking movement behavior of multiple worms on food. *Cold Spring Harb Protoc* doi: 10.1101/pdb.prot067025.