

**LETTER**

**Reduced As components in highly oxidized environments:  
Evidence from full spectral XANES imaging using the Maia massively parallel detector**

**B.E. ETSCHMANN,<sup>1,\*</sup> C.G. RYAN,<sup>2,3,4</sup> J. BRUGGER,<sup>1,5</sup> R. KIRKHAM,<sup>6</sup> R.M. HOUGH,<sup>2</sup> G. MOORHEAD,<sup>3,6</sup>  
D.P. SIDONS,<sup>7</sup> G. DE GERONIMO,<sup>8</sup> A. KUCZEWSKI,<sup>8</sup> P. DUNN,<sup>6</sup> D. PATERSON,<sup>9</sup> M.D. DE JONGE,<sup>9</sup>  
D.L. HOWARD,<sup>9</sup> P. DAVEY,<sup>6</sup> AND M. JENSEN<sup>6</sup>**

<sup>1</sup>Earth and Environmental Sciences, University of Adelaide, South Australia 5005, Australia

<sup>2</sup>CSIRO, Exploration and Mining, Clayton, Victoria 3168, Australia

<sup>3</sup>School of Physics, University of Melbourne, Parkville, Victoria 5010, Australia

<sup>4</sup>CODES Centre of Excellence, University of Tasmania, Hobart, Tasmania 7001, Australia

<sup>5</sup>Mineralogy, South Australian Museum Adelaide, South Australia 5000, Australia

<sup>6</sup>CSIRO Materials Science and Engineering, Clayton, Victoria 3168, Australia

<sup>7</sup>National Synchrotron Light Source, Brookhaven National Laboratory, Upton, New York 11973, U.S.A.

<sup>8</sup>Instrumentation Division, Brookhaven National Laboratory, Upton, New York 11973, U.S.A.

<sup>9</sup>Australian Synchrotron, Clayton, Victoria 3168, Australia

**ABSTRACT**

Synchrotron X-ray fluorescence (SXRF) and X-ray absorption spectroscopy (XAS) have become standard tools to measure element concentration, distribution at micrometer- to nanometer-scale, and speciation (e.g., nature of host phase; oxidation state) in inhomogeneous geomaterials. The new Maia X-ray detector system provides a quantum leap for the method in terms of data acquisition rate. It is now possible to rapidly collect fully quantitative maps of the distribution of major and trace elements at micrometer spatial resolution over areas as large as  $1 \times 5 \text{ cm}^2$ . Fast data acquisition rates also open the way to X-ray absorption near-edge structure (XANES) imaging, in which spectroscopic information is available at each pixel in the map. These capabilities are critical for studying inhomogeneous Earth materials. Using a 96-element prototype Maia detector, we imaged thin sections of an oxidized pisolitic regolith ( $2 \times 4.5 \text{ mm}^2$  at  $2.5 \times 2.5 \text{ }\mu\text{m}^2$  pixel size) and a metamorphosed, sedimentary exhalative Mn-Fe ore ( $3.3 \times 4 \text{ mm}^2$  at  $1.25 \times 5 \text{ }\mu\text{m}^2$ ). In both cases, As *K*-edge XANES imaging reveals localized occurrence of reduced As in parts of these oxidized samples, which would have been difficult to recognize using traditional approaches.

**Keywords:** Arsenic, oxidation state, XANES, element distribution, imaging, X-ray fluorescence