# **Open Access**

# Similarity assessment by multivariate statistics method based on the distance between biosimilar and originator



Jian Xu<sup>1</sup>, Zhihui Shao<sup>2</sup>, Xiaoxiong Han<sup>3</sup>, Yingfeng Huang<sup>4</sup>, Xun Zou<sup>5</sup> and Yaling Shen<sup>1\*</sup>

## Abstract

The development of biosimilar products or follow-on biologics has been flourishing in recent years because of their lower price than the originators. In this study, a multivariate data analysis method based on JMP software was proposed to assess the glycosylation pattern similarity of antibody candidates from different conditions in optimization experiments with a reference. A specific distance was generated by this method and indicated the glycoform similarity between the biosimilar and the reference. This method can be applied to analyze the similarity of other physicochemical and functional characteristics between follow-on biologics and originators. Then, the design of experimental methods can be realized to optimize the conditions of cell culture to attain similar antibody candidates. A higher concentration of GlcNAc added to the basal media made the glycan of the antibody more similar to the glycan of the reference in this study.

Keywords: Biosimilar, Cluster analysis, Glycosylation pattern, Multivariate statistics method, Similarity

## Introduction

During the past two decades, biologic products (also termed "biopharmaceuticals") have been developed intensively by companies for treating cancer and autoimmune diseases such as rheumatoid arthritis. Six of the top pharmaceutical product sales in 2019 were biologic products (Table S1, Additional file 1). Therefore, there has been an increasing trend toward the development of biosimilars, considering market interest and improved access to biologics. Biosimilar products can offer lower prices to optimize efficiencies across healthcare systems. However, biopharmaceuticals typically involve the expression of the gene using living cells followed by purification and formulation to acquire stable drug products. They are large, complex and heterogeneous compared to

<sup>1</sup> State Key Laboratory of Bioreactor Engineering, Shanghai

Collaborative Innovation Center for Biomanufacturing Technology, East China University of Science and Technology, 130 Mei Long Road,

Shanghai 200237, China Full list of author information is available at the end of the article



small molecule products. It is impossible to manufacture identical copies of biologic products (WHO 2016).

To define biosimilars, different regulatory bodies have different criteria. EMA defines it as a biotherapeutic product that is similar in terms of quality, safety, and efficacy to an already licensed reference product (reference medicinal product) in the EMA (EMA 2015). FDA defines it as a biological product that is highly similar to a US-licensed reference product notwithstanding minor differences in clinically inactive components and for which there are no clinically meaningful differences between the biological product and the reference product in terms of safety, purity, and potency of the product (FDA 2015). WHO defines it as a biotherapeutic product that is similar in terms of quality, safety, and efficacy to an already licensed reference product (WHO 2016). In broad terms, a biosimilar is highly similar to a reference product in terms of structure and function. High similarities of physicochemical and functional characteristics are the main aim in the preclinical development phase.

© The Author(s) 2021. This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativecommons.org/licenses/by/4.0/.

<sup>\*</sup>Correspondence: ylshen@ecust.edu.cn

**Table 1** Design of the experiment for investigating the glycan

 profile adjustment in response to supplement addition

Run	GlcNAc (mM) <sup>a</sup>	MnCl2 (µM) <sup>a</sup>	Add in Basal or Feed Media <sup>a</sup>	Distance <sup>b</sup>
SF-1	120	30	Feed	6.71
SF-2	80	30	Feed	10.62
SF-3	40	30	Feed	12.75
SF-4	120	15	Feed	8.94
SF-5	80	15	Feed	11.27
SF-6	40	15	Feed	18.93
SF-7	120	0	Feed	6.44
SF-8	80	0	Feed	6.96
SF-9	40	0	Feed	9.42
SF-10	20	0	Basal	6.06
SF-11	10	0	Basal	5.45
SF-12	5	0	Basal	6.10
SF-13	10	30	Basal	5.42
SF-14	10	15	Basal	4.25
SF-15	120	30	Feed	8.22
SF-16	80	30	Feed	9.74
SF-17	40	30	Feed	18.50
SF-18	120	15	Feed	7.78
SF-19	80	15	Feed	11.20
SF-20	40	15	Feed	20.79
SF-21	120	0	Feed	5.80
SF-22	80	0	Feed	6.40
SF-23	40	0	Feed	9.65
SF-24	20	0	Basal	5.62
SF-25	10	0	Basal	5.64
SF-26	5	0	Basal	5.48
SF-27	10	30	Basal	4.44
SF-28	10	15	Basal	4.51
SF-29	0	0	N/A	16.97
SF-30	0	0	N/A	16.95

 $\mathsf{SF}\xspace{-}\mathsf{SF}\xspace{-}\mathsf{14}$  and  $\mathsf{SF}\xspace{-}\mathsf{15}\xspace{-}\mathsf{SF}\xspace{-}\mathsf{28}$  were duplicate runs, while  $\mathsf{SF}\xspace{-}\mathsf{29}$  and  $\mathsf{SF}\xspace{-}\mathsf{30}$  were controls

<sup>a</sup> Three factors (GlcNAc additive concentration, MnCl<sub>2</sub> additive concentration and the addition mode) were investigated in this experiment. Three levels of GlcNAc concentration in the basal media or feed media were tested. Three MnCl<sub>2</sub> concentrations were tested in the basal or feed media. The experiment was conducted in a shaker flask culture system.

 $^{\rm b}\,$  The column is the distance value from cluster analysis, which are the response of the experiment

The structural and functional elements of therapeutic antibodies include the primary structure, purity, charge heterogeneity, glycosylation and other posttranslational modifications, as well as target and receptor binding activity and bioactivity features (Kirchhoff et al. 2017). A variety of analytic techniques have been developed to demonstrate these elements. However, many analytic results are multivariate data, such as charge heterogeneity, glycosylation and size heterogeneity. It is difficult to assess the similarity between biosimilars and originators (or references) based on these multivariate outcomes.

Identification of the glycosylation pattern is a key consideration during the development of monoclonal antibody (mAb) biosimilars, since the glycan chains in the Fc region can substantially alter protein activity and the PK profile and, in some cases, antigenicity (Kirchhoff et al. 2017). Because many cell culture conditions, such as nutrient availability, pH, dissolved oxygen (DO), ammonia, cell viability, growth phase, and temperature, will affect glycosylation, the culture process parameters should be well monitored and controlled during the manufacturing phase (Patrick et al. 2009). Chemical supplements such as metal ions and substrates for glycan chain synthesis have been described in the literature as effective glycosylation modulators during upstream process development (Crowell et al. 2007; Gramer et al. 2011). To find the optimal additive amount, the design of the experiment (DOE) is used intensively during the screening experiments. However, a specific parameter indicating similarity should be defined as the response before modeling and ANOVA. An effective method that can output a specific parameter to assess the similarity of glycan profiles is needed.

In this work, a multivariate data analysis method was first applied to generate a specific index which can represent the similarity of glycan profiles from different antibody candidates. By this method, similar antibodies were clustered by glycan distribution, and the most similar antibody was easily identified. The effect of supplements could be quantified statistically.

#### Materials and methods

#### **Cell line and reagents**

The cell line expressing recombinant IgG1 was derived from CHO DG44. The basal medium was CD DF1 (Shanghai BasalMedia Technologies Co., LTD., Shanghai, China) supplemented with 6 mM glutamine (Sigma, Shanghai, China). Feed media was Efficient Feed<sup>TM</sup> C+(Thermo Fisher, Shanghai, China), and additives were added as indicated in the design. The additives N-acetyl-D-glucosamine (GlcNAc) and MnCl<sub>2</sub> were purchased from Sigma–Aldrich Shanghai Trading Co., Ltd. (Shanghai, China).

#### Cell culture conditions and process design

The cells were cultured in 125 ml shaker flasks (Corning, New York, USA) with vented caps at 5%  $CO_2$ , 70% humidity and 37 °C (shift to 35 °C from day 4 to harvest). The pH of the basal media was 7.0–7.2. The volume of the initial culture was 30 mL with a shaking rate of 120 rpm and a rotated diameter of 50 mm. The inoculum density was adjusted to 1 million cells/ml. GlcNAc and MnCl<sub>2</sub>

were supplemented into the culture in two ways: as an addition to the feed media or as an addition separately to the basal media on day 4. The supplemented amount was designed, as indicated in Table 1. The first 14 runs and second 14 runs were duplicated, and the last 2 runs were control runs without supplementation. The feed started on day 3, and the feed volume was calculated by Eq. (1):

$$V_{\text{feed}} = V_{\text{current}} \times Q_{\text{feed}} \times \frac{(\text{VCD}_1 + \text{VCD}_2)}{2 \times \text{Glc}_{\text{feed}} \times 1000}$$
(1)

where  $V_{\text{feed}}$  (mL) is the feed volume of the current day;  $V_{\text{current}}$  (mL) is the culture volume before feeding;  $Q_{\text{feed}}$ (g/10<sup>9</sup> cells) is the specific glucose consumption rate in one day, which is 75 in this study; VCD<sub>1</sub> (10<sup>6</sup> cells/mL) is the viable cell density of the current day; VCD<sub>2</sub> (10<sup>6</sup> cells/ mL) is the predicted viable cell density of the next day; and Glc<sub>feed</sub> (g/L) is the glucose concentration in the feed media. The predicted viable cell density is double the current viable cell density in the exponential growth phase and equal to the current viable cell density in the steady growth phase.

Sampling was conducted every day from day 3. The residual glucose was analyzed by a glucose test kit (Beihai, China). Additive glucose was calculated by Eq. (2) and added to maintain the concentration of glucose at 3 g/L when the viable cell density was less than  $1.0 \times 10^7$  cells/mL and shifted to 4 g/L when the viable cell density was higher than  $1.0 \times 10^7$  cells/mL:

- 20 °C. After centrifugation at 13,000 rpm for 15 min, the supernatant was vacuum dried.

#### Fluorescence labeling of glycan chains

5 mg 8-aminopyrene-1,3,6-trisulfonic acid trisodium salt (APTS, AB Sciex Pte. Ltd., Framingham, USA) were added into 0.5 mL aqueous solution containing 15% acetic acid and vortexed it. Then, 15  $\mu$ l APTS solution and 5  $\mu$ L tetrahydrofuran solution containing 1 mol/L sodium cyanoborohydride were added to the vacuum-dried gly-can chains. After fluorescent labeling at 55 °C for 2 h, 400  $\mu$ l ultrapure water was added and the product was analyzed by a PA800 Plus capillary electrophoresis apparatus (AB Sciex Pte. Ltd., Framingham, USA).

#### Analysis by a capillary electrophoresis apparatus

A Beckman N-Cho-coated capillary with a total length of 60.5 cm, effective length of 50 cm and inner diameter of 50  $\mu$ m, and the electrophoresis buffer were purchased from Beckman Coulter Life Sciences (Indianapolis, USA). The capillary temperature was 20 °C. The sample was injected at 2.0 psi for 10 s. Then, the sample was separated at 30 kV for 20 min. Fluorescence detection was implemented at excitation and emission wavelengths of 488 and 520 nm, respectively.

#### **Multivariate statistics**

The peak area percentages from the glycan analysis are

$$V_{\rm Glc} = \left[ \left( {\rm Glc}_{\rm target} - {\rm Glc}_{\rm test} \right) \times V_{\rm current} + \left( {\rm Glc}_{\rm target} - {\rm Glc}_{feed} \right) \times V_{\rm feed} \right] / {\rm Glc}_{\rm solution}$$
(2)

where  $V_{\rm Glc}$  (mL) is the additive volume of glucose solution;  $V_{\rm current}$  (mL) is the culture volume before feeding;  $V_{\rm feed}$  is the volume of feed media;  ${\rm Glc}_{\rm feed}$  (g/L) is the glucose concentration in the feed media;  ${\rm Glc}_{\rm test}$  (g/L) is the glucose concentration in the culture;  ${\rm Glc}_{\rm target}$  (g/L) is the target glucose concentration after feeding; and the  ${\rm Glc}_{\rm solution}$  (g/L) is the concentration of the glucose solution.

After 11 days of culture, the culture was harvested by centrifugation to remove the cell debris and the target antibody was purified by Protein A resin.

#### **Glycan analysis**

#### Enzymolysis of glycan chains from the antibody

First, the buffer of the antibody sample was replaced with 50 mmol/L NH<sub>4</sub>HCO<sub>3</sub> (pH 8.0). Then, 500  $\mu$ g antibody was hydrolyzed with 2000 U PNGase F (New England Biolabs, Beijing, China) at 37 °C for 24 h. Precooled ethanol was added to a final concentration of 75% (v/v). The mixture was blended and allowed to stand for 0.5 h at

used for assigning the similarity in this case. Seven peaks were identified and coded by numbers for each isoform. In this way, all of the data were normalized for each sample. A matrix of the data was built, and hierarchical cluster analysis was performed with JMP software (SAS Institute Inc.). A dendrogram could be generated and the distance matrix could be saved to another data table. The distance can refer the measure of virtual distance between biosimilar to originator. Distance can also be denoted similarity. To analyze the effect of the 3 factors on glycan similarity, the distance was set as the response value and then fit model was run in JMP.

# Results

#### **Glycan analysis**

The cell growth data and residual glucose concentration variation are shown in the Additional file 1. The cells grew well, and glucose was not depleted in all culture runs. The culture was harvested on day 11, and the supernatant was collected by centrifugation. After purification and analysis as described in the methods, the glycan peak distributions of 30 samples from different shaker flask cultures are shown in Fig. 1. All of the samples' glycans showed similar patterns, with peak 3 being the highest portion; however, the percentages of each peak were different. At a glance, supplementation with GlcNAc and  $MnCl_2$  should affect the glycan distribution of antibodies. However, it is difficult to recognize the effect of the supplements, and the effect cannot be tested by a statistical method, because a specific number to represent the similarity is lacking.

#### **Cluster analysis**

A  $31 \times 8$  matrix with a 7 peak area percentage of the reference and each sample was generated. By performing the cluster analysis tool in JMP as described in the Methods, a hierarchical clustering tree is shown in Fig. 2. The single samples are the leaves, and the similar samples are clustered on one branch. We set the duplicate runs to the same color, so it can be seen that almost all duplicate samples are in the third or fourth branch from the main trunk. It is very clear that SF-7, 8, 13, 14, 21, 22, 27, and 28 are located in the same third branch as the reference, which illustrates that the glycan distribution of these samples is very close to that of the reference.

#### Distance

The distance shows the similarity between the samples and the reference. They can be generated by saving the distance matrix option in menu. The distances are plotted to run the label, as shown in Fig. 3. Here, it is very clear that the antibodies from SF-14 and 28 are the nearest to the reference antibody in terms of glycan distribution, while SF-6 and 20 are the farthest. To confirm the real peak distribution similarity, only SF-14, SF-20 and the reference are shown in Fig. 4. SF-14 is similar to the



**Hierarchical Clustering** 

Method = Ward

Dendrogram







**Table 2** ANOVA results of three factors affecting glycan similarity using the distance value as the response.

Source	Nparm	DF	Sum of Squares	F Ratio	Р
GlcNAc	1	1	173.1	27.2	<0.0001*
MnCl <sub>2</sub>	1	1	17.6	2.8	0.1083
Add in Basal or Feed Media	2	2	474.2	37.3	<0.0001*

\*Significant when P is less than 0.05

reference, and the difference between SF-20 and the reference is significant.

#### **Effect Significance Analysis**

Use the "Fit model" in JMP, we set the Distance as the response and put "GlcNAc", "MnCl<sub>2</sub>" and "adding them

in basal or feed media" in the model effect window and then ran the fit program. The ANOVA results are shown in Table 2. In the effect test, the *P* value of GlcNAc additive concentration and addition to basal or feed media were < 0.0001, which means that these 2 factors were significant, while the MnCl<sub>2</sub> additive concentration was 0.1083, which was justified as not significant. A higher concentration of GlcNAc added to the basal media made the glycan of the antibody more similar to the glycan of the reference. In this way, the significance of the factors can be identified statistically, and then the optimal condition can be quickly determined.

#### Discussion

In developing a biosimilar drug, a stepwise approach is needed beginning with chemistry manufacturing controls (CMC) and bioanalytical characterization (Burchiel et al. 2019). In the CMC development phase, because the structure of antibodies is very complex, the quality attributes will be characterized by multiple analysis methods, such as glycan isoforms, capillary isoelectric focusing (cIEF), cation-exchange chromatography (CEX) and peptide maps. We know that these methods output multivariance results rather than single data. Clustering is a ubiquitous data analysis tool to divide complex data into groups of similar items (Andreas et al. 2019). Therefore, it can be utilized to reveal the similarity of these multivariance test results effectively. Kang and Chow (2013) proposed a three-arm parallel design to assess biosimilarity between a biosimilar product and an innovator biological product based on the relative distance of means observed from the test and reference products. In the proposed design, if the relative distance is less than a prespecified margin, they claim that the two products are biosimilar. This method's merit is to provide a specific

standard to access similarity, but it is relatively complex and inexecutable compared to the method we proposed. Beyond this method, we have not found other strategies to assess biosimilar quality by a statistical pathway in the published literature.

The clustering program in JMP can output a hierarchical tree, making the cluster results visualizable and easier to find a similar group to a reference. This can be utilized to control the batch quality in the manufacturing phase. If the third branch is set as a similar margin, the leaves of different batch data can be judged as qualified batches. Once a batch's quality is clustered to another third branch, deviation investigation can be triggered to recall or destroy the batch based on the risk assessment. The distance between samples and the reference can be used as the response value in the optimization experiment. In this way, a quality-by-design (QbD) strategy is feasible to characterize the process effect and to optimize or define the operation space of the critical process parameter (CPP) by DOE. In the future, many batches will be clustered, and a distance can also be set as a margin to select the qualified batch.

GlcNAc is one of the main monosaccharides in conserved N-linked glycan structure. In the medial Golgi, the N-acetylglucosaminyltransferase-I (GnTI) enzyme mediates the transfer of GlcNAc from UDPGlcNAc to the O-2 position of the terminal mannose residue in the  $\alpha 1 \rightarrow 3$  branch of the Man5GlcNAc2 oligosaccharide (Liming 2015). GlcNAc supplement will increase the substrate of glycan synthesis and change the percentage of various glycans. In this study, the additional glycans increased the percentage of peak 3 and decreased the percentage of peak 5, 6 and 7 significantly. Crowell et al. (2007) reported that manganese is the cofactor of galactosyltransferse and increased the galactosylation. For the antibody we studied, effect of  $MnCl_2$  is not significant in our tested range but it may help GlcNAc to tune galacosylation to be more similar to reference.

#### Conclusion

In this paper, a multivariate statistics method is proposed to assess the similarity of antibodies to references from different conditions in optimization experiments. The multivariance quality results can be grouped by this method, and a specific distance can be generated. The distance value indicates the similarity between the biosimilar and the reference and the DOE method can be realized to evaluate the effects of factors and to optimize the culture conditions. In the case study, the highest similar glycoform was easily identified. GlcNAc supplemental and adding mode were significant factors and their impact were clear. The optimal conditions to attain higher similar antibody were 10 mM GlcNAc and 15  $\mu M$  MnCl\_2 supplemented into basal medium.

#### Abbreviations

WHO: World Health Organization; EMA: European Medicines Agency; FDA: U.S. Food and drug administration; mAb: Mono-clonal antibody; PK: Pharmacokinetics; DO: Dissolved oxygen; DOE: Design of experiment; ANOVA: Analysis of Variance; GlcNAc: N-Acetyl-D-glucosamine; CMC: Chemical, manufacture and control; clEF: Capillary isoelectric focusing; CEX: Cation-exchange chromatography; QbD: Quality by design; CPP: Critical process parameter.

#### **Supplementary Information**

The online version contains supplementary material available at https://doi.org/10.1186/s40643-021-00378-2.

Additional file 1: Table S1. The top 10 pharmaceutical products sold in 2019. Fig. S1. Viable cell density variation of 30 shaker flask runs by the culture time. Fig. S2 Viability variation of 30 shake flask runs by the culture time. Fig. S3 Residual glucose variation of 30 shake flask runs by the culture time.

#### Acknowledgements

The authors wish to acknowledge the Quality Analytical Science group at Dragonboat Biopharmaceutical Co. for testing the glycosylation.

#### Authors' contributions

JX wrote the main manuscript text. ZS and XH analyzed the data by statistical tools. YH, XZ and YS conducted the experiment. All authors read and approved the final manuscript.

#### Funding

Not applicable.

#### Availability of data and materials

Not applicable.

#### Declarations

**Ethics approval and consent to participate** Not applicable.

#### **Consent for publication**

Not applicable.

#### **Competing interests**

The authors confirm that there are no conflicts of interests.

#### Author details

<sup>1</sup> State Key Laboratory of Bioreactor Engineering, Shanghai Collaborative Innovation Center for Biomanufacturing Technology, East China University of Science and Technology, 130 Mei Long Road, Shanghai 200237, China. <sup>2</sup> CAS Key Laboratory of Synthetic Biology, CAS Center for Excellence in Molecular Plant Sciences, Institute of Plant Physiology and Ecology, Chinese Academy of Sciences, Shanghai, China. <sup>3</sup> Process Development Department, Dragon Sail Pharmaceutical, Shanghai, China. <sup>5</sup> Shanghai Sanjin Bioscience and Technology, Shanghai, China.

Received: 17 December 2020 Accepted: 16 March 2021 Published online: 29 March 2021

#### References

- Andreas A, Margareta A, Naomi CB (2019) To cluster, or not to cluster: An analysis of clusterability methods. Pattern Recogn 88:13–26. https://doi.org/10. 1016/j.patcog.2018.10.026
- Burchiel SW, Aspbury R, Munday J (2019) The search for biosimilars and biobetters. Drug Discovery Today 24(5):1087–1091. https://doi.org/10.1016/j. drudis.2019.03.016
- Crowell CK, Grampp GE, Rogers GN, Miller J, Scheinman RI (2007) Amino acid and manganese supplementation modulates the glycosylation state of erythropoietin in a CHO culture system. Biotechnol Bioeng 96(3):538–549. https://doi.org/10.1002/bit.21141
- EMA (2015) Guideline on similar biological medicinal products containing biotechnology-derived proteins as active substance: Non-clinical and clinical issues. http://www.ema.europa.eu/docs/en\_GB/document\_libra ry/Scientific\_guideline/2015/01/WC500180219.pdf
- FDA (2015) Guidance for industry: Scientific considerations in demonstrating biosimilarity to a reference product. http://www.fda.gov/downloads/ Drugs/GuidanceComplianceRegulatoryInformation/Guidances/UCM29 1128.pdf
- Gramer MJ, Eckblad JJ, Donahue R, Brown J, Shultz C, Vickerman K, Priem P, Bremer E, Gerritsen J, Berkel P (2011) Modulation of antibody galactosylation through feeding of uridine, manganese chloride, and galactose. Biotechnol Bioeng 108(7):1591–1602. https://doi.org/10.1002/bit.23075

- Kang SH, Chow SC (2013) Statistical assessment of biosimilarity based on relative distance between follow-on biologics. Statist Med 32(3):382–392. https://doi.org/10.1002/sim.5582
- Kirchhoff CF, Wang XZM, Conlon HD, Anderson S, Ryan AM, Bose A (2017) Biosimilars: key regulatory considerations and similarity assessment tools. Biotechnol Bioeng 114(12):2696–2705. https://doi.org/10.1002/bit.26438
- Linming L (2015) Antibody glycosylation and its impact on the pharmacokinetics and pharmacodynamics of monoclonal antibodies and Fc-fusion proteins. J Pharm Sci 104(6):1866–1884. https://doi.org/10.1002/jps.24444
- Patrick H, Khattak SF, Li ZJ (2009) Optimal and consistent protein glycosylation in mammalian cell culture. Glycobiology 19(9):936–949. https://doi.org/ 10.1093/glycob/cwp079
- WHO (2016) Guidelines on evaluation of monoclonal antibodies as similar biotherapeutic products (SBPs). https://www.who.int/biologicals/bioth erapeutics/WHO\_TRS\_1004\_web\_Annex\_2.pdf

#### **Publisher's Note**

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

# Submit your manuscript to a SpringerOpen<sup>®</sup> journal and benefit from:

- Convenient online submission
- ► Rigorous peer review
- Open access: articles freely available online
- High visibility within the field
- Retaining the copyright to your article

Submit your next manuscript at > springeropen.com