

 Open access • Posted Content • DOI:10.1101/2021.02.11.21251605

Viroselect: A novel SARS-CoV-2 detection assay to resolve inconclusive samples

— [Source link](#) 

Jawade K, Sinha Ay, Sharad Bhagat, Shilpa Bhowmick ...+14 more authors

Institutions: Kasturba Medical College, Manipal, Indian Council of Medical Research

Published on: 13 Feb 2021 - medRxiv (Cold Spring Harbor Laboratory Press)

Related papers:

- [High throughput detection and genetic epidemiology of SARS-CoV-2 using COVIDSeq next-generation sequencing.](#)
- [Quantitative detection of SARS-CoV-2 B.1.1.7 variant in wastewater by allele-specific RT-qPCR](#)
- [Rapid base-specific calling of SARS-CoV-2 variants of concern using combined RT-PCR melting curve screening and SIRPH technology](#)
- [Optimizing SARS-CoV-2 molecular diagnostic using N gene target: insights about reinfection](#)
- [Direct Quantitation of SARS-CoV-2 Using Droplet Digital PCR in Suspected Samples With Very Low Viral Load](#)

Share this paper:    

View more about this paper here: <https://typeset.io/papers/viroselect-a-novel-sars-cov-2-detection-assay-to-resolve-1mo3x45yz3>

Viroselect: A novel SARS-CoV-2 detection assay to resolve inconclusive samples

Ketki Jawade^{#1}, Akhauri Yash Sinha^{#2}, Sharad Bhagat¹, Shilpa Bhowmick¹, Bhagyashree Chauhan², Snehal Kaginkar¹, Harsha Palav¹, Nandini Kasarpalkar³, Pratik Devadiga³, Kalyani Karandikar³, Sarchie Agarwal^{4,5}, Jayanthi Shastri^{4,5}, Kiran Munne⁶, Vikrant M Bhor³, Smita D Mahale⁷, Subhanjan Bhowmik², Dhanashree Jagtap^{*7}, Vainav Patel^{*1}

These authors contributed equally

*Corresponding authors

Author Affiliations:

¹ Department of Biochemistry & Virology, ICMR-NRRH, Mumbai

² Bencos Research Solutions Pvt. Ltd., Mumbai

³ Department of Molecular Immunology & Microbiology, ICMR-NIRRH, Mumbai

⁴ Molecular Lab, Kasturba Hospital for Infectious Diseases, Mumbai

⁵ Department of Microbiology, BYL Nair Hospital, Mumbai

⁶ Department of Clinical Research, ICMR-NIRRH, Mumbai

⁷ Department of Structural Biology Division, ICMR-NIRRH, Mumbai

ABSTRACT:

Background: India bears the second largest burden of SARS-CoV-2 infection. A multitude of RT-PCR detection assays with disparate gene targets including automated high throughput platforms are available. Varying concordance and interpretation of diagnostic results in this setting can result in significant reporting delays leading to suboptimal disease management. Here, we report the development of a novel ORF-1a based SARS-CoV-2 RT-PCR assay, Viroselect, showing high concordance with conventional assays and the ability to resolve inconclusive results generated during the peak of the epidemic in Mumbai, India.

NOTE: This preprint reports new research that has not been certified by peer review and should not be used to guide clinical practice.

Methods: We identified a unique target region within SARS-CoV-2 ORF1a, non-structural protein (nsp3), that was used to design and develop our assay. This hypervariable region (1933-3956) between SARS-CoV-2, SARS-CoV, and MERS-COV was utilized to design our primers and probe for RT-PCR assay. We further evaluated concordance of our assay with commonly used EUA (USFDA) manual kits as well as an automated high throughput testing platform. Further, a retrospective analysis using Viroselect on samples reported as ‘inconclusive’ during April-October 2020 was carried out.

Results: A total of 701 samples were tested. Concordance analysis of 477 samples demonstrated high overall agreement of Viroselect assay with both manual (87.6%; 95% CI) as well as automated (84.7%; 95% CI) testing assays. Also, in the retrospective analysis of 224 additional samples reported as ‘inconclusive’, Viroselect was able to resolve 100% (19/19) and 93.7% (192/205) samples which were termed inconclusive by manual and automated high throughput platform respectively.

Conclusion: We show that Viroselect had high concordance with conventional assays, both manual and automated, as well as highlight its potential in resolving inconclusive samples.

Keywords: SARS-CoV 2; RT-PCR; inconclusive samples; Viroselect; high throughput SARS-CoV-2 testing; reporting delay

INTRODUCTION

The striking difference between previous SARS-CoV/MERS-CoV outbreaks and the ongoing SARS-CoV-2 pandemic is the increased contagious nature of the latter (Yang et al., 2020). This has resulted in over 100 million SARS-CoV-2 infections globally and ~2.2 million deaths to date (WHO COVID-19 Dashboard, 2021). India bears the second largest burden of infected individuals with 10 million infections and over 150,000 deaths so far (ICMR COVID-19 Testing Status, 2021). Efficient and rapid diagnosis is critical to optimal management of disease spread as well as clinical management of COVID-19. To ensure maximum coverage (WHO Feature Story, 2020) and accessibility to testing at economically feasible rates, Government of India, has successfully fostered and promoted development of indigenous diagnostic kits as well integrated cutting-edge high throughput testing facilities (Gupta et al., 2020; Poljak et al., 2020). The availability of a range of diagnostic kits with disparate viral targets has also posed a challenge in terms of concordance of results as well as inconclusive results (Bhattacharya et al., 2020; Hur et al., 2020; Pujadas et al., 2020). Also, implementation of high-throughput platform-based testing while rapidly increasing diagnostic capability, also has the potential of generating inconclusive results. As a designated national centre of excellence for COVID-19 diagnostic kit validation as well as a high throughput testing centre and having encountered the aforementioned issues, we developed an in-house RT-PCR assay, Viroselect, which was evaluated for concordance with commonly used EUA (USFDA) manual kits as well as with our high throughput COBAS 6800 testing platform. Further we assessed its utility in resolving inconclusive tests by both testing modalities.

MATERIALS AND METHODS

1. Viroselect Assay Development

1.1 Primer and Probe Design

A novel strategy (described in results) was used to design primers for detection of SARS-CoV-2. These sequence regions were unique to SARS-CoV-2 but had maximum variation with SARS-CoV and MERS-CoV. The primer and probe were optimized for an in-house Viroselect assay described below. All the primer and probes used in the study were synthesized from Macrogen Inc., Korea.

2. SARS-CoV-2 Detection

2.1 RNA Extraction

Oropharyngeal/ nasopharyngeal swab samples collected from suspected COVID-19 patients in Viral Transport Medium (VTM) were used for SARS-CoV-2 detection. Viral RNA Extraction from these VTM samples, were carried out using Maverick Magnetic Bead Based Nucleic Acid Extraction Kit (Cat. No.: EC00001B, MyLabs Discovery Solutions Pvt. Ltd., Lonavla, Maharashtra, India) on Thermo Scientific KingFisher Flex Purification System and MagRNA-II Viral RNA Extraction Kit (Cat. No.: G2M030420, Genes2Me Pvt. Ltd, Gurugram, Haryana, India) on MGISP-100 automated nucleic acid extraction and purification system (MGI Tech Co., Ltd, Yantian, Shenzhen, China) as per the manufacturer's instructions.

2.2 Viroselect Assay

Biplex RT-PCR assay was conducted for Viroselect assay. ORF1a (nsp3) gene for the detection of SARS-CoV-2 Virus along with human RNase P, used as internal control, from human patient samples. The primer and probe concentrations were the following, 0.5 μ mole/L for primers and 0.5 μ mole/L of probe concentration for ORF1a (nsp3) and 0.5 μ mole/L for primers and 0.2 μ mole/L of probe concentration for Human RNaseP. The TaqMan probe was labelled with 6-carboxyfluorescein (6-Fam) 5' and BHQ1 3' for ORF1a and with Cy5 at 5' end and BHQ3 at 3' end for the RNase P probe. Prime Script One-step First strand cDNA synthesis kit (Catalogue number: RR6110A, Takara Bio) along with gene specific primers (ORF1a and RNaseP) was used for cDNA preparation (as per manufacturers protocol). One step RT-PCR condition were 52° C for 5 minutes, 1 cycle (for reverse transcription), initial denaturation 95° C for 10 seconds followed by 40 cycles (for RT-PCR) of denaturation (95° C for 05 sec), (60° C for 90 sec) annealing/extension for PCR reaction.

2.3 Assay Interpretation for Viroselect

The in-house Viroselect assay was considered positive if ORF1a (nsp3) and human RNase P (Internal Control) were detected. Negative, if ORF1a (nsp3) was not detected but RNase P was detected and invalid, if RNase P was not detected. Interpretation of STANDARD M nCoV Real-Time Detection kit (SD Biosensor, Gyeonggi-do, Republic of Korea) and LabGun (LabGenomics Co. Ltd., Gyeonggi-do, Republic of Korea) assays were as per manufacturers' instructions.

2.4 Interpretation for high throughput COBAS 6800 SARS-CoV-2 detection assay

As per the manufacturer, four types of interpretations were possible by the system (positive, negative, invalid, inconclusive). If target 1 (ORF), target 2 (E) and IC (non-infectious RNA in MS2 bacteriophage) were detected, the sample was interpreted as positive and if both targets were not detected but IC was detected, then the sample was interpreted as negative. If the samples showed target 1 negative, target 2 positive & IC positive, then the sample was interpreted as inconclusive.

3. Data analysis:

Contingency Tables was constructed and positive and negative concordance along with over all agreement between various assays was determined.

RESULTS

Development of novel in-house Viroselect assay: Primers were designed using the SARS-CoV-2 reference genome available at NCBI genbank (Accession number: NC 045512). The primers and probe were uniquely designed for Viroselect assay from regions which had maximum variation with the SARS-CoV-2 and MERS-CoV-2. This hypervariable region was identified through NCBI align2 blast (Figure 1a). As shown in Table 1 the sequence targets for designing the primers and probe were present in ORF1a (nsp3) region of the SARS-CoV-2 genome between nucleotide positions 1933-3956. The internal control primer and probe were derived from human ribonuclease (RNase) P gene (GenBank accession number NM_006413) sequence (Emery et al., 2004). As shown in Figure 1a, the unique design of the primer can aid in differentiating between SARS-CoV and SARS-CoV-2 during mixed infection.

NCBI nBlast was performed to analyze the sequence homology of putative amplicon with more than 250 SARS-CoV-2 isolates and 100% nucleotide sequence similarity was observed. Also, primers were designed to ensure that the putative amplicon did not overlap with known major mutation hotspots prevailing globally and in India (Figure 1b). Additionally, care was taken to have minimum homology with the human genome to avoid nonspecific amplification (Table 1). Following design of primers, amplicon was generated from infected samples and verified by agarose gel electrophoresis to be of the expected size (Figure 1c).

Concordance Analysis:

Viroselect with SD Biosensor assay

To assess the concordance of our in-house assay (Viroselect) with established and approved kits used for diagnosis, we evaluated the overall agreement of our assay with an EUA (USFDA) approved kit (SD Biosensor). RNA extraction from samples diagnosed as positive (n = 126) and negative (n = 107) for SARS-CoV-2 infection by SD Biosensor was followed by RT-PCR using Viroselect. Results revealed high concordance of the Viroselect assay with the SD Biosensor Kit as summarized in Table 2. As shown in Figure 2, a total of 286 samples (including 19 inconclusives) were tested of which 34 were invalid (not included in analysis). A concordance of 92.9% among positive samples and (81.3%) among negative samples was observed. An overall agreement of 87.6%, with 95% confidence interval (74.91% to 93.62%) was observed (Table 2). Interestingly a higher rate of discordance was observed in negative samples (20 out of 107; 18.7%) compared to that in positive samples (9 out of 126; 7.1%).

Notably, out of the 19 inconclusives by SD Biosensor kit, 7 (37%) were positive for Viroselect assay and 12 were negative. Of the 7 positives, six were positive when tested with another EUA (USFDA) approved RT-PCR kit, Labgun COVID-19 RT-PCR kit (LabGenomics). Further, clinical data obtained from the ICMR specimen referral form (SRF) revealed that six individuals out of seven that tested positive by Viroselect assay were symptomatic (Figure 2). Thus, Viroselect was able to resolve 100% of inconclusive samples tested by SD Biosensor.

Viroselect assay with COBAS 6800 High Throughput Platform

As our centre transitioned to a high throughput testing facility, we also had the opportunity to establish the concordance of Viroselect with SARS-CoV-2 detection assay on the COBAS 6800 High Throughput Platform (COBAS). Known COBAS positive (n= 102), negative (n=108) were tested by Viroselect. Seven of these samples were found to be invalid and not analysed further. Results revealed high concordance of the Viroselect assay with COBAS as summarized in Table 3. A concordance of 89% among positive samples and 84.5% among negative samples was observed with an overall agreement of 84.7%, with 95% confidence interval (76.00% to 90.85%) (Table 3). Interestingly, a higher rate of discordance was observed in negative samples (16 out of 103; 15.5%) compared to that in positive samples (11 out of 100; 11%).

As shown in Figure 3, a total of 205 inconclusive samples were also tested of which 13 were invalid (not included in analysis). Of these, 175 (91.2%) tested negative and 17 (8.9%) tested positive. This indicates that 192 out of 205 (93.7%) of the COBAS inconclusive samples were resolved.

DISCUSSION

In this study, we demonstrate the utility of using a unique, ORF-1a based SARS-CoV-2 detection assay in resolving inconclusive results obtained from a conventional manual assay as well for those obtained from a high-throughput testing platform.

Real time PCR has been established as a gold standard for detecting SARS-CoV-2 infection. The whole genome sequencing of SARS-CoV-2 in a very large number, and availability of free databases (including, NCBI, GISAID) for accessing the genomes have attributed to development of multitude of RT PCR based kits which can often provide conflicting and inconclusive results (Fang et al., 2020). Thus, a testing strategy that would resolve these issues would greatly enhance efficiency of diagnosis which in turn would significantly impact transmission dynamics and allow efficient epidemiological surveillance. Also, the rapidly accumulating mutations in SARS-CoV-2 genome pose a challenge for accurately diagnosing the infection due to possible mispriming in assays that rely on detection of viral sequences in these regions (Rambaut et al., 2020). Based on hypervariable sequence analysis as well as mutation hotspot enumeration we designed unique primers capable of differentiating infection with SARS-CoV-2, SARS-CoV and MERS-CoV. The hypervariability was observed in the regions (1933-3956) of ORF1a and our primer/probe set spanned within the nsp3 region. This design strategy was implemented to develop an in-house RT-PCR based assay for detection of SARS-CoV-2.

Our unique primer and probe design were subjected to comparative analysis with 1) WHO recommended and EUA (USFDA) approved SD Biosensor kits, as well as with 2) the high throughput, automated SARS-CoV-2 detection assay on the COBAS 6800 platform. Overall high concordance was observed for our inhouse assay with both these detection modalities. Viroselect assay was able to resolve 100% of the inconclusive results obtained by the SD Biosensor kit. Notably, 7 out of 19 (36.8%) of these samples were resolved to be positive by our assay, while 6 out of 7 were resolved to be positive by another EUA (USFDA) approved

kit, LabGun. Additionally, when corroborated clinically, 6 out of 7 inconclusive samples tested positive by Viroselect were documented to have been obtained from symptomatic individuals.

A precise and early diagnosis and treatment is crucial in COVID-19 (Rong et al., 2020). High-throughput diagnostic platforms, such as the COBAS 6800 are critical in addressing diagnostic demand during fast moving epidemics such as COVID-19 and have a direct impact on curtailing transmission as well as clinical management. Considering the similarity in geography, disease prevalence, population dynamics and healthcare spending of all the countries in South Asia, Ministry of Health and Family Welfare, Government of India adopted commendable measures to tackle the COVID-19 pandemic. A major such measure was increased accessibility to high throughput testing facilities. Our laboratory was one of these centres where we have tested close to 100,000 samples to date from Mumbai and surrounding areas. When dealing with such high throughput, a major challenge encountered was the rate of inconclusives. Although the rate was relatively low at 3%, the absolute number translated into approximately 3000 individuals (for every 100,000 individuals). Considering rapid transmission dynamics of SARS-CoV-2, even a day's delay in correctly diagnosing and reporting a false negative or inconclusive sample as positive would result in significant impact on disease transmission in a highly populated and resource setting such as India (Kretzschmar et al., 2020; Larremore et al., 2021; Rong et al., 2020). When Viroselect was used to interrogate a subset of inconclusive samples obtained from our COBAS 6800 platform we found that our assay was able to successfully resolve ~94% of these samples. Most importantly, ~9% of these inconclusives were resolved to be positive. Our assay would thus enable optimal intervention in managing these cases, especially if they are assumed to be potential transmitters.

Consent for publication

Not applicable

Availability of data and materials

Not applicable

Competing interests

The authors declare that they have no competing interests.

Funding sources

This work was supported by Indian Council of Medical Research- National Institute for Research in Reproductive Health (ICMR- NIRRH). The funding body had no role in the design of the study, collection, analysis, and interpretation of data or in writing the manuscript.

Authors' contributions

DJ, SB1, and VP conceptualized the study. KJ, AYS, BC, SK, NK, SB2, SB3, PD, KK, HP performed the experiments. KJ, AYS, BC conducted primary data analysis and prepared the figures. JS, SA, KJ, BC and KM collected the clinical data. DJ, SB1, and VP supervised the study. VP and DJ wrote and edited the manuscript; all authors discussed the results and contributed to the final manuscript.

Acknowledgments

We are grateful to Director General, Indian Council of Medical Research (ICMR), Ministry for Health & Family Welfare, Government of India for encouraging and supporting with resources, the pursuit of important research questions relevant to COVID-19 pathogenesis.

Ethics approval

This study was approved by Institutional Ethics Committee at ICMR-NIRRH (Project No. 399-2020)

REFERENCES:

Bhattacharya S, Vidyadharan A, Joy V. Inconclusive SARS-COV-2 reverse transcription-polymerase chain reaction test reports: Interpretation, clinical and infection control implications. *J Acad Clin Microbiol* 2020;22:59–61.
https://doi.org/10.4103/jacm.jacm_16_20.

Emery SL, Erdman DD, Bowen MD, Newton BR, Winchell JM, Meyer RF, et al. Real-Time Reverse Transcription–Polymerase Chain Reaction Assay for SARS-associated Coronavirus. *Emerg Infect Dis J* 2004;10:311. <https://doi.org/10.3201/eid1002.030759>.

Fang FC, Naccache SN, Greninger AL. The Laboratory Diagnosis of Coronavirus Disease 2019— Frequently Asked Questions. *Clin Infect Dis* 2020;71:2996–3001.
<https://doi.org/10.1093/cid/ciaa742>.

Gupta N, Bhatnagar T, Rade K, Murhekar M, Gangakhedkar RR, Nagar A. Strategic planning to augment the testing capacity for COVID-19 in India. *Indian J Med Res* 2020;151:210–5. https://doi.org/10.4103/ijmr.IJMR_1166_20.

Hur K-H, Park K, Lim Y, Jeong YS, Sung H, Kim M-N. Evaluation of Four Commercial Kits for SARS-CoV-2 Real-Time Reverse-Transcription Polymerase Chain Reaction Approved by Emergency-Use-Authorization in Korea. *Front Med* 2020;7:521.
<https://doi.org/10.3389/fmed.2020.00521>.

ICMR COVID-19 Testing Status. Delhi: Indian Council Of Medical Research, 2021.
Available
online:https://www.icmr.gov.in/pdf/covid/update/ICMR_testing_update_29Jan2021.pdf
(last cited: 29th January 2021)

Kretzschmar ME, Rozhnova G, Bootsma MCJ, van Boven M, van de Wijgert JHHM, Bonten MJM. Impact of delays on effectiveness of contact tracing strategies for COVID-19: a modelling study. *Lancet Public Heal* 2020;5:e452–9. [https://doi.org/10.1016/S2468-2667\(20\)30157-2](https://doi.org/10.1016/S2468-2667(20)30157-2).

Larremore DB, Wilder B, Lester E, Shehata S, Burke JM, Hay JA, et al. Test sensitivity is secondary to frequency and turnaround time for COVID-19 screening. *Sci Adv* 2021;7:eabd5393. <https://doi.org/10.1126/sciadv.abd5393>.

Poljak M, Korva M, Knap Gašper N, Fujs Komloš K, Sagadin M, Uršič T, et al. Clinical Evaluation of the cobas SARS-CoV-2 Test and a Diagnostic Platform Switch during 48 Hours in the Midst of the COVID-19 Pandemic. *J Clin Microbiol* 2020;58:e00599-20. <https://doi.org/10.1128/JCM.00599-20>.

Pujadas E, Ibeh N, Hernandez MM, Waluszko A, Sidorenko T, Flores V, et al. Comparison of SARS-CoV-2 detection from nasopharyngeal swab samples by the Roche cobas 6800 SARS-CoV-2 test and a laboratory-developed real-time RT-PCR test. *J Med Virol* 2020;92:1695–8. <https://doi.org/10.1002/jmv.25988>.

Rambaut A, Loman N, Pybus O, Barclay W, Barrett J, Carabelli A, Connor T, Peacock T, Robertson DL, Volz E. Preliminary genomic characterisation of an emergent SARS-CoV-2 lineage in the UK defined by a novel set of spike mutations. 2020. Available online: <https://virological.org/t/preliminary-genomic-characterisation-of-an-emergent-sars-cov-2-lineage-in-the-uk-defined-by-a-novel-set-of-spike-mutations/563>

Rong X, Yang L, Chu H, Fan M. Effect of delay in diagnosis on transmission of COVID-19. *Math Biosci Eng* 2020;17:2725–40.

WHO coronavirus disease (COVID-19) dashboard. Geneva: World Health Organization, 2021. Available online: <https://covid19.who.int/> (last cited: 28th January 2021)

WHO Feature Story: How India scaled up its laboratory testing capacity for COVID-19. Geneva: World Health Organization, 2020. Available online: <https://www.who.int/india/news/feature-stories/detail/how-india-scaled-up-its-laboratory-testing-capacity-for-covid19>

Yang Y, Peng F, Wang R, Yange M, Guan K, Jiang T, et al. The deadly coronaviruses: The 2003 SARS pandemic and the 2020 novel coronavirus epidemic in China. *J Autoimmun* 2020;109:102434. <https://doi.org/10.1016/j.jaut.2020.102434>.

FIGURES LEGEND

Figure 1a: Depiction of the Hypervariable region between the genomes of SARS-CoV-2, SARS-CoV and MERS-CoV

Figure 1b: Depiction of primer position relative to major mutation hotspots of ORF1ab gene of SARS-CoV-2.

Figure 1c: RT-PCR amplicon separated by 2% agarose gel electrophoresis after RNA extraction and 1st strand cDNA synthesis. RNA extraction was carried out using QIAamp Viral RNA Mini Kit by QIAGEN and 1st strand cDNA synthesis was carried out using Takara Prime Script 1st strand cDNA synthesis kit. PCR of the cDNA product was carried out using GeNei Taq DNA Polymerase and using the designed ORF1a forward and reverse primers. The expected amplicon size is 170bp which is observed in the figure.

Figure 2: Flowchart showing Viroselect test results of the selected sample panel tested on SD Biosensor

Figure 3: Flowchart showing Viroselect test results of the selected sample panel tested on COBAS 6800 platform.

Figure 1a: Depiction of the Hypervariable region between the genomes of SARS-CoV-2, SARS-CoV and MERS-CoV




Figure 1b: Depiction of primer position relative to major mutation hotspots of ORF1ab gene of SARS-CoV-2.





Figure 1c: RT-PCR amplicon separated by 2% agarose gel electrophoresis after RNA extraction and 1st strand cDNA synthesis. RNA extraction was carried out using QIAamp Viral RNA Mini Kit by QIAGEN and 1st strand cDNA synthesis was carried out using Takara Prime Script 1st strand cDNA synthesis kit. The PCR of the cDNA product was carried out using GeNei Taq DNA Polymerase kit and using the designed ORF1a forward and reverse primers.



The expected amplicon size is 170bp which is observed in the figure.

Lane 1: 100bp DNA Ladder (Marker)

Lane 2: Blank Well

Lane 3: Reagent Control

Lane 4: COVID-19 Positive sample

Lane 5: COVID-19 Negative sample

Figure 2: Flowchart showing Viroselect test results of the selected sample panel tested on SD Biosensor.





Figure 3: Flowchart showing Viroselect test results of the selected sample panel tested on COBAS 6800 platform.



TABLES LEGEND

Table 1: Primer sequences for SARS-CoV-2 detection (RT-PCR). Homology of ORF1a (*nsp3*) RT-PCR primers and Probe with Human Genome sprovided using nBLAST ,highly similar sequences (megablast).

Seq ID.	Primer Name	Primer Sequence (5' – 3')	Map Location	Strand binding	Homology with Human (QC+ID)	e Value	Amplicon Size	Description
1	ORF1 (nsp3) QPCR_F P	TTCTGCT GCTCTTC AACCTGA	3160-3180	(-)	80%QC,1 00% ID	1.3	174	Unique to SARS-CoV-2
2	ORF1 (nsp3) QPCR_R P	ATAGTCT GAACAAC TGGTGTA AGT	3310-3330	(+)	79%QC,1 00% ID	2.2		
3	ORF1_P robe (TaqMan)	GTTGGTC AACAAGA CGGCAGT GAGGACA ATC (5' Tamra- 3'BHQ1)	3230-3257	(-)	No significant similarity	NA		
4	Human RNase P QPCR_F P	AGATTTG GACCTGC GAGCG	Chromosome 10: 90872002-90872020	(-)				(Emery et al., 2004)
5	Human RNase P QPCR_R P	GAGCGGC TGTCTCC ACAAGT	Chromosome 10: 908722047-90872066	(+)				

6	Human RNase P_ Probe	TTCTGAC CTGAAGG CTCTGCG CG	Chromosome 10: 90872023- 90872045	(-)				
---	----------------------------	-------------------------------------	-----------------------------------------	-----	--	--	--	--

QC: Query coverage, ID: Percentage Identity

Table 2: Concordance of SD Biosensor positive and negative samples with Viroselect assay

Viroselect results	SD Biosensor			Concordance among positives	Concordance among negatives	Overall agreement (%) (95% CI)
	Positive	Negative	Total			
Positive	117	20	137	92.85%	81.31%	87.55% (74.91 to 93.62)
Negative	9	87	96			
Total	126	107	233			

Table 3: Concordance of COBAS positive and negative samples with Viroselect assay

Viroselect results	COBAS results			Concordance among positives	Concordance among negatives	Overall agreement (%) (95% CI)
	Positive	Negative	Total			
Positive	89	16	105	89%	84.47%	84.7% (76.00% to 90.85%)
Negative	11	87	98			
Total	100	103	203			

Declaration of interests

The authors declare that they have no known competing financial interests or personal relationships that could have appeared to influence the work reported in this paper.

The authors declare the following financial interests/personal relationships which may be considered as potential competing interests:



map2 mutations: T209I, Y378I, G293D, F347R, D448Q, L737V, P756S

map4 mutations: R576I, T330I, T229V, R370R

map8 mutation: F937Y

map5: G378S

map5: G335D

map5: L392I

map12: M455S

map12: T486V

map12: T540R

map13: V556L

map13: P570R

map13: R685V

map13: S702I

Map15: I925I

map5/10r: G668W

● Mutation

○ Promoter

— Amplification

